


## **General Disclaimer**

### **One or more of the Following Statements may affect this Document**

- This document has been reproduced from the best copy furnished by the organizational source. It is being released in the interest of making available as much information as possible.
- This document may contain data, which exceeds the sheet parameters. It was furnished in this condition by the organizational source and is the best copy available.
- This document may contain tone-on-tone or color graphs, charts and/or pictures, which have been reproduced in black and white.
- This document is paginated as submitted by the original source.
- Portions of this document are not fully legible due to the historical nature of some of the material. However, it is the best reproduction available from the original submission.



NASA CR-

147828

# TECHNOLOGY

RESEARCH

and

INCORPORATED

LIFE SCIENCES DIVISION

DEVELOPMENT

SPECIAL REPORT

EMGAN: A Computer Program

for Time and Frequency Domain Reduction of Electromyographic Data

(NASA-CR-147878) EMGAN: A COMPUTER PROGRAM  
FOR TIME AND FREQUENCY DOMAIN REDUCTION OF  
ELECTROMYOGRAPHIC DATA (Technology, Inc.,  
Houston, Tex.) 165 p HC \$6.75 CSCL 09B

N76-31936

Unclas

G3/61 03433

September 5, 1975

CONTRACT NAS 9-13291

National Aeronautics and Space Administration  
Lyndon B. Johnson Space Center  
Houston, Texas 77058





TECHNOLOGY INCORPORATED  
LIFE SCIENCES DIVISION  
HOUSTON, TEXAS

SPECIAL REPORT

EMGAN: A Computer Program


for Time and Frequency Domain Reduction of Electromyographic Data

September 5, 1975


SPECIALIZED CARDIOVASCULAR STUDIES

CONTRACT NAS 9-13291


PREPARED BY:

  
William N. Hursta  
Research Bioengineer

APPROVED BY:

  
Joseph T. Baker  
Project Leader

APPROVED BY:

  
T. Wayne Holt, Manager  
Life Sciences Division

SPECIAL REPORT

EMGAN: A Computer Program  
for Time and Frequency Domain Reduction of Electromyographic Data

September 5, 1975

CONTRACT NAS 9-13291

National Aeronautics and Space Administration  
Lyndon B. Johnson Space Center  
Houston, Texas 77058

## TABLE OF CONTENTS

INTRODUCTION	1
DATA PREPARATION	2
PROGRAM CHARACTERISTICS	3
EMGAN DESCRIPTION	5
SAMPLE EMGAN RUN	13
PROGRAM IMPROVEMENTS	14
REFERENCES	15
APPENDIX A - ASTP Myography Description	
APPENDIX B - Input/Output Magnetic Tape and Data Card Formats	
APPENDIX C - Main Program Flow Chart	
APPENDIX D - Table of Formulas	
APPENDIX E - Example of EMGAN Input Card Deck with Plotted and Printed Output	
APPENDIX F - Sigma Source Listing and Load Map	
APPENDIX G - Nova Source Listing and Load Map	
APPENDIX H - CF-16/A General A/D Program Documentation	



## INTRODUCTION

An experiment in Electromyography (EMG) utilizing surface electrode techniques was developed for the Apollo-Soyuz Test Project. The objective of the study was to obtain quantitative measure of changes in muscle characteristics as a result of exposure to weightlessness (see Appendix A). This report describes the computer program, EMGAN, which was written to provide first order data reduction for the experiment.

EMG signals are produced by the membrane depolarization of muscle fibers during a muscle contraction. Surface electrodes detect a spatially summated signal from a large number of muscle fibers commonly called an interference pattern. An interference pattern is usually so complex that analysis through signal morphology is extremely difficult if not impossible. However, because of the ease of application and noninvasive character of surface electrodes and the often desirable feature of studying large groups of muscle fibers, methods have been sought to evaluate the surface interference pattern.

In the time domain a valuable technique to characterize an interference pattern has been the root-mean-square (RMS) value of the signal. A similar technique has been to integrate the interference pattern over a given time length. A linear relationship is usually found between the integrated EMG or the RMS value and the force exerted by the muscle up to some maximum.

More recently with the increased access to the digital computers and the availability of an algorithm to economically calculate discrete Fourier transforms, it has

become common to process EMG interference patterns in the frequency domain.

Muscle fatigue and certain myopathic conditions can be recognized through changes in muscle frequency spectra.

## DATA PREPARATION

The protocol for the ASTP Myography required that two types of data be recorded, digitized, and processed: the force being exerted by a muscle group and the surface EMG signal. The data was originally recorded on a 4-channel Tandberg analog recorder with one-channel of IRIG B time code, one-channel of force data, and two-channels of simultaneous EMG signal (two different muscles).

A digital tape was prepared from the analog data using the General A/D program (see Appendix H) available on the CF-16A mini-computer located in the Sigma 3 Computer Room. Analog force data was low pass filtered at 10 Hz before being digitized at 20 samples/second. Analog EMG data was low pass filtered at 400 Hz. Consistent with this frequency limit a digitizing rate of 1000 samples/second was chosen as the minimum necessary for power spectral density analysis [2].

The digitized data was organized into half second records with 535 words per record. The first twenty-five words were header information with ten words of digitized force data interspaced into 500 words of EMG data (see data format in Appendix B). Only one-channel of EMG data was digitized with its associated force data at any one time. Although not required by EMGAN, a data tape was usually constructed with the first file containing force calibration and subsequent



files containing EMG data from each muscle with the appropriate EMG calibration signals.

## PROGRAM CHARACTERISTICS

EMGAN exists in two versions. One version was written to execute on the Biomedical Sciences Division Xerox Sigma 3 Computer using the Computek 400/20 CRT terminal for graphical output. A second version was written to execute on the Cardiovascular Laboratory's Data General Nova 820 mini-computer. Both programs are essentially identical to the user. However, there are a number of programming differences made necessary by differences in Xerox and Data General's Fortran IV, supporting software packages (graphical and tape I/O), and peripheral hardware. Whenever these differences are sufficiently important in describing the program they will be noted. Complete program listings of both versions with their associated load maps are given in the Appendices F and G.

An EMG system tape exists for the Sigma 3 Computer which contains an appropriately sized operating system. EMGAN exists in absolute binary as a file in the User Processor area of the Rad. A source listing of EMGAN exists both on magnetic tape stored in the Sigma 3 Computer Room and on punched cards in the Cardiovascular Laboratory. Both source listings and absolute binary files of the Nova version of EMGAN are stored on disk pack and on magnetic tape in the Cardiovascular Laboratory. The absolute binary files are named EMGAN.SV and EMGAN.OL (both are needed for program execution). The source listing is lo-

cated in separate files for the main program and each subroutine with each file named for the subroutine it contains. How to execute a run with either computer version is explained in a later section. The computer peripherals required to fully execute EMGAN are a card reader, line printer, two 9-track tape drives, and the Computek terminal (Sigma) or the Tektronix 4014 terminal (Nova).

EMGAN's execution is totally controlled by the data card information received. Force and EMG data slices are processed with separate data cards. The maximum length of a force data slice is 120 seconds with previous force calibration required. After the raw data is converted into pounds it is averaged over a user selected interval for output to the line printer. The maximum length of an EMG data slice is four seconds with previous EMG calibration required. After scaling into microvolts, an integrated value is found for the EMG data slice. The discrete Fourier transform of the data is taken. The raw power spectral density (PSD) of the data up to 400 Hz is calculated and smoothed over a frequency bandwidth selected by the user. The raw PSD with associated header information is outputted to magnetic tape for use in any subsequent processing program. The smoothed PSD is outputted to the line printer along with accessory information such as a normalized PSD, percent contribution of each bandwidth to the whole spectrum, the mean and standard deviation of the spectrum, and time domain integration and mean square value of the data slice.

EMGAN run time is variable and depends on the number and kind of data cards, the length of the EMG data slices, whether or not the 60 Hz filter is used,

data location on the input tape (processing data in sequential order on the tape is much faster than skipping back and forth on the input tape), the number of data plots, and whether execution is on the Sigma or Nova. For the ASTP EMG experiments a typical run consisting of calibration cards, 16 force data slice cards, and 28 4-second EMG data slice cards (no filtering or plotting) required approximately 25 minutes on the Sigma and 35 minutes on the Nova.

## EMGAN DESCRIPTION

In describing the various subroutines, differences in the Nova version from the Sigma version will be set off with a \*\*\*. A table of formulas used in the program is listed in Appendix D.

### Main Program:

The Main Program's basic task is to load the needed overlays in the proper sequence. A flowchart is given in Appendix C. The program is organized into a large do loop with one pass through the loop for each data card to be processed. A check is made to determine if calibration data has been received before the respective data type, force or EMG, can be processed. If the user should request the processing of data without calibration, the program outputs an error message and halts.

\*\*\* The magnetic tape units are initialized and the output tape file is opened. Upon completion of the program the tape units are released.



### Subroutine CINPUT

Subroutine CINPUT inputs all data card information. Data card format is given in Appendix B. As each card is read it is outputted to the line printer to aid the user in verification. If the number of data cards to be processed is read to be greater than 200, an appropriate error message is given and the program halts. All data card information is stored in a scratch file, ICARD, in the user data area of the RAD. As each data card is to be processed ICARD is read back into main memory. The program variables are then set to the corresponding data card values.

\*\*\* The scratch file, ICARD, is created on the disk pack by the Nova program if it does not already exist. For the Sigma version the user should create ICARD (1 record, 6000 bytes, random format) if it does not already exist before attempting to execute EMGAN.

### Subroutine TINPUT

Subroutine TINPUT acquires the digitized data from the input tape. Tape files are numbered from 00 to 99. The routine keeps track of the current file being accessed. If a data card requests another file, TINPUT skips the tape either forward or backward to the beginning of the requested file. Subroutine FIND is called to position the tape at the start of the data slice within the file. TINPUT reads either force or EMG data according to data card request. After the data has been read, the tape is backed up to the beginning of the data slice.

If an EOT is encountered during a search for a requested file, an error message is given and the program halts. If an EOF or EOT is encountered during the input of a data slice, a error message is given and the program halts.

\*\*\* Nova Fortran calls for magnetic tape I/O are considerably different from the QINOUT package used on the Sigma. Because Nova error checking is more comprehensive, error messages given by EMGAN resulting from tape I/O problems contain much more information.

#### Subroutine FIND

Subroutine FIND searches the currently accessed tape file for the start time of the data slice. If an EOF is encountered before the start time is found, the tape is rewound to the beginning of the file. A second search is then made through the entire file. If the start time is still not found, an error message is given and the program halts.

#### Subroutine OUT60

Subroutine OUT60 is a digital notch filter which can be called to remove 60 Hz and its odd harmonics from EMG data. If at all possible filtering should be avoided, since no distinction can be made between 60 Hz present in the true data and that which may have been picked up from extraneous sources. Since a large number of sines and cosines are generated, considerable time is added to the processing of data slices.



### Subroutine DCAL

Subroutine DCAL processes force and EMG calibration data to obtain the scale factors needed to change the data from raw counts to engineering units. After the scale factors have been calculated, the calibration is converted to engineering units for plotting.

A force calibration data slice consists of two DC levels, the first level corresponding to zero pounds of force, and a second level corresponding to a known amount of force (98 pounds in the ASTP experiment). The data slice is searched for a jump in the data to indicate the beginning of the second level. When the jump is found, the average number of counts three seconds before the jump is subtracted from the average number of counts three seconds after the jump. The difference is divided by the number of pounds represented by the jump to obtain the scale factor. If the jump is not detected, an error message is outputted and the program halts.

The EMG scale factor is calculated using the fact that for a zero mean, stationary signal the standard deviation of the signal is equal to its RMS value. (For the ASTP experiment a nominal 20 Hz, 350 microvolt RMS signal was used for EMG calibration). The calibration data first has any DC mean subtracted from it and the standard deviation of the data in counts is found. The scale factor is computed by dividing the standard deviation by the RMS magnitude of the calibration signal.

### Subroutine GRAPH

Subroutine GRAPH plots on the Computek terminal, if requested by the data card, force and EMG calibration data, force data, and EMG data both in the time and frequency domain. Scaling for force and EMG data is variable according to the maximum amplitude of the data. Force data is plotted without smoothing (such as is done for the printed output). When time domain EMG data is plotted, only 1000 evenly spaced points are plotted no matter what the length of the data slice. The PSD plotted is smoothed and normalized to the maximum value of the spectrum. The return key must be struck for the program to continue whenever the terminal's bell is rung.

\*\*\*As of the writing of this report the graphics package for the Tektronix 4014 terminal is not available. No plotting capability presently exists when executing EMGAN on the Nova. However, if plots are requested when executing on the Nova, no run error occurs.

### Subroutine MODLINE, WORDS, and BELL

These routines refer to the Computek terminal and respectively, plot data arrays, output alphanumeric characters, and ring the Computek bell.

### Subroutine EMG

Subroutine EMG subtracts any DC offset from the unscaled EMG data. The data is then scaled into engineering units and the largest absolute value of the

data array is found for use in plotting of the data. Simpson's rule is used to find the integrated value of the data in units of microvolt-seconds. The mean square value of data is calculated.

#### Subroutine FORCE

Subroutine FORCE scales the force data into engineering units. Compensation is made for force signals recorded at a higher gain than the force calibration. The first second of force data is assumed to be the zero force level. The user should insure that this is so for correct results. The maximum force in the data slice is found for use in plotting the data. Force data is then averaged over intervals requested by the data card for output to the line printer.

#### Subroutine FFTPSD

Subroutine FFTPSD calculates the raw power spectral density of an EMG data slice up to 400 Hz. It is suggested that the user read (2) for information on data analysis using PSD techniques and (1) for an excellent explanation of discrete Fourier transforms and the fast Fourier transform algorithm. A cosine taper is applied to the first and last tenths of the data slice to reduce leakage. The Fourier transform of the data is then taken and the raw power spectral density is computed. The PSD amplitudes are corrected for reduction caused by the cosine taper. The raw PSD is outputted to magnetic tape for use in subsequent processing programs. The PSD is smoothed over the frequency bandwidth requested by the data card and the maximum PSD value found. A second PSD array is computed,



normalized to the maximum PSD value. The area under the PSD is found and used to calculate the percentage each bandwidth contains of the total power and the cumulative of percentage of total power with increasing frequency. The expected value and variance of the PSD are computed.

Viewed statistically the raw PSD has very poor reliability. Each calculated value is an inconsistent estimate of the true value and has a possible random error of 100%. Smoothing the PSD can greatly reduce the random error and provide a better estimate over the smoothed bandwidth (2). The normalized standard error is calculated to provide information on the amount of random error associated with each PSD value.

#### Subroutine WINDOW

Subroutine WINDOW applies a cosine taper to the first and last tenths of the data slice to reduce leakage in the discrete Fourier transform.

#### Subroutines RFORT and FORT

Subroutine RFORT and FORT take the Fourier transform of the EMG data. FORT contains the actual fast Fourier transform algorithm and assumes a transform of  $N$  complex data points. RFORT compensates for the actual case of a transform of  $2N$  real data points. The FORT computes an initial sine/cosine array and refers to this array when computing the transform. Repetitive computation of identical sines and cosines is thereby avoided, substantially increasing the speed of the transform

at the cost of an extra array of size  $N/4$ . Both of these subroutines were received from Mr. Jack McBryde of Lockheed. The only modification was to not divide each data value by the number of data points (an alternate definition of the discrete Fourier transform). EMGAN can transform 1, 2, or 4 seconds of EMG data which corresponds to 1024, 2048, and 4096 data points, respectively. Computation time for the Fourier transform algorithm are as follows:

<u># of Points</u>	<u>Sigma Time</u>	<u>Nova Time</u>
1024	6 seconds	10 seconds
2048	13 seconds	22 seconds
4096	29 seconds	47 seconds

\*\*\* Subroutines RFORT and FORT are named RFFT and FFT to avoid confusion with Nova Command Line Interpreter calls.

#### Subroutine PSDSAVE

Subroutine PSDSAVE outputs to magnetic tape the raw PSD of each EMG data slice. The PSD values are preceded by a 25-word header record. The format of the header record is located in Appendix B. The PSD of a one-second slice of data has 512 real values, a two-second slice has 1024 real values, and a four-second slice has 2048 real values. A PSD record contains 2048 integer words; therefore, a one-second PSD fills  $1/2$  of one record, a two-second PSD fills one full record, and a four-second slice fills two records. Two EOF's are written on the tape after the raw PSD is output and the tape then backs up to just before the



EOF's. At the end of an EMGAN run the output tape consists of a file containing all of raw PSD's from the run terminated by two EOF's.

### Subroutine PRINT

Subroutine PRINT outputs to the line printer all of the processed force and EMG data. A page of header and calibration information is printed whenever the print switch on a calibration data card is set. The length of a force or PSD output is checked and if sufficiently long a two column list of the processed data is printed instead of a one column list. The frequency corresponding to each PSD value is the center of the bandwidth.

\*\*\* Subroutine PRINT is named POUT to avoid confusion with Nova Command Line Interpreter calls.

### SAMPLE EMGAN RUN

To execute an EMGAN run on the Sigma 3 first load the EMG system tape. Place the data cards in the card reader and power up the Computek terminal. Mount the input data tape on unit 0 and the PSD output on unit 1. Assign the user input device to the teletype and type in "!EMGAN".

To execute an EMGAN run on the Nova load the EMG disk pack and bring up the system. Place the data cards in the card reader. Mount the input data tape on unit 0 and the output tape on unit 1. Type in "EMGAN". A message will be outputed to the terminal, "LOAD \$CDR, STRIKE ANY KEY", upon entering any keyboard character, the data cards will be read and the program

executed.

An example of an EMGAN input card deck with the printed and plotted output is in Appendix E.

## PROGRAM IMPROVEMENTS

There are several additions that would be desirable in future work with EMGAN. The major addition would be a subroutine to check the EMG data slice prior to transform for stationarity and normal distribution. Suggestions for the necessary statistical techniques are found in (2).

An alternative FFT algorithm might be used which calculates sines and cosines as needed rather than establishing a table. The core storage saved would allow an 8196 real data point transform to be performed on the Sigma and, with proper sizing of operating system, on the Nova. However a price would be paid in increased processing time considerably beyond twice the transform time of 4096 real data point array.

Of a more minor nature subroutine GRAPH could be modified to output the muscle name with the EMG and PSD plots to allow easier later identification. Although force and EMG data are now processed separately, it could be useful to obtain the average force value over the period of the EMG data slice. The force value could be output with the EMG's PSD for more convenient association between the two.

## REFERENCES

1. Brigham, E. O. The Fast Fourier Transform, Prentice-Hall Inc., Englewood Cliffs, N. J., 1974.
2. Bendat, J. S. and Piersol, A. G. Random Data: Analysis and Measurement Procedures, John Wiley and Sons, New York, N. Y., 1971



APPENDIX A

ASTP MYOGRAPHY DESCRIPTION

## 1.0 General Objective

The general objective of the experiment is to continue the study efforts begun in the Skylab program to identify and describe antigravity muscle dysfunction characteristics and consequences resulting from spaceflight.

## 2.0 Specific Objective

The purpose of the experiment is to assess changes that occur following a period of disuse. Some studies suggest that the first few days of exposure to 0-g may be significant to the ultimate muscle deconditioning resulting from longer missions. The relationship between muscle capability, in terms of strength or tension, fatigability, and muscle electrical activity will be investigated, as well as the differential effects of spaceflight disuse on "fast" and "slow" muscles.

## 3.0 Significance

Muscle function and condition may well be a critical determinant of man's capacity to endure the effects of long duration weightlessness as well as the readaptation to the earth's gravity, or to the gravity of other planets. Data collected in the experiment will aid in quantifying the muscle deconditioning which results from weightlessness.

In addition, an important spinoff will be the extension of an already considerable ground-based body of knowledge about the characteristics and consequences of muscle disuse. Heretofore, research in this area has relied on the "contrived" methods of surgical section and limb or torso fixation to produce the disuse effects. Spaceflight uniquely provides a "pure" form of neuromuscular system disuse.

## 4.0 Method

### A. Concept

The measurement of muscle electrical activity is known as electromyography (EMG). The state of muscle function can be described by EMG measurements combined with knowledge of the force being exerted by the muscle.

Electromyographic studies have shown that skeletal muscle undergoes changes in capability and composition when subjected to periods of disuse, i.e., periods of time when the contractile mechanisms are not subjected to the stresses and forces normally encountered. Also changes in biochemical constituents, as seen in



Skylab, such as calcium, potassium, etc., and alterations in enzymic constituents such as ATP, and acetylcholinesterase have been shown, by ground studies to affect normal muscle function.

The EMG results of Skylab 3 skeletal muscle assessment provide ample evidence that normal muscle function is significantly altered by periods of weightlessness (disuse) of 56 days or more.

To investigate the effects of a shorter period of weightlessness on muscle function, a standardized test protocol and a muscle stress device will be used. The muscle stress apparatus will provide for preplanned isometric muscle forces from the calf muscles and the arm muscles. The standardized test protocol will include measures of muscle strength, muscle endurance, and muscle fatigability. Pre- and postflight measures only will be taken.

#### B. Procedure

The test procedure is identified in Figure 1. Two muscles each from the leg and arm will be instrumented and the muscle action potentials recorded. These muscles are the brachial bicep, brachioradialis, gastrocnemius and soleus. The procedure will require about 11 minutes to complete.

#### C. Data

The EMG data will be recorded on magnetic tape for time and frequency domain analysis at the JSC Cardiovascular Laboratory.

## PROCEDURE

- A. Apply surface electrodes
- B. Seat subject in muscle stress device prepared for muscle stress test
- C. Three short (1 to 2 seconds) efforts to determine the maximum voluntary contraction (MVC)
- D. Remainder of calf procedure
  - 1. 10 seconds at 10% MVC - 20 seconds rest
  - 2. 10 seconds at 20% MVC - 20 seconds rest
  - 3. 10 seconds at 30% MVC - 20 seconds rest
  - 4. 60 seconds at 50% MVC
- E. Reset muscle stress device for arm stress test
- F. Three short (1 to 2 seconds) efforts to determine the maximum voluntary contraction
- G. Remainder of arm procedure
  - 1. 10 seconds at 10% MVC - 20 seconds rest
  - 2. 10 seconds at 20% MVC - 20 seconds rest
  - 3. 10 seconds at 30% MVC - 20 seconds rest
  - 4. 60 seconds at 50% MVC
- H. Remove surface electrodes

APPENDIX B

INPUT/OUTPUT MAGNETIC TAPE AND DATA CARD FORMATS



EMG DATA INPUT TAPE RECORD

1	2	3	4	5	6	7	8	9	10	11	12
Record #	Hours	Minutes	Seconds	Mili-Seconds	Subject No.	Exp. Month	Exp. Day	Exp. Year	Digital Month	Digital Day	Digital Year

13	14	15	16	17	18	19	20	21	22	23	24
Flight Ref. Day	Run No.	Analog Tape No.	EMG Samp Rate (S/S)	Force Samp Rate (S/S)	RMS Amplitude EMG Cal	Not Used	Amount of Force Cal (Lbs)	Not Used	Not Used	Not Used	Not Used

25	26	27	28	29				78	79		
Not Used	EMG Data Samp	Force Data Sample	EMG Data Sample	EMG Data Sample				Force Data Sample	EMG Data Sample		

	129	130			180	181			231	232	
	Force Data Sample	EMG Data Sample			Force Data Sample	EMG Data Sample			Force Data Sample	EMG Data Sample	

	282	283			333	334			334	335	
	Force Data Sample	EMG Data Sample			Force Data Sample	EMG Data Sample			Force Data Sample	EMG Data Sample	

	435	436			486	487				534	535
	Force Data Sample	EMG Data Sample			Force Data Sample	EMG Data Sample				EMG Data Sample	EMG Data Sample

## HEADER RECORD FOR PSD OUTPUT

Word #

- 1 - Subject Number
- 2 - Experiment Month
- 3 - Experiment Day
- 4 - Experiment Year
- 5 - Digitizing Month
- 6 - Digitizing Day
- 7 - Digitizing Year
- 8 - Flight Reference Day
- 9 - Run Number
- 10 - Analog Tape Number
- 11 - Sample Rate for EMG Signal (Samp/sec)
- 12 - Sample Rate for Force Signal (Samp/sec)
- 13 - RMS Amplitude of Sine Wave Cal (in microvolts)
- 14 - Not Used
- 15 - Amount of Force Cal
- 16 - Start Hour
- 17 - Start Minute
- 18 - Start Second
- 19 - Length of Data Slice (sec)
- 20 - Muscle ID Number
- 21 - # of Following Records with PSD Values
- 22 - 25 - Not Used



## DATA CARD FORMAT

### First Card In Data Card Deck

cc: 1-5            Number of data slices (max = 200)

### Data Slice Cards

cc: 1-3            File data slice located in (00 to 99)

cc: 4-6            Data Type  
                    1 = Force Cal  
                    2 = EMG Cal  
                    3 = EMG Data  
                    4 = Force Data

cc: 7-9            HR start time of data slice

cc: 10-12          Min start time of data slice

cc: 13-15          Sec start time of data slice

cc: 16-18          Filter switch  
                    -1 = No filter  
                    0 = Filter 60 Hz  
                    1 = Filter 60, 180 and 300 Hz

cc: 19-21          Length of data slice in seconds  
                    Force data or cal - any integer number up to 120 seconds  
                    EMG Cal - 1, 2, 3 or 4 seconds  
                    EMG Data - 1, 2, or 4 seconds

cc: 22-24          Plot switch  
                    0 = Plot data  
                    1 = Suppress plotting

cc: 25-27          Print switch  
                    0 = Print data  
                    1 = Suppress printing

cc: 28-30

Averaging interval (force and EMG data cards only)

Force Data = over any number of seconds up to the length of data slice

EMG Data = (Frequency smoothing only) any integer number up to 400

cc: 31-33

Force signal gain with respect to force calibration (force cards only)

cc: 34-36

Muscle Type (cal data cards only)

1 = Brachial Bicep

2 = Brachioradialis


3 = Gastrocnemius

4 = Soleus

APPENDIX C


MAIN PROGRAM FLOW CHART






EMGAN Flowchart

Figure 1A



EMGAN Flowchart


Figure 1B



EMGAN Flowchart

Figure 1C





EMGAN Flowchart

Figure 1D

2

3

APPENDIX D

TABLE OF FORMULAS

## TABLE OF FORMULAS

### 1. Force Calibration Variable (counts/pound)

$$\text{FORSLP} = (H/20 - L/20)/\text{FLBCAL}$$

where  $H$  = sum of samples for one second of force calibration level

$L$  = sum of samples for one second of zero force level

$\text{FLBCAL}$  = amount of force calibration in pounds

### 2. EMG Calibration Variable (counts/microvolt)

$$\text{EMGAL} = \sqrt{\sum X_i^2 / (NPTS-1) / \text{MAGTUD}}$$

where  $\sum X_i^2$  = Sum of data values squared

$NPTS$  = number of EMG cal signal samples

$\text{MAGTUD}$  = RMS voltage of EMG cal signal

### 3. Mean Square (Variance)

$$\text{EMGVAR} = \sum X_i^2 / (NPTS-1)$$

where  $\sum X_i^2$  = sum of data values squared

$NPTS$  = number of data points

### 4. Simpson's Rule (without error term)

$$\text{VOLTSEC} = H/3 (2\text{SUMO} + 4 \text{SUME} - \text{Data (1)} + \text{Data (n)})$$

where  $H$  = spacing between the data values



SUMO = sum of the odd numbered data values excluding first data value

SUME = sum of the even numbered data values excluding last data value

Data (1) = first data value

Data (n) = last data value

#### 5. Scaling of EMG Data (microvolts)

$$\text{SEDATA}(i) = (\text{USEDATA}(i) - \text{OFFSET}) / \text{EMGCAL}$$

where  $\text{USEDATA}(i) = i^{\text{th}}$  unscaled data value

OFFSET = DC value of EMG data array

EMGVAL = EMG calibration value

#### 6. Scaling of Force Data (pounds)

$$\text{SFDATA}(i) = (\text{USFDA}(i) - \text{BASE}) * \text{GAIN} / \text{FORSLP}$$

where  $\text{USFDA}(i) = i^{\text{th}}$  unscaled force data value

BASE = zero force level

GAIN = ratio of force calibration gain to force data gain

FORSLP = force calibration variable

#### 7. Discrete Fourier Transform

$$X(n) = \sum_{K=0}^{N-1} x(k) (\text{EXP} (-j2 \pi / N))^{nk} \quad n = 0, 1, \dots, N-1$$

where  $x(k) = k^{\text{th}}$  value of untransformed data

n = number of data points

## 8. Power Spectral Density (magnitude)

$$\text{PSD}(n) = 8H/\text{NPTS} [\text{Re}(X(n))^2 + \text{Im}(X(n))^2]$$

where  $H$  = spacing between data samples

$\text{NPTS}$  = number of data samples

$\text{Re}(X(n))$  = real part of  $n^{\text{th}}$  Fourier transform value

$\text{Im}(X(n))$  = imaginary part of  $n^{\text{th}}$  Fourier transform value

## 9. Expected (Mean) Value of PSD

$$\text{EXPVAL} = \sum F (\text{PSD}(n) - \text{PERCNT}(n)/100)$$

where  $\text{PSD}(n)$  = value of  $n^{\text{th}}$  power spectral density bandwidth

$\text{PERCNT}(n)$  = percent of total power contributed by  $n^{\text{th}}$  bandwidth

$F$  = center frequency of  $n^{\text{th}}$  bandwidth

## 10. Standard Deviation of PSD

$$\text{STDEV} = \sqrt{\sum (F - \text{EXPVAL})^2 * (\text{PERCNT}(n)/100)}$$

where  $\text{EXPVAL}$  = expected value of PSD

$\text{PERCNT}(n)$  = percent of total power contributed by  $n^{\text{th}}$  bandwidth

$F$  = center frequency of  $n^{\text{th}}$  bandwidth

APPENDIX E

EXAMPLE OF EMGAN INPUT CARD DECK WITH PLOTTED  
AND PRINTED OUTPUT



E-1

00 17 41 00 41 30 0 1 1

888-57

NO. OF DATA SLICES= 11

DATA SLICES

00	1	17	41	0	-1	30	0	1	0	1	0	0	0	0
01	2	17	14	22	-1	4	0	0	0	0	3	0	0	0
01	4	17	16	16	-1	20	0	0	1	4	0	0	0	0
01	3	17	16	22	-1	4	0	0	10	0	0	0	0	0
01	4	17	16	55	-1	20	0	0	1	2	0	0	0	0
01	3	17	17	1	-1	4	0	0	10	0	0	0	0	0
01	4	17	17	35	-1	70	0	0	1	1	0	0	0	0
01	3	17	17	41	-1	4	0	0	10	0	0	0	0	0
01	3	17	17	57	-1	4	0	0	10	0	0	0	0	0
01	3	17	18	14	-1	4	0	0	10	0	0	0	0	0
01	3	17	18	31	-1	4	0	0	10	0	0	0	0	0

ORIGINAL PAGE IS  
OF POOR QUALITY

CARDIOVASCULAR LABORATORY  
EMG DATA PROCESSING PROGRAM  
-----

\*\*\* HEADER INFORMATION\*\*\*

SUBJECT NO.:	44	EXPERIMENT DATE:	6/28/75	FLIGHT REFERENCE DAY:	F-15
RUN NO.:	3	ANALOG TAPE NO.:	9	DIGITIZING DATE:	6/29/75
MUSCLE:	GASTROCNEMIUS				
EMG SAMPLE RATE (SAMP/SEC):	1000	FORCE SAMPLE RATE (SAMP/SEC):	20		

\*\* CALIBRATION DATA \*\*

FORCE CAL  
-----

CAL WEIGHT (LBS):	98.0	COUNTS/LB:	39.	AVG. BASELINE COUNT:	-3662.	CAL GAIN:	1.
-------------------	------	------------	-----	----------------------	--------	-----------	----

EMG CAL  
-----

RMS AMPLITUDE (MICROVOLTS):	350.	COUNTS/MICROVOLT:	6.44
-----------------------------	------	-------------------	------



\*\*\* FORCE DATA \*\*\*

TIME- 17:15:16

DATA LENGTH (SECS)- 20

AVERAGE FORCE INTERVAL (SECS)- 1.

INTERVAL	AVERAGE FORCE (LBS)
-----	-----

1	0.0
2	0.0
3	0.0
4	0.1
5	27.7
6	41.1
7	41.9
8	41.5
9	41.4
10	40.9
11	42.5
12	42.2
13	41.3
14	41.6
15	42.3
16	41.6
17	21.7
18	0.0
19	-0.0
20	0.0

\*\*\* POWER SPECTRAL DENSITY OF EMG DATA \*\*\*

START TIME- 17:16:22

DATA LENGTH (SECS)- 4

INTEGRATED EMG (MICROVOLT\*SEC)- 0.2676E 3

BANDWIDTH (HZ)- 9.766

NORMALIZED STANDARD ERROR- 0.158


INTEGRATED PSD (MICROVOLT\*\*2)- 0.7890E 4

MEAN (HZ)- 149.2

STANDARD DEVIATION(HZ)- 80.8

EMG VARIANCE (MICROVOLT\*\*2)- 0.8014E 4

FREQ (HZ)	PSD (MMV**2/HZ)	PSD NORM	% OF TOTAL	CUM % TOTAL
4.88	1.398	0.0306	0.17	0.17
14.65	6.766	0.1482	0.84	1.01
24.41	15.152	0.3319	1.88	2.89
34.18	23.153	0.5072	2.87	5.75
43.95	18.558	0.4065	2.30	8.05
53.71	30.225	0.6621	3.74	11.79
63.48	34.020	0.7452	4.21	16.00
73.24	27.169	0.5951	3.36	19.36
83.01	45.652	1.0000	5.65	25.01
92.77	32.585	0.7138	4.03	29.05
102.54	39.469	0.8646	4.89	33.93
112.30	42.109	0.9224	5.21	39.14
122.07	40.914	0.8962	5.06	44.21
131.84	45.400	0.9945	5.62	49.83
141.60	39.421	0.8635	4.88	54.71
151.37	43.597	0.9550	5.40	60.10
161.13	36.991	0.8103	4.58	64.68
170.90	34.035	0.7455	4.21	68.89
180.66	20.588	0.4510	2.55	71.44
190.43	21.723	0.4758	2.69	74.13
200.20	23.014	0.5041	2.85	76.98
209.96	27.384	0.5998	3.39	80.37
219.73	20.216	0.4428	2.50	82.87
229.49	15.688	0.3436	1.94	84.81
239.26	17.412	0.3814	2.16	86.97
249.02	13.422	0.2940	1.66	88.63
258.79	9.829	0.2153	1.22	89.84
268.55	10.397	0.2277	1.29	91.13
278.32	9.821	0.2151	1.22	92.35
288.09	6.011	0.1317	0.74	93.09
297.85	8.662	0.1897	1.07	94.16
307.62	7.959	0.1743	0.99	95.15
317.38	6.416	0.1405	0.79	95.94
327.15	5.909	0.1294	0.73	96.67
336.91	7.114	0.1558	0.88	97.55
346.68	5.590	0.1224	0.69	98.25
356.45	3.670	0.0803	0.45	98.70
366.21	4.024	0.0881	0.50	99.20
375.98	2.574	0.0564	0.32	99.52
385.74	2.234	0.0489	0.28	99.79
395.51	1.670	0.0366	0.21	100.00



CUBS 44

EXP DATE 6/23/75

RUN# 3

AMOUNT OF FORCE CAL (LBS) 99



EMC CAL

17:14:22 17:14:26

1000

500

MICRO  
VOLTS

-1000

-1000

0 500 1000

1500

2000

2500

3000

3500

4000

TIME IN HILISECONDS


EXP DATE 6/20/75

5000 44

RUN 3

FORCE DATA

17:16:16 17:16:36




SUBJ 44

EXP DATE 6/23/75

RUN# 3

EMG DATA

17:16:22 17:16:26



SU34 44


EXP DATE 6/29/75

RUN 3



## POWER SPECTRAL DENSITY

17:16:22 17:16:26



SUB 44

EXP DATE 6/28/75

RUN 3

BANDWIDTH (HZ) = 9.75

NORMALIZED STANDARD ERROR = .150

APPENDIX F

SIGMA SOURCE LISTING AND LOAD MAP

JJOB  
JASS UO=6  
JEDITOR  
-LIST  
-GEN

C\*\*\*\*\*

C

C

C PROGRAM: EMGAN

C

C AUTHOR: WILLIAM N. HURSTA

C

C PURPOSE: TO PROVIDE FIRST ORDER REDUCTION OF EMG AND FORCE DATA

C

C\*\*\*\*\*

C

C

DIMENSION IVARSW(10)  
INTEGER FILNOW,DTYPE,FILNUM,FILTSW,PLOTSW  
INTEGER FCALSW,ECALSW,WHOA,PRNTSW  
COMMON DATA(8500),IHEAD(20),CAL(2,4)  
EQUIVALENCE (DTYPE,IVARSW(1)),(FILNUM,IVARSW(2))  
EQUIVALENCE (FILTSW,IVARSW(3)),(NSECS,IVARSW(4))  
EQUIVALENCE (PLOTSW,IVARSW(5)),(PRNTSW,IVARSW(6))  
EQUIVALENCE (ISPAN,IVARSW(7)),(IFGAIN,IVARSW(8))  
FILNOW=1  
FCALSW=0  
ECALSW=0  
WHOA=0  
NSLICE=0

C

C

C

GO GET NUMBER OF DATA SLICES FOR PRESENT RUN

C

CALL SEGLD (1)

CALL CINPUT (IVARSW,NSLICE,I,START)

C

C

MAIN PROGRAM DO LOOP



C  
DO 200 I=1,NSLICE

C  
C  
GO GET DATA CARD INFO FOR PRESENT DATA SLICE  
C

CALL SEGLD (1)  
CALL CINPUT (IVARSW,NSLICE,I,START)  
SPAN=FLOAT(ISPAN)  
IF(DTYPE.EQ.4) FDGAIN=FLOAT(IFGAIN)

C  
C  
CHECK TO SEE IF CALIBRATION DATA ACQUIRED BEFORE PROCESSING FIRST DATA  
SLICE  
C

IF(DTYPE.EQ.1) FCALSW=1  
IF(DTYPE.EQ.2) ECALSW=1  
IF(DTYPE.EQ.3.AND.ECALSW.NE.1) WHOA=1  
IF(DTYPE.EQ.4.AND.FCALSW.NE.1) WHOA=1  
IF(WHOA.EQ.1) GO TO 90

C  
C  
GO GET THE DATA FROM TAPE  
C

CALL SEGLD(2)  
CALL TINPUT(DTYPE,FILNOW,FILNUM,START,NSECS)

C  
C  
DECIDE WHERE TO GO FOR EACH DATA TYPE.  
C

GO TO (20,10,30,70),DTYPE

C  
C  
FILTER 60 HZ FOR EMG CAL?  
C

10 IF(FILTSW.EQ.-1) GO TO 20  
CALL OUT60 (DTYPE,NSECS,FILTSW)

C  
C  
CALCULATE CALIBRATION VARIABLES OF EMG OR FORCE DATA.  
C

20 CALL SEGLD(3)  
CALL DCAL (DTYPE,NSECS,IFGAIN)  
C

```

C      PLOT EMG OR FORCE CAL DATA?
C
25  IF(PLOTSW.EQ.1) GO TO 26
    CALL SEGLD(4)
    CALL GRAPH (DTYPE,START,NSECS,SPAN,STDERR)
C
C      PRINT OUT HEADER?
C
26  IF(PRNTSW.EQ.1) GO TO 200
    GO TO 90
C
C      FILTER 60 HZ  FOR EMG DATA SLICE?
C
30  IF(FILTSW.EQ.-1) GO TO 40
    CALL OUT60 (DTYPE,NSECS,FILTSW)
C
C      SCALE EMG DATA, CALCULATE INTEGRATED VALUE AND MEAN SQUARE VALUE
C
40  CALL SEGLD (5)
    CALL EMG (NSECS,VOLTSEC)
C
C      PLOT EMG DATA ON COMPUTEK?
C
    IF(PLOTSW.EQ.1) GO TO 50
    CALL SEGLD(4)
    CALL GRAPH (DTYPE,START,NSECS,SPAN,STDERR)
C
C      TAKE THE FOURIER TRANSFORM OF THE DATA AND FIND PSD UP TO 400 HZ.
C
50  CALL SEGLD (6)
    CALL FFTPSD (NSECS,SPAN,PSDSUM,STDERR)
    DTYPE=5
C
C      PLOT THE PSD ON THE COMPUTEK?
C
    IF(PLOTSW.EQ.1) GO TO 60
    CALL SEGLD(4)
    CALL GRAPH (DTYPE,START,NSECS,SPAN,STDERR)

```

PRINT OUT RESULTS ON LINE PRINTER?

60 IF(PRNTSW.EQ.1) GO TO 200  
GO TO 90

SCALE FORCE DATA FOR PLOTTING AND AVERAGE IT FOR PRINTING.

70 CALL SEGLD (5)  
CALL FORCE (DTYPE,NSECS,SPAN,FDGAIN)

PLOT FORCE DATA ON COMPUTEK?

IF(PLOTSW.EQ.1) GO TO 80  
CALL SEGLD(4)  
CALL GRAPH (DTYPE,START,NSECS,SPAN,STDERR)

PRINT OUT FORCE DATA ON LINE PRINTER?

80 IF(PRNTSW.EQ.1) GO TO 200  
90 CALL SEGLD(7)  
CALL PRINT (DTYPE,START,NSECS,SPAN,VOLTSEC,PSDSUM,STDERR,WHOA)  
200 CONTINUE  
END  
SUBROUTINE CINPUT (IVARSW,NSLICE,I,START)

SUBROUTINE CINPUT ACQUIRES THE CARD DATA AND STORES IT ON A RAD FILE

ARGUMENTS: IVARSW- ARRAY HOLDING THE VARIABLES AND SWITCHES FOR EACH  
DATA SLICE (SEE EQUIVALENCE STATEMENTS IN MAIN  
PROGRAM)

NSLICE- NUMBER OF DATA SLICES IN PRESENT RUN

I - IDENTIFIES PRESENT DATA SLICE BEING PROCESSED

START - BEGIN TIME OF PRESENT DATA SLICE



```

DIMENSION IVARSW(10), ICARD(200,15)
DIMENSION IRDIO(10), IRDF(6)
COMMON DATA(8500), IHEAD(20), CAL(2,4)
EQUIVALENCE (ICARD, DATA(1))
DATA IRDIO(1) /'RD'/
DATA IRDF /Z4000, 'RD', 1,6000,0,0/
CALL DEFINE (IRDIO, IRDF)
IF(NSLICE.EQ.0) GO TO 30
GO TO 60

```

C  
C  
C

INPUT NUMBER OF DATA SLICES AND READ DATA SLICE INFO FOR PRESENT RUN

```

30 READ(105,20) NSLICE
WRITE(108,35) NSLICE
20 FORMAT(I5)
35 FORMAT(1H1, ' **** DATA CARDS ****', ///, ' NO. OF DATA SLICES=', I5, /
@///, 5X, 'DATA SLICES', ///)
IF(NSLICE.GT.200) GO TO 100
DO 50 J=1, NSLICE
READ(105,40) (ICARD(J,K), K=1,15)
40 FORMAT(15I3)
WRITE(108,40) (ICARD(J,K), K=1,15)
50 CONTINUE
CALL QOPEN(IRDIO, ICARD, 3000, 1, 0)
CALL QWRITE(IRDIO, 3000, INDIC)
CALL QCLOSE (IRDIO, 2)
GO TO 200

```

C  
C  
C

INITIALIZE VARIABLES AND SWITCHES FOR EACH DATA SLICE

```

60 CALL QOPEN(IRDIO, ICARD, 3000, 0, 0)
CALL QREAD(IRDIO, 3000, INDIC)
CALL QCLOSE (IRDIO, 2)
IVARSW(1)=ICARD(I,2)
IVARSW(2)=ICARD(I,1)+1
START=3600.*FLOAT(ICARD(I,3))+60.*FLOAT(ICARD(I,4))+FLOAT(ICARD(I,
@5))
IVARSW(3)=ICARD(I,6)

```

```

      IVARSW(4)=ICARD(I,7)
      IVARSW(5)=ICARD(I,8)
      IVARSW(6)=ICARD(I,9)
      IVARSW(7)=ICARD(I,10)
      IVARSW(8)=ICARD(I,11)
      IHEAD(16)=ICARD(I,3)
      IHEAD(17)=ICARD(I,4)
      IHEAD(18)=ICARD(I,5)
      IHEAD(19)=ICARD(I,7)
      IF(IVARSW(1)-2) 90,90,200
90  IHEAD(20)=ICARD(I,12)
      GO TO 200
100 WRITE(108,110) NSLICE
110 FORMAT(///,' ERROR...READ NUMBER OF DATA SLICES TO BE',I6,' CANN
      OT HAVE MORE THAN 200.')
```

STOP

```

200 RETURN
      END
      SUBROUTINE TINPUT(DTYPE,FILNOW,FILNUM,START,NSECS)
```

SUBROUTINE TINPUT READS FORCE AND EMG DATA IN FROM TAPE.

CALLING ARGUMENTS: DTYPE - DATA TYPE

FILNOW- FILE ON TAPE CURRENTLY BEING ACCESSED

FILNUM- FILE ON TAPE WHERE DESIRED DATA IS LOCATED

START - TIME IN TOTALED SECONDS OF BEGINNING OF  
DESIRED DATA SLICE

NSECS - NUMBER OF SECONDS OF DATA TO BE READ .

```

      INTEGER FILMOVE,FILNUM,DTYPE,FILNOW
      DIMENSION INTP(10),IARRAY(535)
      COMMON DATA(8500),IHEAD(20),CAL(2,4)
      EQUIVALENCE (IARRAY,DATA(8230))
```

```
INTP(1)=6  
CALL QOPEN (INTP,IARRAY,535,2,0)
```

```
POSITION TAPE AT THE BEGINNING OF THE DESIRED DATA FILE
```

```
IF(FILMOVE.LT.0) FILMOVE=FILMOVE-1  
CALL QFSKIP (INTP,FILMOVE,INDIC)  
GO TO (300,5),INDIC/4  
IF(FILMOVE.LT.0) CALL QFSKIP (INTP,1,INDIC)  
5 IF(FILMOVE.LT.0) FILMOVE=FILMOVE+1  
FILMOVE=FILNUM-FILNOW  
FILNOW=FILNOW+FILMOVE
```

```
POSITION TAPE AT THE BEGINNING OF THE DATA SLICE.
```

```
CALL FIND (FILNUM,START,INTP)
```

```
DECIDE WHERE TO GO FOR EACH DATA TYPE.
```

```
GO TO (10,50,50,10),DTYPE
```

```
ACQUIRE FORCE CAL OR FORCE DATA
```

```
10 JJ=1  
DO 30 J=1,2*NSECS  
CALL QREAD (INTP,535,INDIC)  
IF(INDIC.NE.0) GO TO 200  
DO 20 K=27,535,51  
DATA(JJ)=FLOAT(IARRAY(K)/8)  
JJ=JJ+1  
20 CONTINUE  
30 CONTINUE  
IF(DTYPE.EQ.1) GO TO 1000  
GO TO 150
```

```
ACQUIRE EMG CAL OR EMG DATA
```

```
50 JJ=1
```



```

DO 80 J=1,2*NSECS+1
CALL QREAD (INTP,535,INDIC)
IF(INDIC.NE.0) GO TO 200
DO 70 K=1,10
L=25+K+50*(K-1)
DO 60 KK=L,L+50
IF(KK.EQ.L+1) GO TO 60
DATA(JJ)=FLOAT(IARRAY(KK)/8)
JJ=JJ+1
60 CONTINUE
70 CONTINUE
80 CONTINUE
IF(DTYPE.EQ.2) GO TO 1000
150 IBACK=-(2*NSECS+2)
CALL QSKIP(INTP,IBACK,INDIC)
GO TO 1000

```

C  
C  
C

# ERROR MESSAGES

```

200 FILNUM=FILNUM-1
WRITE(108,210) INDIC,DTYPE,FILNUM
210 FORMAT(///,' ERROR...INDIC=',I4,' DURING INPUT OF DATA TYPE',I4,' L
eOCATED IN FILE',I4)
STOP
300 FILNUM=FILNUM-1
WRITE(108,310) FILNUM
310 FORMAT(///,' FATAL ERROR...EOT ENCOUNTERED DURING SEARCH FOR FILE'
e,I3)
STOP
1000 CONTINUE
CALL QCLOSE (INTP,0)
RETURN
SUBROUTINE FIND (FILNUM,START,INTP)

```

C  
C  
C  
C  
C

SUBROUTINE FIND POSITIONS THE TAPE AT THE BEGINNING OF A DATA SLICE.

CALLING ARGUMENTS: FILNUM- FILE ON TAPE WHERE DATA SLICE IS LOCATED

START - TIME IN TOTALED SECONDS OF THE BEGINNING OF  
DESIRED DATA SLICE.

INTP - CONTROL BLOCK FOR QINOUT TAPE ROUTINES

DIMENSION INTP(10),IARRAY(535)  
INTEGER FILNUM  
REAL NOW  
COMMON DATA(8500),IHEAD(20),CAL(2,4)  
EQUIVALENCE(IARRAY,DATA(8230))  
IPASS=0

READ EACH RECORD AND COMPARE ITS TIME LABEL WITH THE START TIME.

10 CALL QREAD (INTP,535,INDIC)  
IF(INDIC.NE.0) GO TO 20  
NOW=3600.\*ABS(FLOAT(IARRAY(2)))+60.\*ABS(FLOAT(IARRAY(3)))+ABS(FLOAT(IARRAY(4)))  
IF(NOW.NE.START) GO TO 10

ACQUIRE THE HEADER INFORMATION FOR THE PRESENT DATA SLICE.

DO 15 J=1,15  
IHEAD(J)=IARRAY(J+5)  
15 CONTINUE  
CALL QSKIP(INTP,-1,INDIC)  
GO TO 100

DID NOT FIND THE START TIME BEFORE ENCOUNTERING AN EOF. IF ONLY ONE PASS  
HAS BEEN MADE THRU THE DATA, BACK UP AND TRY AGAIN. OTHERWISE,STOP

20 IF(IPASS.EQ.1) GO TO 40  
CALL QFSKIP(INTP,-2,INDIC)  
IF(INDIC.EQ.8) GO TO 35  
CALL QFSKIP (INTP,1,INDIC)  
35 IPASS=1  
GO TO 10  
40 FILNUM=FILNUM-1

```

WRITE(108,50) FILNUM
50  FORMAT(///,'_CANNOT FIND START TIME IN FILE',I4)
STOP
100 CONTINUE
RETURN
END
SUBROUTINE OUT60 (DTYPE,NSECS,FILTSW)

```

SUBROUTINE OUT60 IS A DIGIT NOTCH FILTER. DEPENDING ON THE DATA CARD REQUEST, ONLY 60 HZ IS REMOVED OR ALL HARMONICS OF 60 HZ UP TO 360 HZ

FILTSW- SWITCH WHICH INDICATES WHETHER ONLY 60 HZ  
IS TO BE REMOVED OR ALL HARMONICS OF 60 HZ

```

      LOOPS=1
      TWOPI=6.28318530717
      IF(FILTSW.EQ.1) LOOPS=5
      NPTS=1024*NSECS
      DO 30 J=1,LOOPS,2
      A=0.0
      B=0.0
      W=(60*J)*TWOPI/FLOAT(IRATE)
      DO 10 K=1,NPTS
      A=A+DATA(K)*COS(K*W)
      B=B+DATA(K)*SIN(K*W)
10    CONTINUE
      A=A*2./FLOAT(NPTS)
      B=B*2./FLOAT(NPTS)
      DO 20 K=1,NPTS
      DATA(K)=DATA(K)-A*COS(K*W)-B*SIN(K*W)
20    CONTINUE

```



30 CONTINUE

100 RETURN

END

SUBROUTINE DCAL(DTYPE,NSECS,IFGAIN)

SUBROUTINE DCAL CALCULATES THE SCALE FACTORS FOR THE FORCE AND EMG DATA  
THE CALIBRATION CONSTANT FOR INTEGRATED EMG AREA IS ALSO FOUND. EMG  
THE CAL DATA IS THEN CONVERTED INTO UNITS OF POUNDS OR MICROVOLTS FOR  
PLOTING ON THE COMPUTEK TERMINAL.

CALLING ARGUMENTS: DTYPE - TYPE OF DATA

NSECS - LENGTH OF CAL DATA

DIMENSION HIGH(10),LOW(10),POUNDS(2600)

INTEGER DTYPE,FLBCAL

REAL LOWSUM

COMMON DATA(8500),IHEAD(20),CAL(2,4)

EQUIVALENCE (FLBCAL,IHEAD(15)),(IRATE,IHEAD(11))

EQUIVALENCE (MAGTUD,IHEAD(13)),(ICPS,IHEAD(14))

EQUIVALENCE (POUNDS,DATA(1))

IF(DTYPE.EQ.2) GO TO 100

FIND FORCE CALIBRATION VARIABLES

CAL(1,1)=FLOAT(FLBCAL)

CAL(1,4)=FLOAT(IFGAIN)

LOOK FOR JUMP IN DATA TO INDICATE FORCE CAL.

BEGIN=DATA(80)

DO 10 J=81,20\*NSECS-80

IF(DATA(J).GT.BEGIN+500.) GO TO 30

10 CONTINUE

WRITE (108,20)

20 FORMAT(///,'COULD NOT FIND THE FORCE CAL. ....CHECK DATA SLICE TIMES

e')

STOP

FOUND THE JUMP, NOW COMPUTE THE AVERAGE BEFORE AND AFTER THE JUMP  
CAL WEIGHT IN COUNTS = AVERAGE AFTER - AVERAGE BEFORE  
FORCE SLOPE = COUNTS/POUND  
ZERO FORCE = AVERAGE COUNTS JUST BEFORE CAL JUMP

30 HISUM=0.0  
LOWSUM=0.0  
DO 40 K=1,20  
HISUM=HISUM+DATA(J+K+60)  
LOWSUM=LOWSUM+DATA(J-K-60)  
40 CONTINUE  
FORSLP= (HISUM/20.-LOWSUM/20.)/FLOAT(FLBCAL)  
CAL(1,2)=FORSLP  
ZEROFOR=LOWSUM/20.  
CAL(1,3)=ZEROFOR

SCALE FORCE CALIBRATION DATA

DO 50 J=1,20\*NSECS  
POUNDS(J)=(DATA(J)-ZEROFOR)/FORSLP  
50 CONTINUE  
GO TO 200

FIND EMG CALIBRATION VARIABLES  
FIRST REMOVE ANY DC OFFSET IN THE CAL DATA

100 SUM=0.0  
INDEX=1000\*NSECS  
DO 110 J=1,INDEX  
SUM=SUM+DATA(J)  
110 CONTINUE  
OFFSET=SUM/FLOAT(INDEX)  
DO 120 J=1,INDEX  
DATA(J)=DATA(J)-OFFSET  
DATA(4000+J)=ABS(DATA(J)-OFFSET)  
120 CONTINUE

C  
C CALCULATE STANDARD DEVIATION AND OBTAIN SCALE FACTOR. FOR ZERO MEAN,  
C STATIONARY SIGNAL STANDARD DEVIATION=RMS VALUE  
C

XSQUAR=0.0  
DO 140 J=1,INDEX  
XSQUAR=XSQUAR+DATA(J)\*\*2  
140 CONTINUE  
EMGCAL=SQRT(XSQUAR/(INDEX-1))/FLOAT(MAGTUD)  
CAL(2,1)=FLOAT(MAGTUD)  
CAL(2,2)=EMGCAL

C  
C SCALE CAL DATA FOR PLOTTING.  
C

DO 180 J=1,INDEX  
DATA(J)=DATA(J)/EMGCAL  
180 CONTINUE  
200 RETURN  
END  
SUBROUTINE GRAPH (DTYPE,START,NSECS,SPAN,STDERR)

C  
C SUBROUTINE GRAPH PLOTS ON THE COMPUTEK TERMINAL THE DATA SLICE BEING  
C PROCESSED ALONG WITH VARIOUS HEADER INFORMATION.  
C

C CALLING ARGUMENTS: DTYPE - TYPE OF DATA  
C

C START - BEGIN TIME OF PRESENT DATA SLICE  
C

C NSECS - LENGTH OF DATA SLICE IN SECONDS  
C

C SPAN - AVERAGING INTERVAL FOR PSD  
C

C STDERR- NORMALIZED STANDARD ERROR FOR PSD ESTIMATE  
C

DIMENSION Z(15),DATE(3),POUNDS(2600),PSDN(410)  
DIMENSION FREQ(410),CALDATA(2000)  
INTEGER DTYPE  
COMMON DATA(8500),IHEAD(20),CAL(2,4)



```

EQUIVALENCE (POUNDS, DATA(1)), (PSDN, DATA(8000))
EQUIVALENCE (FREQ, DATA(5000)), (FLBCAL, CAL(1,1))
EQUIVALENCE (CALDATA, DATA(1)), (EMGMAX, DATA(8500))
EQUIVALENCE (FORMAX, DATA(8500))
SUB=FLOAT(IHEAD(1))
RUN=FLOAT(IHEAD(9))
DATE(1)=FLOAT(IHEAD(2))
DATE(2)=FLOAT(IHEAD(3))
DATE(3)=FLOAT(IHEAD(4))
FLBCAL=FLOAT(IHEAD(15))
SHOUR=FLOAT(IFIX(START/3600.))
SMIN=FLOAT(IFIX((START-SHOUR*3600.)/60.))
SSEC=START-SHOUR*3600.-SMIN*60.
FINISH=START+FLOAT(NSECS)
FHOUR=FLOAT(IFIX(FINISH/3600.))
FMIN=FLOAT(IFIX((FINISH-FHOUR*3600.)/60.))
FSEC=FINISH-FHOUR*3600.-FMIN*60.

```

```

      INITIALIZE PLOT AND OUTPUT HEADER INFORMATION.

```

```

CALL INITAL (Z)
CALL MODSET(Z,4,10.)
CALL MODSET(Z,5,10.)
CALL OBJECT(Z,75.,1008.,100.,770.)
CALL WORDS(Z,660.,780.,16., '      ' : ' ')
CALL WORDS(Z,100.,30.,4, 'SUB#')
CALL WORDS(Z,440.,30.,16, 'EXP DATE / / ')
CALL WORDS(Z,840.,30.,4, 'RUN#')
CALL NUMBER(Z,660.,780.,SHOUR,2,0)
CALL NUMBER(Z,696.,780.,SMIN,2,0)
CALL NUMBER(Z,732.,780.,SSEC,2,0)
CALL NUMBER(Z,780.,780.,FHOUR,2,0)
CALL NUMBER(Z,816.,780.,FMIN,2,0)
CALL NUMBER(Z,852.,780.,FSEC,2,0)
CALL NUMBER(Z,160.,30.,SUB,3,0)
CALL NUMBER(Z,548.,30.,DATE(1),2,0)
CALL NUMBER(Z,584.,30.,DATE(2),2,0)
CALL NUMBER(Z,620.,30.,DATE(3),2,0)

```

CALL NUMBER(Z,880.,30.,RUN,2,0)

DECIDE WHERE TO GO FOR EACH DATA TYPE

GO TO (10,40,40,10,70),DTYPE

LABEL AXES FOR FORCE DATA

10 CALL MODSET(Z,2,1.)

NPTS=20\*NSECS

PTS=FLOAT(NPTS)

CALL WORDS(Z,13.,380.,4,'LBS ')

CALL WORDS(Z,25.,355.,2,'OF')

CALL WORDS(Z,0.,330.,6,'FORCE ')

CALL WORDS(Z,450.,55.,16,'TIME IN SECONDS ')

IF(DTYPE.EQ.4) GO TO 30

FINISH UP PLOT FOR FORCE CALIBRATION DATA

CALL WORDS(Z,100.,780.,16,'FORCE CALIBRATION')

CALL WORDS(Z,350.,0.,26,'AMOUNT OF FORCE CAL (LBS) -')

CALL NUMBER (Z,674.,0.,FLBCAL,3,0)

CALL SUBJEC (Z,0.,FLOAT(NSECS),0.,200.)

CALL GRID(Z,2.,25.)

CALL LABEL(Z,1,2.,2,0)

CALL LABEL (Z,2,50.,3,0)

CALL SUBJEC (Z,0.,PTS,0.,200.)

CALL MODLINE(Z,NPTS,1,POUNDS)

CALL BELL(Z)

CALL PAGE(Z)

GO TO 90

FINISH UP PLOT FOR FORCE DATA

30 CALL WORDS(Z,100.,780.,10,'FORCE DATA')

FMAX=200.

IF(FORMAX.LT.80.) FMAX=100.

CALL SUBJEC (Z,0.,FLOAT(NSECS),0.,FMAX)

```
CALL GRID (Z,2,10.)
TICK=2.+FLOAT(NSECS/31)
CALL LABEL (Z,1,TICK,3,0)
CALL LABEL (Z,2,20,3,0)
CALL SUBJEC (Z,0,PTS,0,FMAX)
CALL MODLINE(Z,NPTS,1,POUNDS)
CALL BELL(Z)
CALL PAGE(Z)
GO TO 90
```

```
C
C      FINISH UP PLOT FOR EMG DATA
C
```

```
40 CALL WORDS(Z,430,55,20,'TIME IN MILLISECONDS')
CALL MODSET(Z,2,1.)
CALL WORDS (Z,0,440,6,'MICRO ')
CALL WORDS (Z,0,410,6,'VOLTS ')
NPTS=1000
SECMIL=FLOAT(1000*NSECS)
IF(DTYPE.EQ.3) GO TO 50
```

```
C
C      EMG CAL DATA
C
```

```
CALL WORDS(Z,100,780,8,'EMG CAL ')
CALL SUBJEC (Z,0,SECMIL,-1000,1000.)
CALL GRID (Z,250,250.)
CALL LABEL(Z,1,500,4,0)
CALL LABEL (Z,2,500,5,0)
CALL MODLINE(Z,NPTS,NSECS,CALDATA)
GO TO 60
```

```
C
C      EMG DATA
C
```

```
50 CALL WORDS(Z,100,780,8,'EMG DATA')
ULIM=(FLOAT(IFIX(EMGMAX/500.))+1.)*500.
BLIM=-ULIM
CALL SUBJEC (Z,0,SECMIL,BLIM,ULIM)
HACK=ULIM/4
CALL LABEL(Z,1,500,4,0)
```



```
CALL GRID (Z,250.,HACK)
WNUM=ULIM/2.
CALL LABEL (Z,2,WNUM,5,0)
CALL MODLINE(Z,NPTS,NSECS,DATA)
60 CALL BELL(Z)
CALL PAGE(Z)
GO TO 90
```

```
C
C
C      FINISH UP PLOT FOR PSD DATA
```

```
70 CALL WORDS(Z,100.,780.,22,'POWER SPECTRAL DENSITY')
CALL WORDS(Z,490.,55.,18,'FREQUENCY IN HERTZ')
CALL WORDS (Z,0.,425.,6,'MICRO ')
CALL WORDS(Z,0.,405.,6,'VOLTS ')
CALL WORDS (Z,60.,410.,2,'2 ')
CALL WORDS (Z,0.0,385.,6,'X1000 ')
CALL WORDS(Z,100.,0.,16,'BANDWIDTH (HZ) =')
CALL WORDS (Z,440.,0.,26,'NORMALIZED STANDARD ERROR=')
CALL NUMBER(Z,232.,0.,SPAN,2,2)
CALL NUMBER (Z,824.,0.,STDERR,1,3)
CALL SUBJEC (Z,0.,400.,0.,1,1)
CALL GRID (Z,25.,.05)
CALL LABEL (Z,1,25.,3,0)
CALL LABEL (Z,2.,10001,1,1)
```

```
C
C
C      CREATE X-AXIS FOR PSD VALUES
```

```
      NPTS=IFIX(400./SPAN)
      DO 80 J=1,NPTS
      FREQ(J)=SPAN*FLOAT(J)
80  CONTINUE
      CALL LINES(Z,NPTS,FREQ,PSDN)
      CALL BELL(Z)
      CALL PAGE(Z)
90  RETURN
      END
      SUBROUTINE MODLINE (Z,NPTS,MODX,Y)
```

SUBROUTINE MODLINE PLOTS A Y-ARRAY OF DATA POINTS WITHOUT REQUIRING A CORRESPONDING X-ARRAY. INSTEAD OF AN X-ARRAY, THE ARRAY POSITION OF A Y-ELEMENT IS USED FOR THE X VALUE.

CALLING ARGUMENTS: Z - THE PLOT ARRAY

SPECIAL CONSIDERATION: NPTS\*MODX SHOULD NOT BE GREATER THAN THE LENGTH  
OF THE Y-ARRAY

```

DIMENSION Z(15),Y(1)
CALL POINTS (Z,1,MODX,Y)
J2=IWRDX (Z(2),1)
DO 70 I=2,NPTS
  IY=I*MODX
  XX=FLOAT (I*MODX)
  GO TO (40,50),J2
40 YY=WORDXY(Y,IY)
  GO TO 60
50 YY=Y(IY)
60 CALL VCTOR(SCALE(Z,XX,0),SCALE(Z,YY,1),1)
70 CONTINUE
  RETURN
END
SUBROUTINE WORDS (Z,X,Y,NL,L)
  DIMENSION Z(15),L(50),M(50)
  NW=NL/2
  DO 10 I=1,NW
10 M(I)=L(I)
  CALL PASTOR (0,M,NW)
  CALL LEGEND (Z,X,Y,M,NL)
  RETURN

```

```

END
SUBROUTINE BELL (Z)
CALL MODE (Z,0)
CALL OUT (7,1)
CALL SCRIBE
DO 10 I=1,300 .
DO 10 J=1,60
10 K=I+J
CALL OUT (Z,1)
CALL SCRIBE
20 CALL IN (I)
IF(I.NE.13) GO TO 20
RETURN
END
SUBROUTINE EMG (NSECS,VOLTSEC)

```

```

C
C
C
C
C
C
C
C
C
SUBROUTINE EMG SCALES EMG DATA FOR PLOTTING ON THE COMPUTER AND FINDS
THE INTEGRATED EMG VALUE FOR THE DATA SLICE.

```

```

CALLING ARGUMENTS: NSECS - LENGTH OF DATA SLICE

```

```

VOLTSEC- INTEGRATED EMG VALUE IN MICROVOLT*SEC.

```

```

C
C
C
C
C
C
C
C
C
C
INTEGER DTYPE
COMMON DATA(8500),IHEAD(20),CAL(2,4)
EQUIVALENCE (IRATE,IHEAD(11)),(EMGCAL,CAL(2,2))
EQUIVALENCE (EMGMAX,DATA(8500)),(EMGVAR,DATA(8499))
SUM=0.0

```

```

C
C
C
CALCULATE OFFSET, SUBTRACT FROM DATA AND APPLY SCALE FACTOR

```

```

DO 10 J=1,1024*NSECS
SUM=SUM+DATA(J)
10 CONTINUE
OFFSET=SUM/(1024.*FLOAT(NSECS))
CAL(2,4)=OFFSET
EMGMAX=0.0
DO 20 J=1,1024*NSECS

```



C  
CALL RFORT (DATA,NEXP,S,-1,IFERR)

C  
C  
CALCULATE RAW POWER SPECTRAL DENSITY

C  
H=1./FLOAT(IRATE)

K=1

DO 10 J=3,N+1,2

PSD(K)=8.\*H/PTS\*(DATA(J)\*\*2+DATA(J+1)\*\*2)

K=K+1

10 CONTINUE

C  
C  
CORRECT PSD AMPLITUDES FOR WINDOW REDUCTION

DO 15 J=1,N/2

PSD(J)=1.1429\*PSD(J)

15 CONTINUE

C  
C  
OUTPUT RAW PSD TO TAPE

CALL PSDSAVE (NSECS)

C  
C  
CALCULATE SMOOTHED PSD

F0=1./((PTS/FLOAT(IRATE))

IOTA=IFIX(SPAN/F0)

SPAN=FLOAT(IOTA)\*F0

PSDMAX=0.0

K=1

DO 30 J=1,N/2,IOTA

TEMP=0.0

DO 20 JJ=1,IOTA

TEMP=TEMP+PSD(J+JJ-1)

20 CONTINUE

PSD(K)=TEMP/FLOAT(IOTA)

IF(FLOAT(J/IOTA)\*SPAN.GT.400.+SPAN) GO TO 25

IF(PSD(K).GT.PSDMAX) PSDMAX=PSD(K)

25 K=K+1

30 CONTINUE

C  
C  
C

CALCULATE THE NORMALIZED PSD

LIM=IFIX(400./SPAN)+1

DO 40 J=1,LIM

PSDN(J)=PSD(J)/PSDMAX

40 CONTINUE

C  
C  
C

INTEGRATE THE SMOOTHED PSD

PSDSUM=0.0

DO 50 J=1,LIM

PSDSUM=PSDSUM+SPAN\*PSD(J)

50 CONTINUE

C  
C  
C

CALCULATE % EACH BANDWIDTH IS OF TOTAL AND FIND THE CUMLATIVE % OF TOTAL

DO 60 J=1,LIM

PERCNT(J)=SPAN\*PSD(J)/PSDSUM\*100.

IF(J.NE.1) GO TO 55

CUMPCT(1)=PERCNT(1)

GO TO 60

55 CUMPCT(J)=CUMPCT(J-1)+PERCNT(J)

60 CONTINUE

C  
C  
C

CALCULATE THE EXPECTED VALUE AND VARIANCE OF THE PSD

F=SPAN/2.

EXPVAL=0.0

DO 70 J=1,LIM

EXPVAL=EXPVAL+(PERCNT(J)/100.)\*F

F=F+SPAN

70 CONTINUE

F=SPAN/2.

VARVAL=0.0

DO 80 J=1,LIM

VARVAL=VARVAL+(PERCNT(J)/100.)\*(F-EXPVAL)\*\*2

F=F+SPAN  
S0 CONTINUE  
STDEV=SQRT(UARVAL)

C  
C  
C  
FIND STANDARDIZED NORMAL ERROR

STDERR=SQRT(1./FLOAT(IOTA))  
RETURN  
END  
SUBROUTINE WINDOW (NSECS)

C  
C  
C  
C  
C  
SUBROUTINE WINDOW APPLIES A COSINE TAPER TO THE FIRST AND LAST TENTHS OF  
THE DATA TO REDUCE LEAKAGE.

CALLING ARGUMENTS: NSECS - LENGTH OF THE DATA

COMMON DATA(8500),IHEAD(20),CAL(2,4)

PI=3.1415927

NPTS=1024\*NSECS

IEDGE=NPTS/10

TTOTAL=1.024\*FLOAT(NSECS)

K=1

DO 10 J=1,IEDGE

T1=.001\*FLOAT(J)

C1=.5\*(1.-COS(PI\*T1/(.1\*TTOTAL)))

DATA(J)=C1\*DATA(J)

DATA(NPTS-J+1)=C1\*DATA(NPTS-J+1)

10 CONTINUE

RETURN

END

SUBROUTINE RFORT(A,M,S,IFS,IFERR)

ONE-DIMENSIONAL REAL FINITE FOURIER TRANSFORM

C  
C  
C  
C  
FOURIER TRANSFORM SUBROUTINE FOR REAL DATA.

PROGRAMMED IN SYSTEM/360, BASIC PROGRAMMING SUPPORT,

C  
C  
FORTAN IV, (SEE FORM C28-6504).

C  
THIS DECK IS SET UP FOR IBSYS ON THE IBM 7094



THIS PROGRAM USES THE SUBROUTINE FORT TO COMPUTE COMPLEX  
FOURIER TRANSFORMS OF REAL DATA. PK FORT S.D.A. NO. 3465 IS  
AVAILABLE THROUGH SHARE.

THE FOURIER SERIES IS

$$X(J) = \text{SUM OVER } K=0 \text{ TO } N, \text{ OF } C(K) * \exp(2 * \pi * I * J * K / N)$$
$$J=0,1,2,\dots,N-1$$

WHERE  $I = \sqrt{-1}$  AND WHERE  $C(K)$  IS COMPLEX.  
SINCE  $X(J)$  IS REAL,  $C(K) = \text{CONJG}(C(N-K))$ . THEREFORE ONLY  
 $C(K), K=0,1,\dots,N/2$  ARE COMPUTED AND/OR USED.

#### ARGUMENTS-

A IS INITIALLY THE INPUT ARRAY, X, WHEN COMPUTING A FOURIER  
TRANSFORM AND C WHEN COMPUTING A FOURIER SERIES. A IS REPLACES BY  
THE OUTPUT ARRAY, C IN THE FORMER CASE, X IN THE LATTER.  
THE X VECTOR CONTAINS THE REAL DATA  $X(0), X(1), \dots, X(N-1)$   
THE C VECTOR CONTAINS THE COMPLEX FOURIER AMPLITUDES  
 $C(0), C(1), \dots, C(N/2)$ . THE COMPLEX VECTOR C IS STORED ACCORDING  
TO THE NORMAL FORTRAN IV CONVENTION FOR STORING COMPLEX NUMBERS.  
I.E., REAL PARTS IN ALTERNATE CELLS STARTING WITH THE FIRST,  
IMAGINARY PARTS IN ALTERNATE CELLS STARTING WITH THE SECOND.  
TO ADHERE TO FORTRAN RULES,  $X(0), X(1), \dots$  ARE REFERRED TO AS  
 $X(1), X(2), \dots$ , RESP. IN THE PROGRAMS. ALSO,  $C(0), C(1), \dots$  ARE  
REFERRED TO AS  $C(1), C(2), \dots$ , RESP., IF C IS DESIGNATED AS  
COMPLEX IN A TYPE STATEMENT.

M GIVES  $N=2*M$

THE ARGUMENTS S, IFS, AND IFERR ARE THE SAME AS IN THE  
SUBROUTINE FORT AND THE USER IS REFERRED TO THE COMMENT CARDS  
IN FORT THEIR EXPLANATION.

DIMENSION STATEMENTS- THE DIMENSIONS OF ARRAYS A AND S SHOULD  
BE  $N+2$  AND  $N/4$ , RESP. FOR THE LARGEST N TO BE USED. FOR  
EXAMPLE, IF THE LARGEST M IS 13, THEN,  $N=8192$  AND ONE SHOULD

```

C      HAVE THE DIMENSION STATEMENT-
C      DIMENSION A(8194), S(2048)
C      IF ONE WISHES TO SPECIFY A TO BE COMPLEX BY A TYPE STATEMENT,
C      ONE SHOULD GIVE IT A DIMENSION OF N/2 +1 , FOR THE LARGEST N.
      DIMENSION A(8500),S(1024)
10  IFERRS = 0
      N=2**M
      NU2 = N / 2
      NU4M1 = N/4 - 1
      MM1 = M - 1
      IF (ABS(IFS)-1) 40,40,20
20  IF (MP-M) 30,50,50
30  IFERRS = 1
40  NP = N
      MP = M
      CALL FORT(A,M,S,0,IFERR1)
      IFERRS = IFERRS + IFERR1
50  KD = NP / N
      KT = KD
      NPU4 = NP / 4
      IF (IFS) 60,80,90
60  CALL FORT(A,MM1,S,-2,IFERR2)
      IFERRS = IFERRS + IFERR2
      DO 70 K=1,NU4M1
          J=NU2-K
          A1R= A(2*K+1) + A(2*J+1)
          A1I=A(2*K+2)-A(2*J+2)
          A2R=A(2*K+2)+A(2*J+2)
          A2I=A(2*J+1)-A(2*K+1)
          KKT = NPU4-KT
          AWR=A2R*S(KKT) + A2I*S(KT)
          AWI = A2I*S(KKT)- A2R*S(KT)
          A(2*K+1)=(A1R+AWR)/4.
          A(2*K+2)=(A1I+AWI)/4.
          A(2*J+1)=(A1R-AWR)/4.
          A(2*J+2)=(AWI-A1I)/4.
70  KT=KT+KD
      T=A(1)

```

```

      A(1)=(T+A(2))/2
      A(N+1) = (T-A(2))/2
      A(2)=0
      A(N+2) = 0
      A(NV2+1) = .5*A(NV2+1)
      A(NV2+2) = -.5 * A(NV2+2)
80  IFERR = IFERRS
      RETURN
90  DO 100 K=1,NV4M1
      J=NV2-K
      A1R=A(2*K+1) + A(2*J+1)
      A1I=A(2*K+2)-A(2*J+2)
      AWR=A(2*K+1)-A(2*J+1)
      AWI=A(2*K+2)+A(2*J+2)
      KKT = NPV4 - KT
      A2R=AWR*S(KKT) - AWI*S(KT)
      A2I=AWR*S(KT) + AWI*S(KKT)
      A(2*K+1) = A1R - A2I
      A(2*K+2) = A1I + A2R
      A(2*J+1) = A1R + A2I
      A(2*J+2) = A2R - A1I
100  KT = KT + KD
      T = A(1)
      A(1) = T + A(N+1)
      A(2) = T - A(N+1)
      A(NV2+1) = 2.*A(NV2+1)
      A(NV2+2) = -2.*A(NV2+2)
      CALL FORT(A,MM1,S,2,IFERR2)
      IFERRS = IFERRS+IFERR2
      GO TO 80
      END
      SUBROUTINE FORT(A,M,S,IFS,IFERR)

```

FORT, ONE-DIMENSIONAL FINITE COMPLEX FOURIER TRANSFORM.

FOURIER TRANSFORM SUBROUTINE, PROGRAMMED IN SYSTEM/360,  
 BASIC PROGRAMMING SUPPORT, FORTRAN IV. FORM C28-6504  
 THIS DECK SET UP FOR IBSYS ON IBM 7094.



DOES EITHER FOURIER SYNTHESIS, I.E., COMPUTES COMPLEX FOURIER SERIES  
GIVEN A VECTOR OF N COMPLEX FOURIER AMPLITUDES, OR, GIVEN A VECTOR  
OF COMPLEX DATA X DOES FOURIER ANALYSIS, COMPUTING AMPLITUDES  
A IS A COMPLEX VECTOR OF LENGTH  $N=2**M$  COMPLEX NOS. OR  $2*N$  REAL  
NUMBERS. A IS TO BE SET BY USER

M IS AN INTEGER 0.LT.M.LE.13, SET BY USER.

S IS A VECTOR  $S(J)=\sin(2*PI*J/NP)$ ,  $J=1,2,\dots,NP/4-1$ ,  
COMPUTED BY PROGRAM.

IFS IS A PARAMETER TO BE SET BY USER AS FOLLOWS-

IFS=0 TO SET  $NP=2**M$  AND SET UP SINE TABLE.

IFS=1 TO SET  $N=NP=2**M$ , SET UP SIN TABLE, AND DO FOURIER  
SYNTHESIS, REPLACING THE VECTOR A BY

$X(J)=\text{SUM OVER } K=0,N-1 \text{ OF } A(K)*\exp(2*PI*I/N)**(J*K)$ ,  
 $J=0,N-1$ , WHERE  $I=\sqrt{-1}$

THE X'S ARE STORED WITH RE  $X(J)$  IN CELL  $2*J+1$   
AND IM  $X(J)$  IN CELL  $2*J+2$  FOR  $J=0,1,2,\dots,N-1$ .  
THE A'S ARE STORED IN THE SAME MANNER.

IFS=-1 TO SET  $N=NP=2**M$ , SET UP SIN TABLE, AND DO FOURIER  
ANALYSIS, TAKING THE INPUT VECTOR A AS X AND  
REPLACING IT BY THE A SATISFYING THE ABOVE FOURIER SERIES.

IFS>0 SET UP SIN TABLE AND RETURN

IFS=+2 TO DO FOURIER SYNTHESIS ONLY, WITH A PRE-COMPUTED S.

IFS=-2 TO DO FOURIER ANALYSIS ONLY, WITH A PRE-COMPUTED S.

IFERR IS SET BY PROGRAM TO-

=0 IF NO ERROR DETECTED.

=1 IF M IS OUT OF RANGE, OR, WHEN IFS=+2,-2, THE  
PRE-COMPUTED S TABLE IS NOT LARGE ENOUGH.

=-1 WHEN IFS =+1,-1, MEANS ONE IS RECOMPUTING S TABLE  
UNNECESSARILY.

NOTE- AS STATED ABOVE, THE MAXIMUM VALUE OF M FOR THIS PROGRAM  
ON THE IBM 7094 IS 13. FOR 360 MACHINES HAVING GREATER STORAGE  
CAPACITY, ONE MAY INCREASE THIS LIMIT BY REPLACING 13 BY

C STATEMENT 3 BELOW BY LOG2 N, WHERE N IS THE MAX NO. OF  
C COMPLEX NUMBERS ONE CAN STORE IN HIGH-SPEED CORE. ONE MUST  
C ALSO ADD MORE DO STATEMENTS TO THE BINARY SORT ROUTINE  
C FOLLOWING STATEMENT 24 AND CHANGE THE EQUIVALENCE STATEMENTS  
C FOR THE K'S.

C DIMENSION A(8500),S(1024),K(15)  
C IF (M) 20,20,10  
10 IF (M-15) 40,40,20  
20 IFERR=1  
30 RETURN  
40 IFERR=0  
N=2\*\*M  
IF (IABS(IFS)-1) 440,440,50  
C WE ARE DOING TRANSFORM ONLY. SEE IF PRE-COMPUTED  
C S TABLE IS SUFFICIENTLY LARGE  
50 IF (N-NP) 70,70,60  
60 IFERR=1  
GO TO 440  
C SCRAMBLE A, BY SANDE'S METHOD  
70 K(1)=2\*N  
DO 80 L=2,M  
80 K(L)=K(L-1)/2  
DO 90 L=M,14  
90 K(L+1)=2  
C THE FOLLOWING 15 STATEMENTS ARE TO COMPENSATE FOR A WEAKNESS IN  
C THE FORTRAN V COMPILER  
K1 =K(15)  
K2 =K(14)  
K3 =K(13)  
K4 =K(12)  
K5 =K(11)  
K6 =K(10)  
K7 =K(9)  
K8 =K(8)  
K9 =K(7)  
K10=K(6)  
K11=K(5)

```

      K12=K(4)
      K13=K(3)
      K14=K(2)
      K15=K(1)
      N2 =K(1)
C     NOTE EQUIVALENCE OF KL AND K(14-L)
C     BINARY SORT-
      IJ =2
100   J1 =2
110   J2 =J1
120   J3 =J2
130   J4 =J3
140   J5 =J4
150   J6 =J5
160   J7 =J6
170   J8 =J7
180   J9 =J8
190   J10=J9
200   J11=J10
210   J12=J11
220   J13=J12
230   J14=J13
240   JI=J14
250   IF (IJ-JI) 260,270,270
260   T=A(IJ-1 )
      A(IJ-1)=A(JI-1)
      A(JI-1)=T
      T=A(IJ)
      A(IJ)=A(JI)
      A(JI)=T
270   IJ=IJ+2
      JI=JI+K14
      IF (JI.LE.K15) GO TO 250
      J14=J14+K13
      IF (J14.LE.K14) GO TO 240
      J13=J13+K12
      IF (J13.LE.K13) GO TO 230
      J12=J12+K11

```



```

IF (J12.LE.K12) GO TO 220
J11=J11+K10
IF (J11.LE.K11) GO TO 210
J10=J10+K9
IF (J10.LE.K10) GO TO 200
J9 =J9+K8
IF (J9.LE.K9) GO TO 190
J8=J8+K7
IF (J8.LE.K8) GO TO 180
J7=J7+K6
IF (J7.LE.K7) GO TO 170
J6=J6+K5
IF (J6.LE.K6) GO TO 160
J5=J5+K4
IF (J5.LE.K5) GO TO 150
J4=J4+K3
IF (J4.LE.K4) GO TO 140
J3=J3+K2
IF (J3.LE.K3) GO TO 130
J2=J2+K1
IF (J2.LE.K2) GO TO 120
J1=J1+2
IF (J1.LE.K1) GO TO 110
IF (IFS) 280,20,300

```

C DOING FOURIER ANALYSIS,SO DIV. BY N AND CONJUGATE.

280 FN = N

DO 290 I=1,N

A(2\*I-1)=A(2\*I-1)

290 A(2\*I)=-A(2\*I)

C SPECIAL CASE- L=1

300 DO 310 I=1,N,2

T = A(2\*I-1)

A(2\*I-1) = T + A(2\*I+1)

A(2\*I+1)=T-A(2\*I+1)

T=A(2\*I)

A(2\*I) = T + A(2\*I+2)

310 A(2\*I+2)= T - A(2\*I+2)

IF (M-1) 20,30,320

```

C      SET FOR L=2
320    LEXP1=2
C      LEXP1=2***(L-1)
      LEXP=8
C      LEXP=2***(L+1)
      NPL= 2***MT
C      NPL = NP* 2**-L
330    DO 390 L=2,M
C      SPECIAL CASE- J=0
      DO 340 I=2,N2,LEXP
      I1=I + LEXP1
      I2=I1+ LEXP1
      I3 =I2+LEXP1
      T=A(I-1)
      A(I-1) = T +A(I2-1)
      A(I2-1) = T-A(I2-1)
      T =A(I)
      A(I) = T+A(I2)
      A(I2) = T-A(I2)
      T= -A(I3)
      TI = A(I3-1)
      A(I3-1) = A(I1-1) - T
      A(I3 ) = A(I1 ) - TI
      A(I1-1) = A(I1-1) +T
340    A(I1) = A(I1 ) +TI
      IF (L-2) 380,380,350
350    KLAST=N2-LEXP
      JJ=NPL
      DO 370 J=4,LEXP1,2
      NPJJ=NT-JJ
      UR=S(NPJJ)
      UI=S(JJ)
      ILAST=J+KLAST
      DO 360 I=J,ILAST,LEXP
      I1=I+LEXP1
      I2=I1+LEXP1
      I3=I2+LEXP1
      T=A(I2-1)*UR-A(I2)*UI

```

```

      IARRAY(J)=0
10  CONTINUE
      DO 20 J=1,20
      IARRAY(J)=IHEAD(J)
20  CONTINUE
      NRECS=NSECS/2
      IF(NRECS.EQ.0) NRECS=1
      IARRAY(21)=NRECS
      CALL QWRITE(IOUTP,25,INDIC)

```

```

      OUTPUT FOLLOWING DATA RECORDS CONTAINING THE RAW PSD VALUES

```

```

      K=0
      DO 50 J=1,NRECS
      DO 40 JJ=1,1024
      DATA(4099+JJ)=DATA(K+JJ)
40  CONTINUE
      CALL QWRITE (IOUTP,2048,INDIC)
      K=K+1024
50  CONTINUE

```

```

      WRITE TWO EOFs AND BACK UP TO JUST BEFORE THEM

```

```

      CALL QWEOF (IOUTP,1,INDIC)
      CALL QCLOSE (IOUTP,0)
      RETURN
      END

```

```

      SUBROUTINE PRINT (DTYPE,START,NSECS,SPAN,VOLTSEC,PSDSUM,STDERR,WHO
      @A)

```

```

      SUBROUTINE PRINT OUTPUTS TO THE LINE PRINTER HEADER AND CALIBRATION DATA,
      AVERAGED FORCE DATA, AND SMOOTHED PSD SPECTRUM.

```

```

      CALLING ARGUMENTS:  DTYPE - TYPE OF DATA

```

```

                        START - BEGIN TIME OF PRESENT DATA SLICE

```

```

                        NSECS - LENGTH OF DATA SLICE

```



ORIGINAL PAGE IS  
OF POOR QUALITY

SPAN - FOR FORCE DATA NUMBER OF SECS FORCE DATA IS  
AVERAGED OVER; FOR PSD DATA FREQUENCY  
INTERVAL FOR SMOOTHING ESTIMATES

VOLTSEC- INTEGRATED EMG VALUE (TIME DOMAIN)

PSDSUM- INTEGRATED EMG VALUE (FREQUENCY DOMAIN)

STDERR- NORMALIZED STANDARD ERROR OF ESTIMATE FOR  
PSD VALUES

WHOA - FLAG SET WHEN DATA PROCESSING ATTEMPTED  
WITHOUT CALIBRATION

INTEGER DTYPE,WHOA  
DIMENSION IEDATE(3),IDDATE(3),MUSCLE(7,4),IFORR(2)  
DIMENSION AUGLBS(130),PSD(410),ITIME(3),PSDN(410),PERCNT(410)  
DIMENSION CUMPCT(410)  
COMMON DATA(8500),IHEAD(20),CAL(2,4)  
EQUIVALENCE (ISUBNO,IHEAD(1)),(IEDATE,IHEAD(2))  
EQUIVALENCE (IDDATE,IHEAD(5)),(IFLITE,IHEAD(8))  
EQUIVALENCE (IRUN,IHEAD(9)),(MUS,IHEAD(20))  
EQUIVALENCE (IATAPE,IHEAD(10)),(ISAMPE,IHEAD(11))  
EQUIVALENCE (ISAMPH,IHEAD(12)),(AUGLBS,DATA(5000))  
EQUIVALENCE (PSD,DATA(1)),(PSDN,DATA(8000))  
EQUIVALENCE (PERCNT,DATA(7000)),(CUMPCT,DATA(7500))  
EQUIVALENCE (EXPVAL,DATA(6500)),(STDEV,DATA(6501))  
EQUIVALENCE (EMGUAR,DATA(8499))  
DATA MUSCLE // 'BR','AC','HI','AL','B','IC','EP','B',' ','RA','DI  
e','AL','IS',' ','G','AS','TR','OC','NE','HI','US','S','OL','EU'  
e','S',' ',' ',' ',' '  
DATA IFORR // 'F','R+'  
IF(WHOA.EQ.1) GO TO 150  
ITIME(1)=IFIX(START/3600.)  
ITIME(2)=IFIX((START-FLOAT(ITIME(1))\*3600.)/60.)  
ITIME(3)=IFIX(START-FLOAT(ITIME(1))\*3600.-FLOAT(ITIME(2))\*60.)

BRANCH TO PROPER PRINT OUT SECTION FOR THE DATA TYPE

IF(DTYPE-4) 5,50,100

PRINT OUT HEADER

5 IP=1

IF (IFLITE.GE.0) IP=2

WRITE(108,10)

10 FORMAT(1H1,/,48X,'CARDIOVASCULAR LABORATORY',/,47X,'EMG DATA PRO  
@CESSING PROGRAM',/,47X,'-----')

WRITE(108,20)ISUBNO,IEDATE,IFORR(IP),IFLITE,IRUN,IATAPE,IDDATE,(MU  
@SCLE(J,MUS),J=1,7),ISAMPE,ISAMPM

20 FORMAT(////,47X,'\*\*\* HEADER INFORMATION\*\*\*',////,' SUBJECT NO.: ', 7  
@X,I3,18X,'EXPERIMENT DATE: ',6X,I2,'/',I2,'/',I2,10X,'FLIGHT REFERE  
@NCE DAY: ',A3,I3,/,', RUN NO.: ',11X,I2,19X,'ANALOG TAPE NO.: ', 7X,I  
@3,14X,'DIGITIZING DATE: ',6X,I2,'/',I2,'/',I2,/,', MUSCLE: ', 9X,7A2  
@,/,', EMG SAMPLE RATE (SAMP/SEC): ',I5,8X,'FORCE SAMPLE RATE (SAMP/  
@SEC): ',I5)

WRITE(108,30)(CAL(1,J),J=1,4)

30 FORMAT(////,49X,'\*\* CALIBRATION DATA \*\*',////,' FORCE CAL',/,', ----  
@-----',/,', CAL WEIGHT (LBS): ',F7.1,16X,'COUNTS/LB: ',2X,F6.0,5X,'A  
@UG. BASELINE COUNT: ',2X,F6.0,5X,'CAL GAIN: ',2X,F3.0)

WRITE(108,40) (CAL(2,J),J=1,2)

40 FORMAT(////,' EMG CAL',/,', -----',/,', RMS AMPLITUDE (MICROVOLTS)  
@: ',F6.0,7X,' COUNTS/MICROVOLT: ',F5.2,/,')

GO TO 200

PRINT OUT FORCE DATA

50 WRITE(108,60) ITIME,NSECS,SPAN

60 FORMAT(1H1,51X,'\*\*\* FORCE DATA \*\*\*'///5X,'START TIME-',I3,' ',I2,'  
@-',I2,15X,'DATA LENGTH (SECS)-',I4,15X,'AVERAGE FORCE INTERVAL (SE  
@CS)-',F3.0,/,')

LOOPS=NSECS/IFIX(SPAN)

IF(LOOPS.GT.49) GO TO 66

WRITE(108,62)

```

62 FORMAT(47X,'INTERVAL',5X,'AVERAGE FORCE (LBS)'/,47X,'-----',5X,
@'-----',//)
DO 65 J=1,LOOPS
WRITE(108,64) J,AUGLBS(J)
64 FORMAT(49X,I3,15X,F5.1)
65 CONTINUE
GO TO 200

```

```

PRINT OUT PSD DATA

```

```

66 IF((LOOPS/2)*2.NE.LOOPS) LOOPS=LOOPS-1
WRITE(108,68)
68 FORMAT(15X,'INTERVAL',5X,'AVERAGE FORCE (LBS)',20X,'INTERVAL',5X,
@AVERAGE FORCE (LBS)'/,15X,'-----',5X,'-----',20X
@,'-----',5X,'-----',//)
DO 75 J=1,LOOPS/2
JJ=J+LOOPS/2
WRITE(108,73) J,AUGLBS(J),JJ,AUGLBS(JJ)
73 FORMAT(17X,I3,15X,F5.1,29X,I3,15X,F5.1)
75 CONTINUE
GO TO 200
100 WRITE(108,110) ITIME,NSECS,VOLTSEC,SPAN,STDERR,PSDSUM,EXPVAL,STDEV
@,EMGUAR
110 FORMAT(1H1,39X,'*** POWER SPECTRAL DENSITY OF EMG DATA ***'///4X,
@START TIME-',I3,',',I2,',',I2,11X,'DATA LENGTH (SECS)-',I4,14X,'IN
@TEGRATED EMG (MICROVOLT*SEC)-',E12.4//4X,'BANDWIDTH (HZ)-',F6.3,10
@X,'NORMALIZED STANDARD ERROR-',F6.3,5X,'INTEGRATED PSD (MICROVOLT*
@*2)-',E13.4,/,4X,'MEAN (HZ)-',F6.1,15X,'STANDARD DEVIATION
@ (HZ)-',F6.1,8X,'EMG VARIANCE (MICROVOLT**2)-',E15.4,/)
LOOPS=IFIX(400./SPAN)+1
FREQ=SPAN/2
IF(LOOPS.GT.46) GO TO 130
WRITE(108,115)
115 FORMAT(30X,' FREQ',13X,'PSD',11X,'PSD',10X,'% OF',10X,'CUM %'//,30X
@,' (HZ)',8X,'(MMU**2/HZ)',8X,'NORM',9X,'TOTAL',9X,'TOTAL'//,31X,
@'-----',8X,'-----',8X,'-----',9X,'-----',9X,'-----',/)
DO 125 J=1,LOOPS
WRITE(108,120) FREQ,PSD(J),PSDN(J),PERCNT(J),CUMPCT(J)

```



```

120 FORMAT(30X,F6.2,7X,G10.5,8X,F6.4,7X,F6.2,8X,F6.2)
    FREQ=FREQ+SPAN
125 CONTINUE
    GO TO 200
130 IF((LOOPS/2)*2.NE.LOOPS) LOOPS=LOOPS-1
    FREQ2=FLOAT(LOOPS/2)*SPAN+FREQ
    WRITE(108,135)
135 FORMAT(' FREQ', 8X,'PSD', 9X,'PSD',8X,'% OF',6X,'CUM %',22X,'FREQ
@', 9X,'PSD', 8X,'PSD',8X,'% OF',6X,'CUM %',2X,'(HZ)',4X,'(MMU**2/
@HZ)',5X,'NORM',7X,'TOTAL',5X,'TOTAL',22X,'(HZ)',5X,'(MMU**2/HZ)',5
@X,'NORM',6X,'TOTAL',6X,'TOTAL',/, ' ----',4X,'-----',5X,'---
@-',7X,'-----',5X,'-----',22X,'----',5X,'-----',5X,'----',6X,
A'-----',6X,'-----',/)
    DO 145 J=1,LOOPS/2
    JJ=J+LOOPS/2
    WRITE(108,140) FREQ,PSD(J),PSDN(J),PERCNT(J),CUMPCT(J),FREQ2,PSD(J
@J),PSDN(JJ),PERCNT(JJ),CUMPCT(JJ)
140 FORMAT(F7.2,3X,G10.5,5X,F6.4,5X,F6.2,5X,F6.2,20X,F6.2,4X,G10.5,5X,
@F6.4,5X,F6.2,5X,F6.2)
    FREQ=FREQ+SPAN
    FREQ2=FREQ2+SPAN
145 CONTINUE
    GO TO 200

C
C      ERROR MESSAGE
C
150 WRITE(108,160)
160 FORMAT(' ERROR...DATA ANALYSIS ATTEMPTED WITHOUT CALIBRAT\ -NN \I:ST
@ CHECK DATA CARD ORDER.')
    STOP
200 RETURN
    END

```

R

!JOB EH

!ASS GO-EMGGU,UD

!ASS CY-EMGAN,UP

!LOAD 7,B

!\$R001 400,,G0,1

!\$SEG 1,0,G0,1

!\$SEG 2,0,G0,3

!\$SEG 3,0,G0,1

!\$SEG 4,0,G0,4

!\$SEG 5,0,G0,2

!\$SEG 6,0,G0,5

!\$SEG 7,0,G0,1

!\$HL

END

ORIGINAL PAGE IS  
OF POOR QUALITY

MAP

OVERLAY TASK BA ORG=2100 HLUC=3AD4 CBAS=3574 CS17=428C UPEN=029F

DEF	M:FSAVE	01B3
DEF	D:KEY	1758
DEF	D:CARD	1759
DEF	D:SNAP	175A
DEF	M:SAVE	1753
DEF	M:EXIT	1754
DEF	M:IOEX	175B
DEF	M:READ	175D
DEF	M:WRITE	175E
DEF	M:CTRL	175F
DEF	M:TERM	1761
DEF	M:DATIME	1760
DEF	M:ABORT	1762
DEF	M:HEXIN	1763
DEF	M:INHEX	1764
DEF	M:CKREST	1765
DEF	M:LOAD	1755
DEF	M:OPEN	1766
DEF	M:CLOSE	1767
DEF	M:DKEYS	1768
DEF	M:WAIT	1769
DEF	M:SEGLD	176A
DEF	M:DEFINE	176B
DEF	M:ASSIGN	176C
DEF	M:BPFILE	176D
DEF	M:PUP	176E
DEF	M:RES	176F
DEF	M:DYN	1770
DEF	M:RSVP	1756
DEF	M:DOW	1757
DEF	M:COC	175C

R00T ORG=2290 LWA=264A LEN=03BB TRA=2290 SEV=0000 OV:LOAD=2415

DEF	L:88Z	L S M	263F
DEF	M:PSHC	L S M	25EC
DEF	M:PUSH	L S M	25EF
DEF	L:ERROR	L S M	257C
DEF	SEGLDX	L S M	2546
DEF	SEGLD	L S M	2544
DEF	L:88X1	L S M	24A9
DEF	L:88X	L S M	24A7
DEF	L:33R3	L S M	248F
DEF	L:33R2	L S M	2493
DEF	L:33R1	L S M	2497
DEF	M:SR	L S M	2630
DEF	L:32R4	L S M	2444
DEF	L:32L4	L S M	247D
DEF	L:32L3	L S M	2477
DEF	L:32L2	L S M	2473
DEF	L:32L1	L S M	247B
DEF	L:32R3	L S M	243E
DEF	L:32R2	L S M	243A
DEF	L:32R1	L S M	2442
DEF	L:32C	L S M	2483



DEF	::FLUAT	L S M	2480
DEF	:FLUAT	L S M	2480
DEF	OV:LOAD	I	2415
P REF	PRINT	I	0007
P REF	FORCE	I	0005
P REF	FFTPSD	I	0006
P REF	ENG	I	0005
P REF	GRAPH	I	0004
P REF	DCAL	I	0003
P REF	OUT60	I	0002
P REF	TINPUT	I	0002
P REF	CINPUT	I	0001
DEF	F:MAIN	I	2290

SEGMENT	IDENT	NODE	ORG	LWA	LEN	TRA	SEV
	0001	0000	264B	3610	0FC6	0000	0000

DEF	L:3N	L S M	35FE
DEF	L:44R2	L S M	35E4
DEF	L:44R1	L S M	35E8
DEF	L:44N2	L S M	3527
DEF	L:44M1	L S M	352F
DEF	L:44L4	L S M	3517
DEF	L:44L1	L S M	3515
DEF	L:44T4	L S M	34EE
DEF	L:44T2	L S M	34E4
DEF	L:44T1	L S M	34EC
DEF	X:3NORM	L S M	35A4
DEF	L:44S4	L S M	347A
DEF	L:44S3	L S M	3470
DEF	L:44S2	L S M	3474
DEF	L:44S1	L S M	3478
DEF	L:44A4	L S M	3486
DEF	L:44A3	L S M	347C
DEF	L:44A2	L S M	3480
DEF	L:43R3	L S M	3410
DEF	L:43R2	L S M	3419
DEF	L:43R1	L S M	3421
DEF	L:43R4	L S M	3423
DEF	L:43L4	L S M	3465
DEF	::DBLE	L S M	3416
DEF	L:43L3	L S M	345F
DEF	L:43L1	L S M	3463
DEF	L:43C	L S M	3469
DEF	:DBLE	L S M	3416
DEF	L:34R4	L S M	33C3
DEF	L:34R3	L S M	338D
DEF	L:34R1	L S M	33C1
DEF	::SNGL	L S M	3386
DEF	L:34L4	L S M	3407
DEF	L:34L3	L S M	3401
DEF	L:34L2	L S M	33FD
DEF	L:34L1	L S M	3405
DEF	L:34C	L S M	340D
DEF	:SNGL	L S M	3386
P REF	M:DEFINE	L S M	0000
P REF	M:OPEN	L S M	0000
DEF	QREW	L S M	3342
DEF	QFSKIP	L S M	3329
DEF	QWE0F	L S M	3354
DEF	QSKIP	L S M	32FA
DEF	L:43L2	L S M	345B

	DEF	L:4H	L S M	35FE
	DEF	L:44R3	L S M	35E0
	DEF	L:44A1	L S M	3484
	DEF	L:44L3	L S M	350D
	DEF	L:44L2	L S M	3511
	DEF	L:44H3	L S M	3528
	DEF	L:34R2	L S M	33B9
	DEF	L:FMT0	L S M	2F38
	DEF	L:FMT1	L S M	2DA4
	DEF	L:44M4	L S M	3531
	DEF	F:STCHR	L S M	2D6E
	DEF	F:CPWRT	L S M	2C9D
	DEF	F:FMTERR	L S M	2D30
	DEF	F:FIURET	L S M	2C13
	DEF	L:FI0	L S M	2D2B
P	REF	M:SR	L S M	0000
	DEF	L:XMT1	L S M	2AC0
	DEF	L:XMT	L S M	2AB9
	DEF	L:CLEAR	L S M	3606
	DEF	L:WSEQ	L S M	2A11
	DEF	L:IUSNER	L S M	2A36
P	REF	L:ERROR	L S M	0000
P	REF	M:PUSH	L S M	0000
	DEF	L:FI00P	L S M	2D28
	DEF	L:WSEQFN	L S M	2A9F
	DEF	L:WSEQI	L S M	2A13
	DEF	L:SEQRST	L S M	2978
	DEF	L:RSEQFN	L S M	2A9D
	DEF	L:FI0IN	L S M	2B25
	DEF	L:RSEQ	L S M	2A00
	DEF	L:SEQST	L S M	2947
	DEF	L:33M4	L S M	28C1
	DEF	L:33M2	L S M	28B7
	DEF	L:33M1	L S M	28BF
	DEF	L:33L3	L S M	28A7
	DEF	L:33L2	L S M	28AB
	DEF	L:33S3	L S M	280E
	DEF	L:33S2	L S M	2815
	DEF	L:33S1	L S M	281C
	DEF	L:33A2	L S M	2817
	DEF	L:33A1	L S M	281E
	DEF	L:33S4	L S M	280C
	DEF	L:33A4	L S M	280A
P	REF	M:POP	I	0000
P	REF	L:88X	I	0000
P	REF	L:33R2	I	0000
	DEF	L:33A3	L S M	2810
	DEF	L:33M3	L S M	28BB
	DEF	L:33L1	L S M	28AF
P	REF	L:33R3	I	0000
P	REF	::FL0AT	I	0000
	DEF	QREAD	L S M	32AC
	DEF	QCLOSE	L S M	3283
	DEF	QWRITE	L S M	325C
	DEF	QOPEN	L S M	31D8
	DEF	L:88W	L S M	2AB9
	DEF	L:88C2	L S M	2C86
	DEF	L:88S2	L S M	293E
	DEF	L:88C	L S M	2C86
	DEF	L:88S	L S M	2935
	DEF	DEFINE	L S M	3388
P	REF	M:PSHC	I	0000



DEF CINPUT I 2648

SEGMENT IDENT CODE ORG LWA LEN TRA REV  
0002 0000 2648 344E 1294 0000 0000

DEF	L:88W	L S M	3873
DEF	L:XMT	L S M	3873
DEF	L:44T4	L S M	3849
DEF	L:44T2	L S M	383F
DEF	L:44T1	L S M	3847
DEF	L:43R3	L S M	37EE
DEF	L:43R2	L S M	37EA
DEF	L:43R1	L S M	37F2
DEF	L:43R4	L S M	37F4
DEF	::DBLE	L S M	37E7
DEF	L:43L3	L S M	3830
DEF	L:43L1	L S M	3834
DEF	L:43C	L S M	383A
DEF	::DBLE	L S M	37E7
DEF	L:34R4	L S M	3794
DEF	L:34R3	L S M	378E
DEF	L:34R1	L S M	3792
DEF	::SNGL	L S M	3787
DEF	L:34L4	L S M	37D8
DEF	L:34L3	L S M	37D2
DEF	L:34L2	L S M	37CE
DEF	L:34L1	L S M	37D6
DEF	::SNGL	L S M	3787
DEF	QREW	L S M	373B
DEF	QWEUF	L S M	374D
DEF	QWRITE	L S M	3655
DEF	L:43L2	L S M	382C
DEF	L:44T3	L S M	3843
DEF	L:34R2	L S M	378A
DEF	L:FMT0	L S M	3331
DEF	L:FMT1	L S M	319D
DEF	L:XMT1	L S M	387A
DEF	F:STCHR	L S M	3167
DEF	F:CPWRT	L S M	3096
DEF	L:88C	L S M	307F
DEF	F:FNTRR	L S M	3129
DEF	F:FIURET	L S M	300C
DEF	L:FI0IN	L S M	2F1E
DEF	L:FI0	L S M	2F24
DEF	L:RSEQFN	L S M	2F02
DEF	L:CLEAR	L S M	3868
DEF	L:WSEQ	L S M	2E76
DEF	L:IUSQER	L S M	2E9B
DEF	L:SEQRST	L S M	2DEE
DEF	L:FI00P	L S M	2F21
DEF	L:WSEQFN	L S M	2F04
DEF	L:WSEQI	L S M	2E78
DEF	L:SEQST	L S M	2D8D
DEF	L:44R2	L S M	2D92
DEF	L:44R1	L S M	2D96
DEF	L:44M4	L S M	2CDF
DEF	L:44M2	L S M	2CD5
DEF	L:44L4	L S M	2CC5
DEF	L:44L3	L S M	2CBB



	DEF	L:44L2	L S M	2C8F
	DEF	L:44S4	L S M	2C51
	DEF	L:44S3	L S M	2C47
	DEF	L:44S2	L S M	2C4B
	DEF	L:44S1	L S M	2C4F
	DEF	L:44A4	L S M	2C5D
	DEF	L:44A3	L S M	2C53
	DEF	L:44A2	L S M	2C57
	DEF	L:34C	L S M	37DE
	DEF	L:4N	L S M	2D4C
	DEF	X:3NORM	L S M	2D52
P	REF	L:ERROR	L S M	0000
	DEF	L:44M3	L S M	2CD9
	DEF	L:44M1	L S M	2CD0
	DEF	L:44A1	L S M	2C5B
	DEF	L:44R3	L S M	2D8E
	DEF	L:44L1	L S M	2CC3
	DEF	L:33M4	L S M	2B21
	DEF	L:33L2	L S M	2B0B
	DEF	L:33D4	L S M	2A93
	DEF	L:33D2	L S M	2A89
	DEF	L:33D1	L S M	2A91
	DEF	L:33T3	L S M	2A76
	DEF	L:33T1	L S M	2A7A
	DEF	L:33S3	L S M	29D9
	DEF	L:33S1	L S M	29E7
	DEF	L:33A2	L S M	29E2
	DEF	L:33A1	L S M	29F9
	DEF	L:33S4	L S M	29D7
	DEF	L:33A4	L S M	29D5
	DEF	:DSIN:	L S M	2BBA
	DEF	L:43L4	L S M	3B36
	DEF	:RECNVRT	L S M	2C42
P	REF	M:PUSH	L S M	0000
P	REF	M:SR	L S M	0000
	DEF	:ABS	L S M	29AE
	DEF	L:3N	L S M	2D4C
	DEF	L:33S2	L S M	29E0
	DEF	L:33M1	L S M	2B1F
	DEF	L:33L3	L S M	2B07
	DEF	SIN	L S M	29B9
	DEF	L:33M2	L S M	2B17
	DEF	CBS.	L S M	29BB
P	REF	L:32L3	I	0000
	DEF	L:33D3	L S M	2A8D
P	REF	L:32C	I	0000
	DEF	OUT60	I	2B83
	DEF	L:33T2	L S M	2A72
	DEF	L:33A3	L S M	29DB
	DEF	L:33M3	L S M	2B1B
	DEF	L:33L1	L S M	2B0F
	DEF	::ABS	L S M	29AE
P	REF	L:33R3	I	0000
P	REF	M:POP	I	0000
	DEF	GCLUSE	L S M	367C
P	REF	L:33X	I	0000
	DEF	L:88C2	L S M	307F
	DEF	L:88S2	L S M	2DD4
	DEF	NSKIP	L S M	36F3

P	REF	L:33R2	I	0000
P	REF	::FLBAT	I	0000
	DEF	QREAD	L S M	36A5
	DEF	FIND	I	27ED
P	REF	L:88Z	I	0000
	DEF	QFSKIP	L S M	3722
	DEF	QOPEN	L S M	35D1
P	REF	M:PSHC	I	0000
	DEF	TINPUT	I	264B

SEGMENT	IDENT	NODE	ORG	LWA	LEN	TRA	SEV
	0003	0000	264B	34F9	0EAF	0000	0000

DEF	L:88W	L S M	348E
DEF	L:XMT	L S M	348E
DEF	L:44R2	L S M	3469
DEF	L:44R1	L S M	346D
DEF	L:44M2	L S M	33AC
DEF	L:44M1	L S M	33B4
DEF	L:44L4	L S M	339C
DEF	L:44L1	L S M	339A
DEF	L:44T4	L S M	3373
DEF	L:44T2	L S M	3369
DEF	L:44T1	L S M	3371
DEF	X:3NORM	L S M	3429
DEF	L:44S4	L S M	32FF
DEF	L:44S3	L S M	32F5
DEF	L:44S2	L S M	32F9
DEF	L:44S1	L S M	32F0
DEF	L:44A4	L S M	330B
DEF	L:44A3	L S M	3301
DEF	L:44A2	L S M	3305
DEF	L:43R3	L S M	32A2
DEF	L:43R2	L S M	329E
DEF	L:43R1	L S M	32A6
DEF	L:43R4	L S M	32A8
DEF	L:43L4	L S M	32EA
DEF	::DBLE	L S M	329B
DEF	L:43L3	L S M	32E4
DEF	L:43L1	L S M	32E8
DEF	L:43C	L S M	32EE
DEF	:DBLE	L S M	329B
DEF	L:34R4	L S M	3248
DEF	L:34R3	L S M	3242
DEF	L:34R1	L S M	3246
DEF	::SNGL	L S M	323B
DEF	L:34L4	L S M	328C
DEF	L:34L3	L S M	3286
DEF	L:34L2	L S M	3282
DEF	L:34L1	L S M	328A
DEF	L:34C	L S M	3292
DEF	:SNGL	L S M	323B
DEF	L:43L2	L S M	32E0
DEF	L:44T3	L S M	336D
DEF	L:44R3	L S M	3465
DEF	L:44A1	L S M	3309
DEF	L:44L3	L S M	3392
DEF	L:44L2	L S M	3396
DEF	L:44M3	L S M	33B0
DEF	L:34R2	L S M	323E
DEF	L:FMT0	L S M	2FC7
DEF	L:FMT1	L S M	2FC7

	DEF	L:43D9	L S M	3306
	DEF	L:43D1	L S M	3495
	DEF	F:STCHR	L S M	20FD
	DEF	F:CPWRT	L S M	202C
	DEF	L:88C	L S M	2015
	DEF	F:FIICPR	L S M	209F
	DEF	F:FIICST	L S M	20A0
	DEF	L:FI0IN	L S M	28B4
	DEF	L:FI0	L S M	28BA
	DEF	L:RSEQEN	L S M	2898
	DEF	L:ICLEAR	L S M	3483
	DEF	L:WSEQ	L S M	280C
	DEF	L:IOSQER	L S M	2831
P	REF	M:IPUSH	L S M	0000
	DEF	L:SEQRST	L S M	2A84
	DEF	L:FI0OP	L S M	28B7
	DEF	L:WSEQEN	L S M	289A
	DEF	L:WSEQI	L S M	280E
	DEF	L:SEQST	L S M	2A53
	DEF	L:4N	L S M	2A42
	DEF	L:33M4	L S M	29CE
	DEF	L:33M3	L S M	29C8
	DEF	L:33M1	L S M	29CC
	DEF	L:33D4	L S M	2940
	DEF	L:33D2	L S M	2936
	DEF	L:33T3	L S M	2923
	DEF	L:33T1	L S M	2927
	DEF	L:33S2	L S M	2880
	DEF	L:33S1	L S M	2894
	DEF	L:33S4	L S M	2884
	DEF	L:33A4	L S M	2882
P	REF	L:ERROR	L S M	0000
P	REF	M:SR	L S M	0000
P	REF	M:POP	I	0000
	DEF	SQRT	L S M	2826
P	REF	L:32R3	I	0000
	DEF	L:33A3	L S M	2888
	DEF	L:33M2	L S M	29C4
P	REF	L:33R2	I	0000
	DEF	L:33L2	L S M	29B8
	DEF	L:3N	L S M	2A42
	DEF	L:33D3	L S M	293A
	DEF	L:33S3	L S M	2886
	DEF	L:33D1	L S M	293E
	DEF	L:33A2	L S M	288F
P	REF	L:88X	I	0000
	DEF	L:88C2	L S M	2015
	DEF	L:88S2	L S M	2A4A
	DEF	L:33T2	L S M	291F
	DEF	L:33A1	L S M	2896
	DEF	L:33L3	L S M	29B4
P	REF	L:33R3	I	0000
	DEF	L:33L1	L S M	29BC
P	REF	L:33R1	I	0000
P	REF	::FLOAT	I	0000
P	REF	M:PSHC	I	0000
	DEF	DCAL	I	264B

SEGMENT	IDENT	MODE	ORG	LWA	LEN	TRA	SEV
	0004	0000	264B	3AD4	148A	0000	0000

DEF	M:PU5HK	L S M	3AA4
DEF	L:23E3	L S M	3A4D



	DEF	L:23E1	L S M	3A49
	DEF	:AINT	L S M	3A21
	DEF	L:AINT	L S M	3A21
	DEF	C:HARC	L S M	33CF
P	REF	M:WAIT	L S M	0000
P	REF	M:DOW	L S M	0000
	DEF	UTBM1	L S M	37C3
	DEF	BYTP	L S M	37C1
	DEF	UTB	L S M	37C4
P	REF	M:IOEX	L S M	0000
	DEF	L:23E2	L S M	3A4B
	DEF	::AINT	L S M	3A21
P	REF	M:PUSH	L S M	0000
	DEF	VCT01	L S M	35EE
	DEF	OUTNUM	L S M	3631
	DEF	MODRER	L S M	39F2
	DEF	L:4H	L S M	2EAF
	DEF	L:33M4	L S M	2E3B
	DEF	L:33M3	L S M	2E35
	DEF	L:33D4	L S M	2DAD
	DEF	L:33D3	L S M	2DA7
	DEF	L:33T3	L S M	2D90
	DEF	L:33T2	L S M	2D8C
	DEF	L:33S1	L S M	2D01
	DEF	L:33S4	L S M	2CF1
	DEF	L:33A4	L S M	2CEF
P	REF	M:SR	L S M	0000
	DEF	L:23R4	L S M	2C9F
	DEF	L:23L4	L S M	2CB3
	DEF	L:23R3	L S M	2C99
	DEF	::INT	L S M	2C92
	DEF	L:23R2	L S M	2C95
	DEF	L:23R1	L S M	2C9D
	DEF	L:23L3	L S M	2CAD
	DEF	L:23L2	L S M	2CA9
	DEF	L:23L1	L S M	2CB1
	DEF	L:23C	L S M	2CB7
	DEF	:IFIX	L S M	2C92
	DEF	:INT	L S M	2C92
	DEF	IN	L S M	374B
	DEF	SCRIBE	L S M	383F
	DEF	OUT	L S M	378B
	DEF	MODE	L S M	35BF
	DEF	LEGEND	L S M	32CE
	DEF	PASTOR	L S M	38B3
	DEF	VCTOR	L S M	35D8
	DEF	SCALE	L S M	3511
	DEF	WORDXY	L S M	3A02
	DEF	IWRDXY	L S M	39E4
	DEF	POINTS	L S M	332A
P	REF	M:POP	I	0000
	DEF	LINES	L S M	33B9
	DEF	L:33M2	L S M	2E31
	DEF	L:33D2	L S M	2DA3
	DEF	L:33A1	L S M	2D03
	DEF	L:33T1	L S M	2D94
	DEF	L:33L1	L S M	2E29
	DEF	PAGE	L S M	3447
	DEF	BELL	I	2C4B

	DEF	MODLIN	I	2B69
	DEF	LABEL	L S M	3195
	DEF	GR10	L S M	3016
	DEF	SUBJEC	L S M	2F93
P	REF	L:188Z	I	0000
	DEF	NUMBER	L S M	32FC
	DEF	WORDS	I	2BE2
	DEF	OBJECT	L S M	2F20
	DEF	MODSET	L S M	347E
	DEF	INITAL	L S M	2EB7
	DEF	L:133S3	L S M	2CF3
	DEF	L:133A2	L S M	2CFC
	DEF	L:133A3	L S M	2CF5
	DEF	L:133L3	L S M	2E21
	DEF	L:13N	L S M	2EAF
	DEF	L:133S2	L S M	2CFA
	DEF	L:133M1	L S M	2E39
	DEF	::IFIX	L S M	2C92
	DEF	L:133D1	L S M	2DAB
	DEF	L:133L2	L S M	2E25
P	REF	L:133R1	I	0000
P	REF	L:133R2	I	0000
P	REF	L:133R3	I	0000
P	REF	::FL0AT	I	0000
P	REF	M:PSHC	I	0000
	DEF	GRAPH	I	264B

SEQUENT	IDENT	NODE	ORG	LWA	LEN	TRA	SEV
	0005	0000	264B	2A88	043E	0000	0000

	DEF	L:4N	L S M	2A81
	DEF	L:133M4	L S M	2A0D
	DEF	L:133D4	L S M	297F
	DEF	L:133T1	L S M	2966
	DEF	L:133S2	L S M	28CC
	DEF	L:133S4	L S M	28C3
	DEF	L:133A4	L S M	28C1
	DEF	L:123R4	L S M	2871
	DEF	L:123L4	L S M	2885
	DEF	L:123R3	L S M	286B
	DEF	::INT	L S M	2864
	DEF	L:123R2	L S M	2867
	DEF	L:123R1	L S M	286F
	DEF	L:123L3	L S M	287F
	DEF	L:123L2	L S M	287B
	DEF	L:123L1	L S M	2883
	DEF	L:123C	L S M	2889
	DEF	::IFIX	L S M	2864
	DEF	::INT	L S M	2864
	DEF	L:13N	L S M	2A81
P	REF	M:SR	L S M	0000
	DEF	:ABS	L S M	2857
	DEF	::IFIX	L S M	2864
	DEF	L:133T2	L S M	295E
	DEF	L:133D2	L S M	2975
	DEF	FORCE	I	2787
P	REF	M:POP	I	0000
	DEF	L:133S1	L S M	28D3
	DEF	L:133A1	L S M	28D5

	DEF	L:33M2	L S M	2A03
	DEF	L:33M3	L S M	2A07
	DEF	L:33A3	L S M	28C7
	DEF	L:33T3	L S M	2962
	DEF	::ABS	L S M	2857
P	REF	L:33R2	I	0000
	DEF	L:33D1	L S M	2970
	DEF	L:33S3	L S M	28C5
	DEF	L:33L2	L S M	29F7
P	REF	L:33R1	I	0000
	DEF	L:33D3	L S M	2979
	DEF	L:33M1	L S M	2A08
P	REF	::FLBAT	I	0000
	DEF	L:33A2	L S M	28CE
	DEF	L:33L3	L S M	29F3
P	REF	L:33R3	I	0000
	DEF	L:33L1	L S M	29F8
P	REF	N:PSHC	I	0000
	DEF	ENG	I	2648

SEGMENT	IDENT	NODE	ORG	LWA	LEN	TRA	SEV
	0006	0000	2648	3978	1331	0000	0000

	DEF	L:43R3	L S M	3928
	DEF	L:43R2	L S M	3927
	DEF	L:43R1	L S M	392F
	DEF	L:43R4	L S M	3931
	DEF	::DBLE	L S M	3924
	DEF	L:43L3	L S M	3960
	DEF	L:43L2	L S M	3969
	DEF	L:43L1	L S M	3971
	DEF	L:43C	L S M	3977
	DEF	::DBLE	L S M	3924
	DEF	L:34R4	L S M	38D1
	DEF	L:34R3	L S M	38CB
	DEF	L:34R2	L S M	38C7
	DEF	L:34R1	L S M	38CF
	DEF	::SNGL	L S M	38C4
	DEF	L:34L4	L S M	3915
	DEF	L:34L3	L S M	390F
	DEF	L:34L2	L S M	3908
	DEF	L:34L1	L S M	3913
	DEF	::SNGL	L S M	38C4
	DEF	QREW	L S M	3878
	DEF	QSKIP	L S M	3830
	DEF	QREAD	L S M	37E2
	DEF	L:22E1	L S M	36AD
	DEF	::IA:5	L S M	3690
	DEF	L:44R2	L S M	3683
	DEF	L:44R1	L S M	3687
	DEF	L:44M4	L S M	35D0
	DEF	L:44M2	L S M	35C6
	DEF	L:44L4	L S M	35B6
	DEF	L:44L3	L S M	35AC
	DEF	L:44L2	L S M	35B0
	DEF	L:44S4	L S M	3542
	DEF	L:44S3	L S M	3538
	DEF	L:44S2	L S M	353C
	DEF	L:44S1	L S M	3540
	DEF	L:44A4	L S M	354E
	DEF	L:44A3	L S M	3544



DEF	L:34C	L S M	331B
DEF	L:41	L S M	36D4
DEF	X:3NORM	L S M	3643
DEF	L:44H3	L S M	35CA
DEF	L:44M1	L S M	35CE
DEF	L:44A1	L S M	354C
DEF	L:44C3	L S M	357F
DEF	L:44L1	L S M	35B4
DEF	L:33M4	L S M	3412
DEF	L:33D4	L S M	3384
DEF	L:33T1	L S M	3368
DEF	L:33S4	L S M	32C8
DEF	L:33A4	L S M	32C6
P REF	L:ERRHR	L S M	0000
DEF	:DSIN:	L S M	34AB
DEF	L:43L4	L S M	3973
DEF	:RECNVRT	L S M	3533
P REF	M:PUSH	L S M	0000
DEF	L:23R4	L S M	31F4
DEF	L:23L4	L S M	3208
DEF	L:23R3	L S M	31EE
DEF	::INT	L S M	31E7
DEF	L:23R2	L S M	31EA
DEF	L:23R1	L S M	31F2
DEF	L:23L3	L S M	3202
DEF	L:23L2	L S M	31FE
DEF	L:23L1	L S M	3206
DEF	L:23C	L S M	320C
DEF	:IFIX	L S M	31E7
DEF	:INT	L S M	31E7
P REF	M:SR	L S M	0000
DEF	:ABS	L S M	31DA
DEF	QCLOSE	L S M	37B9
DEF	QWEOF	L S M	388A
DEF	QWRITE	L S M	3792
DEF	QFSKIP	L S M	385F
DEF	QOPEN	L S M	370E
DEF	SIN	L S M	3245
P REF	L:32R3	I	0000
DEF	::IABS	L S M	369D
DEF	L:33S3	L S M	32CA
DEF	L:33S2	L S M	32D1
DEF	FORT	I	28F2
P REF	L:32L1	I	0000
DEF	::ABS	L S M	31DA
DEF	L:22E2	L S M	36A5
DEF	L:3N	L S M	36D4
DEF	COS	L S M	3247
P REF	M:POP	I	0000
DEF	SQRT	L S M	326A
DEF	L:33S1	L S M	32D8
DEF	L:33A1	L S M	32DA
DEF	L:33D1	L S M	3382
P REF	L:33R1	I	0000
DEF	L:33M1	L S M	3410
DEF	L:33D2	L S M	337A
DEF	L:33T2	L S M	3363
DEF	L:33T3	L S M	3367
DEF	L:33A2	L S M	32D3
DEF	::IFIX	L S M	31E7
DEF	L:33L3	L S M	33F8
DEF	PSDSAV	I	3143
P REF	L:33R2	I	0000

DEF	L:33A3	L S M	32CC
DEF	L:33M2	L S M	3408
DEF	L:33L2	L S M	33FC
DEF	L:33H3	L S M	340C
DEF	L:33D3	L S M	337E
DEF	L:33L1	L S M	3400
DEF	RFORT	I	294E
P REF	L:33R3	I	0000
P REF	::FLUAT	I	0000
DEF	L:22E3	L S M	36A9
DEF	WINDOW	I	28B2
P REF	M:PSHC	I	0000
DEF	FFTPSD	I	264B

SEGMENT	IDENT	NODE	ORG	LWA	LEN	TRA	SEV
	0007	0000	264B	3ABA	1470	0000	0000

DEF	L:44R2	L S M	3A96
DEF	L:44R1	L S M	3A9A
DEF	L:44M2	L S M	39D9
DEF	L:44M1	L S M	39E1
DEF	L:44L4	L S M	39C9
DEF	L:44L1	L S M	39C7
DEF	L:44T4	L S M	39A0
DEF	L:44T2	L S M	3996
DEF	L:44T1	L S M	399E
DEF	X:3NORM	L S M	3A56
DEF	L:44S4	L S M	392C
DEF	L:44S3	L S M	3922
DEF	L:44S2	L S M	3926
DEF	L:44S1	L S M	392A
DEF	L:44A4	L S M	3938
DEF	L:44A3	L S M	392E
DEF	L:44A2	L S M	3932
DEF	L:43R3	L S M	38CF
DEF	L:43R2	L S M	38CB
DEF	L:43R1	L S M	38D3
DEF	L:43R4	L S M	38D5
DEF	L:43L4	L S M	3917
DEF	::DBLE	L S M	38C8
DEF	L:43L3	L S M	3911
DEF	L:43L1	L S M	3915
DEF	L:43C	L S M	391B
DEF	:DBLE	L S M	38C8
DEF	L:34R4	L S M	3875
DEF	L:34R3	L S M	386F
DEF	L:34R1	L S M	3873
DEF	::SNGL	L S M	3868
DEF	L:34L4	L S M	38B9
DEF	L:34L3	L S M	38B3
DEF	L:34L2	L S M	38AF
DEF	L:34L1	L S M	38B7
DEF	L:34C	L S M	38BF
DEF	:SNGL	L S M	3868
DEF	L:43L2	L S M	390D
DEF	L:44T3	L S M	399A
DEF	L:44R3	L S M	3A92
DEF	L:44A1	L S M	3936
DEF	L:44L3	L S M	39BF



DEF	L:44L2	L S M	39C3
DEF	L:44L3	L S M	39D0
DEF	L:33R2	L S M	33A0
DEF	L:FHT0	L S M	35F4
DEF	L:FHT1	L S M	3460
DEF	L:44M4	L S M	39E3
DEF	F:STCHR	L S M	342A
DEF	F:CPWRT	L S M	3359
DEF	L:88C	L S M	3342
DEF	F:ENTERR	L S M	33EC
DEF	F:FI0RET	L S M	32CF
DEF	L:FI0IN	L S M	31E1
DEF	L:FI0	L S M	31E7
DEF	L:XMT1	L S M	317C
DEF	L:XMT	L S M	3175
DEF	L:RSEQFN	L S M	3159
DEF	L:CLEAR	L S M	3AB0
DEF	L:WSEQ	L S M	30C0
DEF	L:I0SQER	L S M	30F2
P REF	L:ERROR	L S M	0000
P REF	M:PUSH	L S M	0000
DEF	L:SEQRST	L S M	3045
DEF	L:FI0BP	L S M	31E4
DEF	L:WSEQFN	L S M	3158
DEF	L:WSEQI	L S M	30CF
DEF	L:SEQST	L S M	3014
DEF	L:4N	L S M	3003
DEF	L:33N4	L S M	2F8F
DEF	L:33M3	L S M	2F89
DEF	L:33D4	L S M	2F01
DEF	L:33D3	L S M	2EFB
DEF	L:33S3	L S M	2E5E
DEF	L:33S1	L S M	2E6C
DEF	L:33A1	L S M	2E6E
DEF	L:33S4	L S M	2E5C
DEF	L:33A4	L S M	2E5A
P REF	M:SR	L S M	0000
DEF	L:23R4	L S M	2E0A
DEF	L:23L4	L S M	2E1E
DEF	L:23R3	L S M	2E04
DEF	::INT	L S M	2DFD
DEF	L:23R2	L S M	2E00
DEF	L:23R1	L S M	2E08
DEF	L:23L3	L S M	2E18
DEF	L:23L2	L S M	2E14
DEF	L:23L1	L S M	2E1C
DEF	L:23C	L S M	2E22
DEF	:IFIX	L S M	2DFD
DEF	:INT	L S M	2DFD
P REF	M:P0P	I	0000
P REF	L:88X	I	0000
DEF	L:33M2	L S M	2F85
DEF	L:33A2	L S M	2E67
DEF	L:33D2	L S M	2EF7
DEF	L:33L1	L S M	2F7D
DEF	L:88W	L S M	3175
DEF	L:88C2	L S M	3342
DEF	L:88S2	L S M	300B
DEF	L:33A3	L S M	2E60



	DEF	L:33L3	L S M	2F75
	DEF	L:3N	L S M	3003
	DEF	L:33S2	L S H	2E65
	DEF	L:3301	L S H	2F80
P	REF	:::FLOAT	I	0000
	DEF	:::IFIX	L S M	2DFD
P	REF	L:33R3	I	0000
	DEF	L:33D1	L S M	2EFF
	DEF	L:33L2	L S M	2F79
P	REF	M:PSHC	I	0000
	DEF	PRINT	I	264B

ERRSEV= 0000

END MAP

APPENDIX G

NOVA SOURCE LISTING AND LOAD MAP

PROGRAM: EMGAN

AUTHOR: WILLIAM L. HURSTA

PURPOSE: TO PROVIDE FIRST ORDER REDUCTION OF EMG AND FORCE DATA.

```
EXTERNAL OV1,OV2,OV3,OV4,OV5,OV6,OV7
INTEGER DTYPE,FILNUM,FILT$W,PLOT$W
INTEGER FCAL$W,ECAL$W,WHOA,PRNT$W,OVL
DIMENSION ITFILE(4)
COMMON DATUMS(4500),IHEAD(20),CAL(2,4),IVAR$W(10)
COMMON/IGD/ MUSCLE(7,4),IFORR(2)
EQUIVALENCE (DTYPE,IVAR$W(1));(FILNUM,IVAR$W(2))
EQUIVALENCE (FILT$W,IVAR$W(3)),(NSEC$S,IVAR$W(4))
EQUIVALENCE (PLOT$W,IVAR$W(5)),(PRNT$W,IVAR$W(6))
EQUIVALENCE (ISPAN,IVAR$W(7)),(IFGAIN,IVAR$W(8))
MTO=2
OVL=0
FCAL$W=0
ECAL$W=0
WHOA=0
NSLICE=0
```

INITIALIZE MAG TAPE UNITS AND OPEN PSD OUTPUT TAPE FILE

```
CALL INIT ("MTO",0,IER)
CALL INIT ("MT1",0,IER)
CALL MTOPD (MTO,"MT1:0",0,IER)
```

OPEN OVERLAY FILE

```
CALL OVOPN (OVL,"EMGAN.OL",IER)
```

READ IN DATA CARDS

```
CALL OVLOD (OVL,OV1,0,IER)
CALL CINPT (NSLICE,I,START)
DO 200 I=1,NSLICE
```

GO GET DATA CARD INFO FOR PRESENT DATA SLICE

```
CALL OVLOD (OVL,OV1,0,IER)
CALL CINPT (NSLICE,I,START)
SPAN=FLOAT(ISPAN)
IF(DTYPE.EQ.4) FDGAIN=FLOAT(IFGAIN)
```

CHECK TO SEE IF CALIBRATION DATA ACQUIRED BEFORE PROCESSING FIRST DATA SLICE

```
IF(DTYPE.EQ.1) FCAL$W=1
IF(DTYPE.EQ.2) ECAL$W=1
IF(DTYPE.EQ.3.AND).ECAL$W.NE.1) WHOA=1
IF(DTYPE.EQ.4.AND).FCAL$W.NE.1) WHOA=1
IF(WHOA.EQ.1) GO TO 90
```

GO GET THE DATA FROM TAPE



CALL OVLOD (OVL,OV2,0,IER)  
CALL TINPUT(DTYPE,ITFILE,FILNUM,START,NSECS,I)

DECIDE WHERE TO GO FOR EACH DATA TYPE.

GO TO (20,10,30,70),DTYPE

FILTER 60 HZ FOR EMG CAL?

10 IF(FILT SW.EQ.-1) GO TO 20  
CALL OUT60 (DTYPE,NSECS,FILT SW)

CALCULATE CALIBRATION VARIABLES OF EMG OR FORCE DATA.

20 CALL OVLOD (OVL,OV3,0,IER)  
CALL DCAL (DTYPE,NSECS,IFGAIN)

PLOT EMG OR FORCE CAL DATA?

IF(PLOT SW.EQ.1) GO TO 25  
CALL OVLOD (OVL,OV4,0,IER)  
CALL GRAPH (DTYPE,START,NSECS,SPAN,STDERR)

PRINT OUT HEADER?

25 IF(PRNT SW.EQ.1) GO TO 200  
GO TO 90

FILTER 60 HZ FOR EMG DATA SLICE?

30 IF(FILT SW.EQ.-1) GO TO 40  
CALL OUT60 (DTYPE,NSECS,FILT SW)

SCALE EMG DATA, CALCULATE INTEGRATED VALUE AND MEAN SQUARE VALUE

40 CALL OVLOD (OVL,OV5,0,IER)  
CALL EMG (NSECS,VOLTSEC)

PLOT EMG DATA ?

IF(PLOT SW.EQ.1) GO TO 50  
CALL OVLOD (OVL,OV4,0,IER)  
CALL GRAPH (DTYPE,START,NSECS,SPAN,STDERR)

TAKE THE FOURIER TRANSFORM OF THE DATA AND FIND PSD UP TO 400 HZ.

50 CALL OVLOD (OVL,OV6,0,IER)  
CALL FFTPSD (NSECS,SPAN,PSDSUM,STDERR)  
DTYPE=5

PLOT THE PSD ?

IF(PLOT SW.EQ.1) GO TO 60  
CALL OVLOD (OVL,OV4,0,IER)  
CALL GRAPH (DTYPE,START,NSECS,SPAN,STDERR)

PRINT OUT RESULTS ON LINE PRINTER?

60 IF(PRNT SW.EQ.1) GO TO 200  
GO TO 90

SCALE FORCE DATA FOR PLOTTING AND AVERAGE IT FOR PRINTING.

70 CALL OVLOD (OVL,OVS,0,IER)  
CALL FORCE (DTYPE,NSECS,SPAN,FDGAIN)

PL0T FORCE DATA ?

IF (PLOTSW.EQ.1) GO TO 80  
CALL OVLOD (OVL,OVS,0,IER)  
CALL GRAPH (DTYPE,START,NSECS,SPAN,STDERR)

PRINT OUT FORCE DATA ON LINE PRINTER?

80 IF (PRNTSW.EQ.1) GO TO 200  
90 CALL OVLOD (OVL,OVS,0,IER)  
CALL POUT (DTYPE,START,NSECS,SPAN,VOLTSEC,PSDSUM,STDERR,WHOA)  
200 CONTINUE  
CALL RLSE ("MT0",IER)  
CALL RLSE ("MT1",IER)  
CALL CLOSE (OVL,IER)  
END



OVERLAY OV1  
SUBROUTINE CINPUT (NSLICE,I,START)

SUBROUTINE CINPUT ACQUIRES THE CARD DATA AND STORES IT ON A RAD FILE

ARGUMENTS: NSLICE= NUMBER OF DATA SLICES

I - IDENTIFIES PRESENT DATA SLICE BEING PROCESSED

START = BEGIN TIME OF PRESENT DATA SLICE

DIMENSION ICARD(200,15)  
COMMON DATUMS(4500),IHEAD(20),CAL(2,4),IVARSW(10)  
COMMON/IGD/ MUSCLE(7,4),IFORR(2)  
EQUIVALENCE (ICARD,DATUMS(1))  
INTEGER CDR,DSK  
LPT=12  
CDR=9  
DSK=1  
IF(NSLICE) 120,20,60

INPUT NUMBER OF DATA SLICES AND READ DATA SLICE INFO FOR PRESENT RUN

20 READ(CDR,30) NSLICE  
30 FORMAT(I5)  
WRITE(LPT,35) NSLICE  
35 FORMAT(1H1,' \*\*\*\* DATA CARDS \*\*\*\*',///,' NO. OF DATA SLICES=',I5,/  
@//,5X,'DATA SLICES',//)  
IF(NSLICE.GT.200) GO TO 100  
IF(NSLICE.LT.1) GO TO 120  
DO 50 J=1,NSLICE  
READ(CDR,40) (ICARD(J,K),K=1,15)  
40 FORMAT(1X,A2,14I3)  
WRITE(LPT,40) (ICARD(J,K),K=1,15)  
50 CONTINUE  
CALL FOPEN (DSK,"ICARD",6000)  
WRITE BINARY (DSK) ICARD  
CALL CLOSE (DSK,IER)  
GO TO 200

INITIALIZE VARIABLES AND SWITCHES FOR EACH DATA SLICE

60 CALL FOPEN (DSK,"ICARD",6000)  
READ BINARY (DSK) ICARD  
CALL CLOSE (DSK,IER)  
IVARSW(1)=ICARD(I,2)  
IVARSW(2)=ICARD(I,1)  
START=3600.\*FLOAT(ICARD(I,3))+60.\*FLOAT(ICARD(I,4))+FLOAT(ICARD(I,  
@5))  
IVARSW(3)=ICARD(I,6)  
IVARSW(4)=ICARD(I,7)  
IVARSW(5)=ICARD(I,8)  
IVARSW(6)=ICARD(I,9)  
IVARSW(7)=ICARD(I,10)  
IVARSW(8)=ICARD(I,11)  
IHEAD(16)=ICARD(I,3)  
IHEAD(17)=ICARD(I,4)  
IHEAD(18)=ICARD(I,5)  
IHEAD(19)=ICARD(I,7)  
IF(IVARSW(1)=2) 90,90,200  
90 IHEAD(20)=ICARD(I,12)  
GO TO 200



C  
C  
ERROR MESSAGES

```
100 WRITE(LPT,110) NSLICE
110 FORMAT(///,' FATAL ERROR...READ NUMBER OF DATA SLICES TO BE',I6,'.
      & CANNOT HAVE MORE THAN 200.')
      STOP
120 WRITE(LPT,130) NSLICE
130 FORMAT(///,' FATAL ERROR...READ # OF DATA CARDS TO BE NEGATIVE.
      &# OF SLICES= ',I6)
      STOP
200 IF(I.EQ.NSLICE) CALL DFILW ("ICARD",IER)
      RETURN

      END
```

ORIGINAL PAGE IS  
OF POOR QUALITY

OVERLAY 002

SUBROUTINE TINPUT(DTYPE,ITFILE,FILNUM,START,NSECS,I)

SUBROUTINE TINPUT READS FORCE AND EMG DATA IN FROM TAPE.

CALLING ARGUMENTS: DTYPE - DATA TYPE

ITFILE- FILE CURRENTLY BEING ACCESSED ON TAPE

FILNUM- FILE ON TAPE WHERE DESIRED DATA IS LOCATED

START - TIME IN TOTALED SECONDS OF BEGINNING OF  
DESIRED DATA SLICE

NSECS - NUMBER OF SECONDS OF DATA TO BE READ .

I - DATA SLICE NUMBER

INTEGER FILNUM,DTYPE

DIMENSION IARRAY(536),ITFILE(4),ITIME(3)

COMMON DATUMS(4500),IHEAD(20),CAL(2,4),IVARSW(10)

COMMON/IGD/ MUSCLE(7,4),IFORR(2)

EQUIVALENCE (IARRAY,DATUMS(4230)),(ITIME,IARRAY(2))

LPT=12

MTI=3

IREAD=1030K

IGOOD=4105K

IBREC=40000K

CHECK TO SEE IF DESIRED FILE IS CURRENTLY BEING ACCESSED.

IF NOT, EXCEPT ON FIRST DATA SLICE CLOSE PRESENT FILE BEFORE  
OPENING NEW ONE.

IF(I.EQ.1) GO TO 1

IF(FILNUM.EQ.ITFILE(3)) GO TO 5

CALL CLOSE (MTI,IER)

1 ITFILE(1)="MT"

ITFILE(2)="0:"

ITFILE(3)=FILNUM

ITFILE(4)="<><>"

POSITION TAPE AT THE BEGINNING OF THE DESIRED DATA FILE.

CALL MTOPD (MTI,ITFILE,0,IER)

IF(IER.NE.1) GO TO 300

POSITION TAPE AT THE BEGINNING OF THE DATA SLICE.

5 CALL FIND (FILNUM,START)

DECIDE WHERE TO GO FOR EACH DATA TYPE.

GO TO (10,50,50,10),DTYPE

ACQUIRE FORCE CAL OR FORCE DATA

10 JJ=1

NRECS=2\*NSECS

DO 30 J=1,NRECS

CALL MTDIO (MTI,IREAD,IARRAY,ISTAT,IER,NWRDS)

IF(IER.NE.1.OR.ISTAT.NE.IGOOD.OR.NWRDS.NE.535) GO TO 200

DO 20 K=27,535,51

DATUMS(JJ)=FLOAT(IARRAY(K))/8.



```
JJ=JJ+1
20 CONTINUE
50 CONTINUE
  IF(DTYPE.EQ.1) GO TO 1000
  GO TO 150
```

```
C
C
C      ACQUIRE EMG CAL OR EMG DATA
```

```
50 JJ=1
  NRECS=2*NSECS+1
  DO 80 J=1,NRECS
    CALL MTDIO (MTI,IREAD,IARRAY,ISTAT,IER,NWRDS)
    IF(IER.NE.1.OR.ISTAT.NE.IGOOD.OR.NWRDS.NE.535) GO TO 200
    DO 70 K=1,10
      L=25+K+50*(K-1)
      L50=L+50
      DO 60 KK=L,L50
        IF(KK.EQ.L+1) GO TO 60
        DATUMS(JJ)=FLOAT(IARRAY(KK))/8.
        JJ=JJ+1
      60 CONTINUE
    70 CONTINUE
  80 CONTINUE
    IF(DTYPE.EQ.2) GO TO 1000
150 IBACK=IBKREC+2*NSECS+2
    CALL MTDIO (MTI,IBACK,IARRAY,ISTAT,IER)
    GO TO 1000
```

```
C
C
C      ERROR MESSAGES
```

```
200 WRITE(LPT,210) DTYPE,FILNUM,ISTAT,IER,NWRDS,ITIME
210 FORMAT(///," FATAL ERROR...DURING INPUT OF DATA TYPE",I2,"IN FILE
@",A2,/,10X,"TAPE STATUS(8)=",O18,/,10X,"FORTRAN ERROR=",I7,/,10X,"
@# OF WORDS READ=",I6,/,10X,"DATA TIME=",4X,3I3)
  STOP
300 WRITE(LPT,310) IER,FILNUM
310 FORMAT(///," FATAL ERROR...FORTRAN ERROR CODE=",I3," DURING SEARCH
@ FOR FILE ",A2)
  STOP
1000 CONTINUE
  RETURN
  END
```



SUBROUTINE FIND (FILNUM,START)

SUBROUTINE FIND POSITIONS THE TAPE AT THE BEGINNING OF A DATA SLICE.

CALLING ARGUMENTS: FILNUM- FILE ON TAPE WHERE DATA SLICE IS LOCATED

START - TIME IN TOTALED SECONDS OF THE BEGINNING OF THE  
DESIRED DATA SLICE.

DIMENSION IARRAY(536)

INTEGER FILNUM

REAL NOW

COMMON DATUMS(4500),IHEAD(20),CAL(2,4),IVARSW(10)

COMMON/IGD/ MUSCLE(7,4),IFORR(2)

EQUIVALENCE (IARRAY,DATUMS(4230))

IPASS=0

LPT=12

MTI=3

IEOF=-73273K

IFFILE=30000K

IBFILE=40000K

IREAD=1030K

READ EACH RECORD AND COMPARE ITS TIME LABEL WITH THE START TIME.

10 CALL MTDIO (MTI,IREAD,IARRAY,ISTAT,IER)

IF(ISTAT.EQ.IEOF) GO TO 20

NOW=3600.\*ABS(FLOAT(IARRAY(2)))+60.\*ABS(FLOAT(IARRAY(3)))+ABS(FLOAT(IARRAY(4)))

IF(NOW.LE.START) GO TO 10

ACQUIRE THE HEADER INFORMATION FOR THE PRESENT DATA SLICE.

DO 15 J=1,15

IHEAD(J)=IARRAY(J+5)

15 CONTINUE

IBREC=IBFILE+1

CALL MTDIO (MTI,IBREC,IARRAY,ISTAT,IER)

GO TO 100

DID NOT FIND THE START TIME BEFORE ENCOUNTERING AN EOF. IF ONLY ONE PASS  
HAS BEEN MADE THRU THE DATA, BACK UP AND TRY AGAIN. OTHERWISE,STOP]

20 IF(IPASS.EQ.1) GO TO 40

DO 30 J=1,2

30 CALL MTDIO (MTI,IBFILE,IARRAY,ISTAT,IER)

CALL MTDIO (MTI,IFFILE,IARRAY,ISTAT,IER)

IPASS=1

GO TO 10

40 WRITE(LPT,50) FILNUM

50 FORMAT(//,, ' CANNOT FIND START TIME IN FILE ',A2)

STOP

100 CONTINUE

RETURN

END

SUBROUTINE OUT60 (DTYPE,NSECS,FILT5W)

SUBROUTINE OUT60 IS A DIGITAL NOTCH FILTER. DEPENDING ON THE DATA CARD REQUEST, ONLY 60 HZ IS REMOVED OR ALL HARMONICS OF 60 HZ UP TO 360 HZ

CALLING ARGUMENTS: \* DTYPE = TYPE OF DATA

NSECS = NUMBER OF SECS OF DATA TO BE FILTERED

FILT5W= SWITCH WHICH INDICATES WHETHER ONLY 60 HZ IS TO BE REMOVED OR ALL HARMONICS OF 60 HZ

```
INTEGER DTYPE, FILT5W
COMMON DATUMS(4500), IHEAD(20), CAL(2,4), IVAR5W(10)
COMMON/IGD/ MUSCLE(7,4), IFORR(2)
EQUIVALENCE (IRATE, IHEAD(11))
LOOPS=1
TWOPI=6.28318530717
IF(FILT5W.EQ.1) LOOPS=5
NPTS=1024*NSECS
DO 30 J=1, LOOPS, 2
  A=0.0
  B=0.0
  W=(60*J)*TWOPI/FLOAT(IRATE)
  DO 10 K=1, NPTS
    A=A+DATUMS(K)*COS(K*W)
    B=B+DATUMS(K)*SIN(K*W)
10  CONTINUE
  A=A*2./FLOAT(NPTS)
  B=B*2./FLOAT(NPTS)
  DO 20 K=1, NPTS
    DATUMS(K)=DATUMS(K)-A*COS(K*W)-B*SIN(K*W)
20  CONTINUE
30  CONTINUE
  RETURN
END
```

ORIGINAL PAGE IS  
OF POOR QUALITY



OVERLAY UV3  
SUBROUTINE DCAL(DTYPE,NSECS,IFGAIN)

SUBROUTINE DCAL CALCULATES THE SCALE FACTORS FOR THE FORCE AND EMG DATA  
THE CALIBRATION CONSTANT FOR INTEGRATED EMG AREA IS ALSO FOUND. EMG  
THE CAL DATA IS THEN CONVERTED INTO UNITS OF POUNDS OR MICROVOLTS FOR  
PLOTING ON THE COMPUTEK TERMINAL.

CALLING ARGUMENTS: DTYPE - TYPE OF DATA

NSECS - LENGTH OF CAL DATA

IFGAIN- GAIN APPLIED TO FORCE CAL SIGNAL

DIMENSION HIGH(10),LOW(10),POUNDS(2600)  
INTEGER DTYPE,FLBCAL  
REAL LOWSUM  
COMMON DATUMS(4500),IHEAD(20),CAL(2,4),IVARSW(10)  
COMMON/IGD/ MUSCLE(7,4),IFORR(2)  
EQUIVALENCE (FLBCAL,IHEAD(15)),(IRATE,IHEAD(11))  
EQUIVALENCE (MAGTUD,IHEAD(13)),(POUNDS,DATUMS(1))  
LPT=12  
IF(DTYPE.EQ.2) GO TO 100

FIND FORCE CALIBRATION VARIABLES

CAL(1,1)=FLOAT(FLBCAL)  
CAL(1,4)=FLOAT(IFGAIN)

LOOK FOR JUMP IN DATA TO INDICATE FORCE CAL.

BEGIN=DATUMS(80)  
NPTS=20\*NSECS-80  
DO 10 J=81,NPTS  
IF(DATUMS(J).GT.BEGIN+500.) GO TO 30  
10 CONTINUE  
WRITE (LPT,20)  
20 FORMAT(///,'COULD NOT FIND THE FORCE CAL.---CHECK DATA SLICE TIMES  
@')  
STOP

FOUND THE JUMP, NOW COMPUTE THE AVERAGE BEFORE AND AFTER THE JUMP.  
CAL WEIGHT IN COUNTS = AVERAGE AFTER - AVERAGE BEFORE  
FORCE SLOPE = COUNTS/POUND  
ZERO FORCE = AVERAGE COUNTS JUST BEFORE CAL JUMP

30 HISUM=0.0  
LOWSUM=0.0  
DO 40 K=1,20  
HISUM=HISUM+DATUMS(J+K+60)  
LOWSUM=LOWSUM+DATUMS(J-K-60)  
40 CONTINUE  
FORSLP=(HISUM/20.-LOWSUM/20.)/FLOAT(FLBCAL)  
CAL(1,2)=FORSLP  
ZEROFOR=LOWSUM/20.  
CAL(1,3)=ZEROFOR

SCALE FORCE CALIBRATION DATA

NPTS=20\*NSECS  
DO 50 J=1,NPTS  
POUNDS(J)=(DATUMS(J)-ZEROFOR)/FORSLP



```

50 CONTINUE
GO TO 200

C
C   FIND EMG CALIBRATION VARIABLES
C   FIRST REMOVE ANY DC OFFSET IN THE CAL DATA
C

100 SUM=0.0
    INDEX=1000*NSECS
    DO 110 J=1,INDEX
        SUM=SUM+DATUMS(J)
110 CONTINUE
    OFFSET=SUM/FLOAT(INDEX)
    DO 120 J=1,INDEX
        DATUMS(J)=DATUMS(J)-OFFSET
120 CONTINUE

C
C   CALCULATE STANDARD DEVIATION AND OBTAIN SCALE FACTOR. FOR ZERO MEAN,
C   STATIONARY SIGNAL STANDARD DEVIATION=RMS VALUE
C

    XSQUAR=0.0
    DO 140 J=1,INDEX
        XSQUAR=XSQUAR+DATUMS(J)**2
140 CONTINUE
    EMGCAL=SQRT(XSQUAR/(INDEX-1))/FLOAT(MAGTUD)
    CAL(2,1)=FLOAT(MAGTUD)
    CAL(2,2)=EMGCAL

C
C   SCALE CAL DATA FOR PLOTTING.
C

    DO 180 J=1,INDEX
        DATUMS(J)=DATUMS(J)/EMGCAL
180 CONTINUE
200 RETURN
END

```

OVERLAY OV4

SUBROUTINE GRAPH (DTYPE, START, NSECS, SPAN, STDERR)

SUBROUTINE GRAPH PLOTS ON THE COMPUTEK TERMINAL THE DATA SLICE BEING  
PROCESSED ALONG WITH VARIOUS HEADER INFORMATION.

CALLING ARGUMENTS: DTYPE - TYPE OF DATA

START - BEGIN TIME OF PRESENT DATA SLICE

NSECS - LENGTH OF DATA SLICE IN SECONDS

SPAN - AVERAGING INTERVAL FOR PSD

STDERR- NORMALIZED STANDARD ERROR FOR PSD ESTIMA

DIMENSION DATE(3), POUNDS(2600), PSDN(410)

DIMENSION FREQ(410), CALDATA(2000)

INTEGER DTYPE

COMMON DATUMS(4500), IHEAD(20), CAL(2,4), IVARSW(10)

COMMON/IGD/ MUSCLE(7,4), IFORR(2)

EQUIVALENCE (POUNDS, DATUMS(1)), (PSDN, DATUMS(4000))

EQUIVALENCE (FREQ, DATUMS(3500)), (FLBCAL, CAL(1,1))

EQUIVALENCE (CALDATA, DATUMS(1)), (EMGMAX, DATUMS(4499))

EQUIVALENCE (FORMAX, DATUMS(4500))

RETURN

END



OVERLAY OV5  
SUBROUTINE EMG (NSECS,VOLTSEC)

SUBROUTINE EMG SCALES EMG DATA FOR PLOTTING ON THE COMPUTER AND FINDS  
THE INTEGRATED EMG VALUE FOR THE DATA SLICE.

CALLING ARGUMENTS: NSECS - LENGTH OF DATA SLICE

VOLTSEC- INTEGRATED EMG VALUE IN MICROVOLT\*SEC.

COMMON DATUMS(4500),IHEAD(20),CAL(2,4),IVARSW(10)  
COMMON/IGD/ MUSCLE(7,4),IFORR(2)  
EQUIVALENCE (IRATE,IHEAD(11)),(EMGCAL,CAL(2,2))  
EQUIVALENCE (EMGMAX,DATUMS(4499)),(EMGVAR,DATUMS(4500))

CALCULATE OFFSET, SUBTRACT FROM DATA AND APPLY SCALE FACTOR

SUM=0.0  
NPTS=1024\*NSECS  
DO 10 J=1,NPTS  
SUM=SUM+DATUMS(J)  
10 CONTINUE  
OFFSET=SUM/FLOAT(NPTS)  
CAL(2,4)=OFFSET  
EMGMAX=0.0  
DO 20 J=1,NPTS  
DATUMS(J)=(DATUMS(J)-OFFSET)/EMGCAL  
IF(ABS(DATUMS(J)).GT.EMGMAX) EMGMAX=DATUMS(J)  
20 CONTINUE

FIND INTEGRATED EMG VALUE

NPTS=NSECS\*1024-1  
H=1./FLOAT(IRATE)  
SUMO=0.0  
SUME=0.0  
NPTS2=NPTS-2  
DO 40 I=1,NPTS2,2  
SUME=SUME+ABS(DATUMS(I+1))  
SUMO=SUMO+ABS(DATUMS(I))  
40 CONTINUE  
VOLTSEC=(2.0\*SUMO+4.0\*SUME-ABS(DATUMS(1))+ABS(DATUMS(NPTS)))\*H/3.0

CALCULATE THE VARIANCE OF THE DATA

EMGVAR=0.0  
DO 50 J=1,NPTS  
EMGVAR=EMGVAR+DATUMS(J)\*\*2  
50 CONTINUE  
EMGVAR=EMGVAR/FLOAT(NPTS-1)  
RETURN  
END



SUBROUTINE FORCE (DTYPE,NSECS,SPAN,FDGAIN)

SUBROUTINE FORCE SCALES ALL FORCE DATA FOR PLOTTING ON THE COMPUTER AND  
AVERAGES FORCE DATA (OVER INTERVALS DETERMINED BY DATA CARD REQUEST)  
FOR PRINTED OUTPUT.

CALLING ARGUMENTS: DTYPE - TYPE OF DATA

NSECS - LENGTH OF DATA SLICE

SPAN - AVERAGING INTERVAL FOR PRINTED OUTPUT

FDGAIN - AMOUNT OF GAIN WHICH WAS APPLIED TO FORCE  
DATA

DIMENSION POUNDS(2600),AVGLBS(130)  
INTEGER DTYPE  
COMMON DATUMS(4500),IHEAD(20),CAL(2,4),IVARSW(10)  
COMMON/IGD/ MUSCLE(7,4),IFORR(2)  
EQUIVALENCE (POUNDS,DATUMS(1)),(AVGLBS,DATUMS(4000))  
EQUIVALENCE (FCGAIN,CAL(1,4)),(ZEROFOR,CAL(1,3))  
EQUIVALENCE (FORSLP,CAL(1,2)),(FORMAX,DATUMS(4500))  
GAIN=FCGAIN/FDGAIN

FIND ZERO BASELINE BEFORE BEGINNING OF FORCE DATA

SUM=0.0  
DO 30 J=1,20  
SUM=SUM+POUNDS(J)  
30 CONTINUE  
BASE=SUM/20.

SCALE FORCE DATA INTO POUNDS

FORMAX=0.0  
NPTS=20\*NSECS  
DO 40 J=1,NPTS  
POUNDS(J)=(POUNDS(J)-BASE)\*GAIN/FORSLP  
IF(POUNDS(J).GT.FORMAX) FORMAX=POUNDS(J)  
40 CONTINUE

AVERAGE FORCE DATA OVER SPAN INTERVAL

LOOPS=NSECS/IFIX(SPAN)  
NPTS=20\*IFIX(SPAN)  
DO 60 J=1,LOOPS  
SUM=0.0  
DO 50 K=1,NPTS  
SUM=SUM+POUNDS((J-1)\*NPTS+K)  
50 CONTINUE  
AVGLBS(J)=SUM/FLOAT(NPTS)  
60 CONTINUE  
RETURN  
END

ORIGINAL PAGE IS  
OF POOR QUALITY

OVERLAY 0V6  
SUBROUTINE FFIPSD (NSECS,SPAN,PSDSUM,STDERR)

SUBROUTINE FFIPSD TAKES THE FOURIER TRANSFORM OF EMG DATA AND CALCULATES THE RAW POWER SPECTRAL DENSITY UP TO 400 HZ. ACCORDING TO THE DATA REQUEST A SMOOTHED PSD IS CALCULATED WHICH HAS IMPROVED STATISTICAL PROPERTIES. TO REFLECT THOSE STATISTICAL PROPERTIES THE NORMALIZED STANDARD ERROR IS ALSO CALCULATED.

CALLING ARGUMENTS: NSECS - LENGTH OF DATA SLICE

SPAN - BANDWIDTH OF SMOOTHED PSD

PSDSUM- AREA UNDER POWER SPECTRUM

STDERR- NORMALIZED STANDARD ERROR OF PSD

DIMENSION PSD(2048),PSDN(410),PERCNT(410)  
DIMENSION CUMPCT(410),S(1024)  
COMMON DATUMS(4500),IHEAD(20),CAL(2,4),IVARSW(10)  
COMMON/IGD/ MUSCLE(7,4),IFORR(2)  
EQUIVALENCE (PSD,DATUMS(1))  
EQUIVALENCE (IRATE,IHEAD(11)),(PSDN,DATUMS(4000))  
EQUIVALENCE (PERCNT,DATUMS(3000)),(CUMPCT,DATUMS(3500))  
EQUIVALENCE (EXPVAL,DATUMS(2500)),(STDEV,DATUMS(2501))  
CONTINUE

APPLY COSINE WINDOW TO DATA

CALL WINDOW (NSECS)

DETERMINE THE POWER OF 2 AND CORRESPONDING NUMBER OF DATA POINTS

NEXP=10+NSECS/2  
IF(NEXP.EQ.14) NEXP=13  
N=2\*\*NEXP  
PTS=FLOAT(N)

TAKE FFT OF THE DATA

CALL RFFT (DATUMS,NEXP,S,-1,IFERR)

CALCULATE RAW POWER SPECTRAL DENSITY

H=1./FLOAT(IRATE)  
K=1  
N1=N+1  
DO 10 J=3,N1,2  
PSD(K)=8.\*H/PTS\*(DATUMS(J)\*\*2+DATUMS(J+1)\*\*2)  
K=K+1

10 CONTINUE

CORRECT PSD AMPLITUDES FOR WINDOW REDUCTION

N2=N/2  
DO 15 J=1,N2  
PSD(J)=1.1429\*PSD(J)

15 CONTINUE

OUTPUT RAW PSD TO TAPE

CALL PSDSAVE (NSECS)



# CALCULATE SMOOTHED PSD

```

C
C
F0=1./(PTS/FLOAT(IRATE))
IOTA=IFIX(SPAN/F0)
SPAN=FLOAT(IOTA)*F0
PSDMAX=0.0
K=1
N2=N2/IOTA*IOTA
DO 30 J=1,N2,IOTA
TEMP=0.0
DO 20 JJ=1,IOTA
TEMP=TEMP+PSD(J+JJ-1)
20 CONTINUE
PSD(K)=TEMP/FLOAT(IOTA)
IF(FLOAT(J/IOTA)*SPAN.GT.400.+SPAN) GO TO 25
IF(PSD(K).GT.PSDMAX) PSDMAX=PSD(K)
25 K=K+1
30 CONTINUE

```

## CALCULATE THE NORMALIZED PSD

```

C
C
LIM=IFIX(400./SPAN)+1
DO 40 J=1,LIM
PSDN(J)=PSD(J)/PSDMAX
40 CONTINUE

```

## INTEGRATE THE SMOOTHED PSD

```

C
C
PSDSUM=0.0
DO 50 J=1,LIM
PSDSUM=PSDSUM+SPAN*PSD(J)
50 CONTINUE

```

## CALCULATE % EACH BANDWIDTH IS OF TOTAL AND FIND THE CUMLATIVE % OF TOT

```

C
C
DO 60 J=1,LIM
PERCNT(J)=SPAN*PSD(J)/PSDSUM*100.
IF(J.NE.1) GO TO 55
CUMPCT(1)=PERCNT(1)
GO TO 60
55 CUMPCT(J)=CUMPCT(J-1)+PERCNT(J)
60 CONTINUE

```

## CALCULATE THE EXPECTED VALUE AND VARIANCE OF THE PSD

```

C
C
F=SPAN/2.
EXPVAL=0.0
DO 70 J=1,LIM
EXPVAL=EXPVAL+(PERCNT(J)/100.)*F
F=F+SPAN
70 CONTINUE
F=SPAN/2.
VARVAL=0.0
DO 80 J=1,LIM
VARVAL=VARVAL+(PERCNT(J)/100.)*(F-EXPVAL)**2
F=F+SPAN
80 CONTINUE
SIDEV=SQRT(VARVAL)

```

## FIND THE NORMALIZED STANDARD ERROR

```

C
C
STDERR=SQRT(1./FLOAT(IOTA))
RETURN

```



END

ORIGINAL PAGE IS  
OF POOR QUALITY

SUBROUTINE WINDOW (NSECS)

SUBROUTINE WINDOW APPLIES A COSINE TAPER TO THE FIRST AND LAST TENTHS  
THE DATA TO REDUCE LEAKAGE.

CALLING ARGUMENTS: NSECS = LENGTH OF THE DATA

COMMON DATUMS(4500), IHEAD(20), CAL(2,4), IVARSW(10)

COMMON/IGD/ MUSCLE(7,4), IFORR(2)

PI=3.1415927

NPTS=1024\*NSECS

IEDGE=NPTS/10

TTOTAL=1.024\*FLOAT(NSECS)

K=1

DO 10 J=1, IEDGE

T1=.001\*FLOAT(J)

C1=.5\*(1.-COS(PI\*T1/(.1\*TTOTAL)))

DATUMS(J)=C1\*DATUMS(J)

DATUMS(NPTS-J+1)=C1\*DATUMS(NPTS-J+1)

10 CONTINUE

RETURN

END



SUBROUTINE RFFT(A,M,S,IFS,IFERR)  
ONE-DIMENSIONAL REAL FINITE FOURIER TRANSFORM

FOURIER TRANSFORM SUBROUTINE FOR REAL DATA,  
PROGRAMMED IN SYSTEM/360, BASIC PROGRAMMING SUPPORT,  
FORTRAN IV, (SEE FORM C28-6504).  
THIS DECK IS SET UP FOR IBSYS ON THE IBM 7094

THIS PROGRAM USES THE SUBROUTINE FFT TO COMPUTE COMPLEX  
FOURIER TRANSFORMS OF REAL DATA. PK FORT S.D.A. NO. 3465 IS  
AVAILABLE THROUGH SHARE.

THE FOURIER SERIES IS

$$X(J) = \text{SUM OVER } K=0 \text{ TO } N, \text{ OF } C(K) * \exp(2 * \pi * I * J * K / N)$$
  
$$J=0,1,2,\dots,N-1$$

WHERE  $I = \sqrt{-1}$  AND WHERE  $C(K)$  IS COMPLEX.  
SINCE  $X(J)$  IS REAL,  $C(K) = \text{CONJG}(C(N-K))$ . THEREFORE ONLY  
 $C(K), K=0,1,\dots,N/2$  ARE COMPUTED AND/OR USED.

ARGUMENTS-

A IS INITIALLY THE INPUT ARRAY, X, WHEN COMPUTING A FOURIER  
TRANSFORM AND C WHEN COMPUTING A FOURIER SERIES. A IS REPLACES BY  
THE OUTPUT ARRAY, C IN THE FORMER CASE, X IN THE LATTER.  
THE X VECTOR CONTAINS THE REAL DATA  $X(0), X(1), \dots, X(N-1)$   
THE C VECTOR CONTAINS THE COMPLEX FOURIER AMPLITUDES  
 $C(0), C(1), \dots, C(N/2)$ . THE COMPLEX VECTOR C IS STORED ACCORDING  
TO THE NORMAL FORTRAN IV CONVENTION FOR STORING COMPLEX NUMBERS.  
I.E., REAL PARTS IN ALTERNATE CELLS STARTING WITH THE FIRST,  
IMAGINARY PARTS IN ALTERNATE CELLS STARTING WITH THE SECOND.  
TO ADHERE TO FORTRAN RULES,  $X(0), X(1), \dots$ , ARE REFERRED TO AS  
 $X(1), X(2), \dots$ , RESP. IN THE PROGRAMS. ALSO,  $C(0), C(1), \dots$  ARE  
REFERRED TO AS  $C(1), C(2), \dots$ , RESP., IF C IS DESIGNATED AS  
COMPLEX IN A TYPE STATEMENT.

M GIVES  $N=2**M$

THE ARGUMENTS S, IFS, AND IFERR ARE THE SAME AS IN THE  
SUBROUTINE FFT AND THE USER IS REFERRED TO THE COMMENT CARDS  
IN FFT FOR THEIR EXPLANATION.  
DIMENSION STATEMENTS- THE DIMENSIONS OF ARRAYS A AND S SHOULD  
BE  $N+2$  AND  $N/4$ , RESP. FOR THE LARGEST N TO BE USED. FOR  
EXAMPLE, IF THE LARGEST M IS 13, THEN,  $N=8192$  AND ONE SHOULD  
HAVE THE DIMENSION STATEMENT-  
DIMENSION A(8194), S(2048)  
IF ONE WISHES TO SPECIFY A TO BE COMPLEX BY A TYPE STATEMENT,  
ONE SHOULD GIVE IT A DIMENSION OF  $N/2 + 1$ , FOR THE LARGEST N.

DIMENSION A(4500),S(1024)

IFERRS = 0

N=2\*\*M

NV2 = N / 2

NV4M1 = N/4 - 1

MM1 = M - 1

IF (IABS(IFS)-1) 40,40,20

20 IF (MP-M) 30,50,50

30 IFERRS = 1

40 NP = N

MP = M



```

CALL FFT (A,M,S,0,IFERR1)
IFERRS = IFERRS + IFERR1
50 KD = NP / N
KT = KD
NPV4 = NP / 4
IF (IFS) 60,80,90
60 CALL FFT(A,MM1,S,-2,IFERR2)
IFERRS = IFERRS + IFERR2
DO 70 K=1,NV4M1
J=NV2-K
A1R=A(2*K+1) + A(2*J+1)
A1I=A(2*K+2)-A(2*J+2)
A2R=A(2*K+2)+A(2*J+2)
A2I=A(2*J+1)-A(2*K+1)
KKT = NPV4-KT
AWR=A2R*S(KKT) + A2I*S(KT)
AWI = A2I*S(KKT)-A2R*S(KT)
A(2*K+1)=(A1R+AWR)/4.
A(2*K+2)=(A1I+AWI)/4.
A(2*J+1)=(A1R-AWR)/4.
A(2*J+2)=(AWI-A1I)/4.
70 KT=KT+KD
T=A(1)
A(1)=(T+A(2))/2.
A(N+1) = (T-A(2))/2.
A(2)=0.
A(N+2) = 0.
A(NV2+1) = .5*A(NV2+1)
A(NV2+2) = -.5 * A(NV2+2)
80 IFERR = IFERRS
RETURN
90 DO 100 K=1,NV4M1
J=NV2-K
A1R=A(2*K+1) + A(2*J+1)
A1I=A(2*K+2)-A(2*J+2)
AWR=A(2*K+1)-A(2*J+1)
AWI=A(2*K+2)+A(2*J+2)
KKT = NPV4 - KT
A2R=AWR*S(KKT) - AWI*S(KT)
A2I=AWR*S(KT) + AWI*S(KKT)
A(2*K+1) = A1R - A2I
A(2*K+2) = A1I + A2R
A(2*J+1) = A1R + A2I
A(2*J+2) = A2R - A1I
100 KT = KT + KD
T = A(1)
A(1) = T + A(N+1)
A(2) = T - A(N+1)
A(NV2+1) = 2.*A(NV2+1)
A(NV2+2) = -2.*A(NV2+2)
CALL FFT(A,MM1,S,2,IFERR2)
IFERRS = IFERRS+IFERR2
GO TO 80
END

```

SUBROUTINE FFT(A,M,S,IFS,IFERR)  
FORT, ONE-DIMENSIONAL FINITE COMPLEX FOURIER TRANSFORM.

FOURIER TRANSFORM SUBROUTINE, PROGRAMMED IN SYSTEM/360,  
BASIC PROGRAMMING SUPPORT, FORTRAN IV, FORT 020-6504  
THIS DECK SET UP FOR IBSYS ON IBM 7094.

DOES EITHER FOURIER SYNTHESIS, I.E., COMPUTES COMPLEX FOURIER SERIES  
GIVEN A VECTOR OF N COMPLEX FOURIER AMPLITUDES, OR, GIVEN A VECTOR  
OF COMPLEX DATA X DOES FOURIER ANALYSIS, COMPUTING AMPLITUDES.  
A IS A COMPLEX VECTOR OF LENGTH  $N=2**M$  COMPLEX NOS. OR  $2*N$  REAL  
NUMBERS. A IS TO BE SET BY USER.

M IS AN INTEGER 0.LT.M.LE.13, SET BY USER.

S IS A VECTOR  $S(J)=\sin(2*PI*J/NP)$ ,  $J=1,2,\dots,NP/4-1$ ,  
COMPUTED BY PROGRAM.

IFS IS A PARAMETER TO BE SET BY USER AS FOLLOWS-

IFS=0 TO SET  $NP=2**M$  AND SET UP SINE TABLE.

IFS=1 TO SET  $N=NP=2**M$ , SET UP SIN TABLE, AND DO FOURIER  
SYNTHESIS, REPLACING THE VECTOR A BY

$X(J)=\text{SUM OVER } K=0,N-1 \text{ OF } A(K)*\exp(2*PI*I/N)**(J*K)$ ,  
 $J=0,N-1$ , WHERE  $I=\sqrt{-1}$

THE X'S ARE STORED WITH RE  $X(J)$  IN CELL  $2*J+1$   
AND IM  $X(J)$  IN CELL  $2*J+2$  FOR  $J=0,1,2,\dots,N-1$ .  
THE A'S ARE STORED IN THE SAME MANNER.

IFS=-1 TO SET  $N=NP=2**M$ , SET UP SIN TABLE, AND DO FOURIER  
ANALYSIS, TAKING THE INPUT VECTOR A AS X AND  
REPLACING IT BY THE A SATISFYING THE ABOVE FOURIER SERIES.

IFS=0 SET UP SIN TABLE AND RETURN

IFS=+2 TO DO FOURIER SYNTHESIS ONLY, WITH A PRE-COMPUTED S.

IFS=-2 TO DO FOURIER ANALYSIS ONLY, WITH A PRE-COMPUTED S.

IFERR IS SET BY PROGRAM TO-

=0 IF NO ERROR DETECTED.

=1 IF M IS OUT OF RANGE., OR, WHEN IFS=+2,-2, THE  
PRE-COMPUTED S TABLE IS NOT LARGE ENOUGH.

=-1 WHEN IFS =+1,-1, MEANS ONE IS RECOMPUTING S TABLE  
UNNECESSARILY.

NOTE- AS STATED ABOVE, THE MAXIMUM VALUE OF M FOR THIS PROGRAM  
ON THE IBM 7094 IS 13. FOR 360 MACHINES HAVING GREATER STORAGE  
CAPACITY, ONE MAY INCREASE THIS LIMIT BY REPLACING 13 IN  
STATEMENT 3 BELOW BY  $\log_2 N$ , WHERE N IS THE MAX. NO. OF  
COMPLEX NUMBERS ONE CAN STORE IN HIGH-SPEED CORE. ONE MUST  
ALSO ADD MORE DO STATEMENTS TO THE BINARY SORT ROUTINE  
FOLLOWING STATEMENT 24 AND CHANGE THE EQUIVALENCE STATEMENTS  
FOR THE K'S.

DIMENSION A(4500),S(1024),K(15)

IF (M) 20,20,10

10 IF (M-15) 40,40,20

20 IFERR=1

30 RETURN

40 IFERR=0

$N=2**M$

IF (IABS(IFS)-1) 440,440,50

WE ARE DOING TRANSFORM ONLY. SEE IF PRE-COMPUTED  
S TABLE IS SUFFICIENTLY LARGE

50 IF (N-NP) 70,70,60

60 IFERR=1



```

GO TO 440
C   SCRAMBLE A, BY SANDE'S METHOD
70  K(1)=2*N
    DO 80 L=2,M
80  K(L)=K(L-1)/2
    DO 90 L=M,14
90  K(L+1)=2

C   THE FOLLOWING 15 STATEMENTS ARE TO COMPENSATE FOR A WEAKNESS IN
C   THE FORTRAN V COMPILER
    K1 =K(15)
    K2 =K(14)
    K3 =K(13)
    K4 =K(12)
    K5 =K(11)
    K6 =K(10)
    K7 =K(9)
    K8 =K(8)
    K9 =K(7)
    K10=K(6)
    K11=K(5)
    K12=K(4)
    K13=K(3)
    K14=K(2)
    K15=K(1)
    N2 =K(1)

C   NOTE EQUIVALENCE OF KL AND K(14-L)
C   BINARY SORT-
    IJ =2
    J1 =2
110  J2 =J1
120  J3 =J2
130  J4 =J3
140  J5 =J4
150  J6 =J5
160  J7 =J6
170  J8 =J7
180  J9 =J8
190  J10=J9
200  J11=J10
210  J12=J11
220  J13=J12
230  J14=J13
240  J1=J14
250  IF (IJ-J1) 260,270,270
260  T=A(IJ-1 )
    A(IJ-1)=A(J1-1)
    A(J1-1)=T
    T=A(IJ)
    A(IJ)=A(J1)
    A(J1)=T
270  IJ=IJ+2
    J1=J1+K14
    IF (J1.LE.K15) GO TO 250
    J14=J14+K13
    IF (J14.LE.K14) GO TO 240
    J13=J13+K12
    IF (J13.LE.K13) GO TO 230
    J12=J12+K11
    IF (J12.LE.K12) GO TO 220
    J11=J11+K10
    IF (J11.LE.K11) GO TO 210
    J10=J10+K9
    IF (J10.LE.K10) GO TO 200

```



```

J9=J9+K8
IF (J9.LE.K9) GO TO 190
J8=J8+K7
IF (J8.LE.K8) GO TO 180
J7=J7+K6
IF (J7.LE.K7) GO TO 170
J6=J6+K5
IF (J6.LE.K6) GO TO 160
J5=J5+K4
IF (J5.LE.K5) GO TO 150
J4=J4+K3
IF (J4.LE.K4) GO TO 140
J3=J3+K2
IF (J3.LE.K3) GO TO 130
J2=J2+K1
IF (J2.LE.K2) GO TO 120
J1=J1+2
IF (J1.LE.K1) GO TO 110
IF (IFS) 280,20,300
C DOING FOURIER ANALYSIS,SO DIV. BY N AND CONJUGATE.
280 FN = N
DO 290 I=1,N
A(2*I-1)=A(2*I-1)
290 A(2*I)=-A(2*I)
C SPECIAL CASE= L=1
300 DO 310 I=1,N,2
T = A(2*I-1)
A(2*I-1) = T + A(2*I+1)
A(2*I+1)=T-A(2*I+1)
T=A(2*I)
A(2*I) = T + A(2*I+2)
310 A(2*I+2)= T - A(2*I+2)
IF (M-1) 20,30,320
C SET FOR L=2
320 LEXP1=2
C LEXP1=2**(L-1)
LEXP=8
C LEXP=2**(L+1)
NPL= 2**MT
C NPL = NP* 2**L
DO 390 L=2,M
C SPECIAL CASE= J=0
DO 340 I=2,N2,LEXP
I1=I + LEXP1
I2=I1+ LEXP1
I3 =I2+LEXP1
T=A(I-1)
A(I-1) = T +A(I2-1)
A(I2-1) = T-A(I2-1)
T =A(I)
A(I) = T+A(I2)
A(I2) = T-A(I2)
T= -A(I3)
TI = A(I3-1)
A(I3-1) = A(I1-1) - T
A(I3 ) = A(I1 ) - TI
A(I1-1) = A(I1-1) +T
340 A(I1) = A(I1 ) +TI
IF (L-2) 380,380,350
350 KLAST=N2-LEXP
JJ=NPL
DO 370 J=4,LEXP1,2
NPJJ=N1-JJ

```

```

UR=S(NPJJ)
UI=S(JJ)
ILAST=J+KLAST
DO 360 I=J,ILAST,LEXP
  I1=I+LEXP1
  I2=I1+LEXP1
  I3=I2+LEXP1
  T=A(I2-1)*UR-A(I2)*UI
  TI=A(I2-1)*UI+A(I2)*UR
  A(I2-1)=A(I1-1)-T
  A(I2)=A(I1)-TI
  A(I1-1)=A(I1-1)+T
  A(I1)=A(I1)+TI
  T=-A(I3-1)*UI-A(I3)*UR
  TI=A(I3-1)*UR-A(I3)*UI
  A(I3-1)=A(I1-1)-T
  A(I3)=A(I1)-TI
  A(I1-1)=A(I1-1)+T
360 A(I1)=A(I1)+TI
C   END OF I LOOP
370 JJ=JJ+NPL
C   END OF J LOOP
380 LEXP1=2*LEXP1
  LEXP = 2*LEXP
390 NPL=NPL/2
C   END OF L LOOP
  IF (IFS) 410,20,30
C   DOING FOURIER ANALYSIS. REPLACE A BY CONJUGATE.
410 DO 420 I=1,N
420 A(2*I)=-A(2*I)
  GO TO 30
C   RETURN
C   MAKE TABLE OF S(J)=SIN(2*PI*J/NP), J=1,2,...,NT-1, NT=NP/4
440 NP=N
  MP=M
  NT=N/4
  MT=M-2
  IF (MT) 510,510,450
450 THETA=.7853981634
C   THETA=PI/2*(L+1)   FOR L=1
  JSTEP = NT
C   JSTEP = 2*(MT-L+1) FOR L=1
  JDIF = NT/2
C   JDIF = 2*(MT-L)   FOR L=1
  S(JDIF) = SIN(THETA)
  IF (MT-2) 510,470,470
470 DO 500 L=2,MT
  THETA = THETA/2.
  JSTEP2 = JSTEP
  JSTEP = JDIF
  JDIF = JDIF/2
  S(JDIF)=SIN(THETA)
  JC1=NT-JDIF
  S(JC1)=COS(THETA)
  JLAST=NT-JSTEP2
  IF (JLAST-JSTEP) 500,480,480
480 DO 490 J=JSTEP,JLAST,JSTEP
  JC=NT-J
  JD=J+JDIF
490 S(JD)=S(J)*S(JC1)+S(JDIF)*S(JC)
500 CONTINUE
510 IF (IFS) 70,30,70
  END

```



SUBROUTINE PSDSAVE (NSECS)

SUBROUTINE PSDSAVE OUTPUTS TO TAPE THE RAW PSD (1024 REAL NUMBERS PER  
2048 WORD RECORD) PRECEDED BY A 25 WORD HEADER RECORD.

CALLING ARGUMENT: NSECS = LENGTH OF TIME SERIES DATA SLICE

DIMENSION IARRAY(2048)

COMMON DATUMS(4500), IHEAD(20), CAL(2,4), IVARSW(10)

COMMON/IGD/ MUSCLE(7,4), IFORR(2)

EQUIVALENCE (IARRAY,DATUMS(2100))

MTD=2

IEOF=60000K

IRITE=50000K

IBFILE=40000K

CREATE AN OUTPUT HEADER RECORD

DO 10 J=1,25

IARRAY(J)=0

10 CONTINUE

DO 20 J=1,20

IARRAY(J)=IHEAD(J)

20 CONTINUE

NRECS=NSECS/2

IF(NRECS.EQ.0) NRECS=1

IARRAY(21)=NRECS

IWRD=IRITE+25

CALL MTDIO (MTD,IWRD,IARRAY,ISTAT,IER)

OUTPUT FOLLOWING DATA RECORDS CONTAINING THE RAW PSD VALUES

K=0

DO 50 J=1,NRECS

DO 40 JJ=1,1024

DATUMS(2099+JJ)=DATUMS(K+JJ)

40 CONTINUE

IWRD=IRITE+2048

CALL MTDIO (MTD,IWRD,IARRAY,ISTAT,IER)

K=K+1024

50 CONTINUE

WRITE TWO EOFs AND BACK UP TO JUST BEFORE THEM

DO 60 J=1,2

60 CALL MTDIO (MTD,IEOF,IARRAY,ISTAT,IER)

IBACK=IBFILE+1

DO 70 J=1,2

70 CALL MTDIO (MTD,IBACK,IARRAY,ISTAT,IER)

RETURN

END



OVERLAY OV7  
SUBROUTINE POUT(DTYPE,START,NSECS,SPAN,VOLTSEC,PSDSUM,STDERR,WHOA)

SUBROUTINE POUT OUTPUTS TO THE LINE PRINTER HEADER AND CALIBRATION DATA  
AVERAGED FORCE DATA, AND SMOOTHED PSD SPECTRUM.

CALLING ARGUMENTS: DTYPE - TYPE OF DATA

START - BEGIN TIME OF PRESENT DATA SLICE

NSECS - LENGTH OF DATA SLICE

SPAN - FOR FORCE DATA NUMBER OF SECS FORCE DATA  
AVERAGED OVER; FOR PSD DATA FREQUENCY  
INTERVAL FOR SMOOTHING ESTIMATES

VOLTSEC - INTEGRATED EMG VALUE (TIME DOMAIN)

PSDSUM - INTEGRATED EMG VALUE (FREQUENCY DOMAIN)

STDERR - NORMALIZED STANDARD ERROR OF ESTIMATE FOR  
PSD VALUES

WHOA - FLAG SET WHEN DATA PROCESSING ATTEMPTED  
WITHOUT CALIBRATION

INTEGER DTYPE,WHOA

DIMENSION IEDATE(3),IDDATE(3),CUMPCT(410)

DIMENSION AVGLBS(130),PSD(410),ITIME(3),PSDN(410),PERCNT(410)

COMMON DATUMS(4500),IHEAD(20),CAL(2,4),IVARSW(10)

COMMON/IGD/ MUSCLE(7,4),IFORR(2)

EQUIVALENCE (ISUBNO,IHEAD(1)),(IEDATE,IHEAD(2))

EQUIVALENCE (IDDATE,IHEAD(5)),(IFLITE,IHEAD(8))

EQUIVALENCE (IRUN,IHEAD(9)),(MUS,IHEAD(20))

EQUIVALENCE (IATAPE,IHEAD(10)),(ISAMPE,IHEAD(11))

EQUIVALENCE (ISAMPM,IHEAD(12)),(AVGLBS,DATUMS(4000))

EQUIVALENCE (PSD,DATUMS(1)),(PSDN,DATUMS(4000))

EQUIVALENCE (PERCNT,DATUMS(3000)),(CUMPCT,DATUMS(3500))

EQUIVALENCE (EXPVAL,DATUMS(2500)),(STDEV,DATUMS(2501))

EQUIVALENCE (EMGVAR,DATUMS(4500)),(ITIME,IHEAD(16))

DATA MUSCLE /'BR','AC','HI','AL','B','IC','EP','B',' ','RA','DI  
'AL','IS',' ','G','AS','TR','OC','NE','MI','US','S','OL','EU'  
'S',' ',' ',' '

DATA IFORR /'F','R+'/'

LPT=12

IF(WHOA.EQ.1) GO TO 150

BRANCH TO THE PROPER PRINT OUT SECTION FOR THE DATA TYPE

IF(DTYPE=4) 5,50,100

PRINT OUT HEADER INFORMATION

5 IP=1

IF (IFLITE.GE.0) IP=2

WRITE(LPT,10)

10 FORMAT(1H1,/,48X,'CARDIOVASCULAR LABORATORY',/,47X,'EMG DATA PRO  
CESSING PROGRAM',/,47X,'-----')

WRITE(LPT,20)ISUBNO,IEDATE,IFORR(IP),IFLITE,IRUN,IATAPE,IDDATE,(MU  
SCLE(J,MUS),J=1,7),ISAMPE,ISAMPM

20 FORMAT(////,47X,'\*\* HEADER INFORMATION\*\*',/,/, ' SUBJECT NO.:', 7  
@X,13,18X,'EXPERIMENT DATE:',6X,12,/,/,12,/,/,12,10X,'FLIGHT REFERENCE

```

@NCE DAY:',A3,I3,/,/, ' RUN NO.:',I1X,I2,I9X,'ANALOG TAPE NO.:',7X,I
@3,I4X,'DIGITIZING DATE:',6X,I2,/,/,I2,/,/,I2,/,/, ' MUSCLE:',9X,I7A2
@,/,/, ' EMG SAMPLE RATE (SAMP/SEC):',I5,8X,'FORCE SAMPLE RATE (SAMP/
@SEC):',I5)

```

```

WRITE(LPT,30) (CAL(1,J),J=1,4)
30 FORMAT(////,49X,'** CALIBRATION DATA **',/,/, ' FORCE CAL',/, ' ----
@-----',/,/, ' CAL WEIGHT (LBS):',F7.1,I6X,'COUNTS/LB:',2X,F6.0,5X,'A
@VCG. BASELINE COUNT:',2X,F6.0,5X,'CAL GAIN:',2X,F3.0)
WRITE(LPT,40) (CAL(2,J),J=1,2)
40 FORMAT(///, ' EMG CAL',/, ' -----',/,/, ' RMS AMPLITUDE (MICROVOLTS)
@:',F6.0,7X,' COUNTS/MICROVOLT:',F5.2,/,/)
GO TO 200

```

# PRINT OUT FORCE DATA

```

50 WRITE(LPT,60) ITIME,NSECS,SPAN
60 FORMAT(1H1,51X,'*** FORCE DATA ***'///5X,'START TIME-',I3,':',I2,'
@:',I2,I5X,'DATA LENGTH (SECS)-',I4,I5X,'AVERAGE FORCE INTERVAL (SE
@CS)-',F3.0,/,/)
LOOPS=NSECS/IFIX(SPAN)
IF (LOOPS.GT.49) GO TO 66
WRITE(LPT,62)
62 FORMAT(47X,'INTERVAL',5X,'AVERAGE FORCE (LBS)'/,47X,'-----',5X,
@,'-----',/,/)
DO 65 J=1,LOOPS
WRITE(LPT,64) J,AVGLBS(J)
64 FORMAT(49X,I3,I5X,F5.1)
65 CONTINUE
GO TO 200
66 IF ((LOOPS/2)*2.NE.LOOPS) LOOPS=LOOPS-1
WRITE(LPT,68)
68 FORMAT(15X,'INTERVAL',5X,'AVERAGE FORCE (LBS)',20X,'INTERVAL',5X,'
@AVERAGE FORCE (LBS)'/,15X,'-----',5X,'-----',20X
@,'-----',5X,'-----',/,/)
LOOPS2=LOOPS/2
DO 75 J=1,LOOPS2
JJ=J+LOOPS2
WRITE(LPT,73) J,AVGLBS(J),JJ,AVGLBS(JJ)
73 FORMAT(17X,I3,I5X,F5.1,29X,I3,I5X,F5.1)
75 CONTINUE
GO TO 200

```

# PRINT OUT PSD DATA

```

100 WRITE(LPT,110) ITIME,NSECS,VOLTSEC,SPAN,STDERR,PSDSUM,EXPVAL,STDEV
@,EMGVAR
110 FORMAT(1H1,39X,'*** POWER SPECTRAL DENSITY OF EMG DATA ***'///4X,'
@START TIME-',I3,':',I2,':',I2,I1X,'DATA LENGTH (SECS)-',I4,I4X,'IN
@TEGRATED EMG (MICROVOLT*SEC)-',E12.4//4X,'BANDWIDTH (HZ)-',F6.3,I10
@X,'NORMALIZED STANDARD ERROR-',F6.3,5X,'INTEGRATED PSD (MICROVOLT*
@*2)-',E13.4,/,/,4X,'MEAN (HZ)-',F6.1,I5X,'STANDARD DEVIATION
@ (HZ)-',F6.1,8X,'EMG VARIANCE (MICROVOLT**2)-',E15.4,/)
LOOPS=IFIX(400./SPAN)+1
FREQ=SPAN/2.
IF (LOOPS.GT.46) GO TO 130
WRITE(LPT,115)
115 FORMAT(30X,' FREQ',13X,'PSD',11X,'PSD',10X,'% OF',10X,'CUM %'/,30X
@,' (HZ)',8X,'(MMV**2/HZ)',8X,'NORM',9X,'TOTAL',9X,'TOTAL',/,31X,'
@-----',8X,'-----',8X,'-----',9X,'-----',9X,'-----',/,)
DO 125 J=1,LOOPS
WRITE(LPT,120) FREQ,PSD(J),PSDN(J),PERCNT(J),CUMPC(J)
120 FORMAT(30X,F6.2,7X,F10.3,8X,F6.4,7X,F6.2,8X,F6.2)
FREQ=FREQ+SPAN

```



```

125 CONTINUE
GO TO 200
130 IF((LOOPS/2)*2.NE.LOOPS) LOOPS=LOOPS-1
FREQ2=FLOAT(LOOPS/2)*SPAN+FREQ
WRITE(LPT,135)
135 FORMAT(' FREQ', 8X,'PSD', 9X,'PSD',8X,'% OF',6X,'CUM %',22X,'FREQ
@', 9X,'PSD' , 8X,'PSD',8X,'% OF',6X,'CUM %'/2X,'(HZ)',4X,'(MMV**2/
@HZ)',5X,'NORM',7X,'TOTAL',5X,'TOTAL',22X,'(HZ)',5X,'(MMV**2/HZ)',5
@X,'NORM',6X,'TOTAL',6X,'TOTAL',/, ' ----',4X,'-----',5X,'---
@-',7X,'-----',5X,'-----',22X,'----',5X,'-----',5X,'----',6X,
@'-----',6X,'-----',/)
LOOPS2=LOOPS/2
DO 145 J=1,LOOPS2
JJ=J+LOOPS2
WRITE(LPT,140) FREQ,PSD(J),PSDN(J),PERCNT(J),CUMPC(J),FREQ2,PSD(J
@J),PSDN(JJ),PERCNT(JJ),CUMPC(JJ)
140 FORMAT(F7.2,3X,F10.3,5X,F6.4,5X,F6.2,5X,F6.2,20X,F6.2,4X,F10.3,5X,
@F6.4,5X,F6.2,5X,F6.2)
FREQ=FREQ+SPAN
FREQ2=FREQ2+SPAN
145 CONTINUE
GO TO 200

C
C
C
150 WRITE(LPT,160)
160 FORMAT(' ERROR...DATA ANALYSIS ATTEMPTED WITHOUT CALIBRATION FIRST
@. CHECK DATA CARD ORDER.')
```

```

STOP
200 RETURN
END
```



.MAIN

001620  
000,000 POUT 003477  
000,001 T1NPU  
FIND  
OUT60 002304  
000,002 DCAL 001035  
000,003 GRAPH 000165  
000,004 EMG  
FORCE 001167  
000,005 FFIPS  
WINDO  
RFFT  
FFT  
PSDSA 007540  
000,006 CINPU 001217  
011620

I  
FREAD  
THREA  
RDFLD  
READL  
OPEN  
FOVLD  
FOPEN  
DFILW  
CLOSE  
FREDI  
FALOC  
ARYSZ  
FSBR  
CGT  
IABS  
SMPY  
SDVD  
IFIX  
ABSLT  
RLSE  
INIT  
MTDIO  
IPWER  
RIPWR  
CUS  
SQRT  
PLY1  
BREAK  
FLIP  
ARGUM  
FRGLD  
FARGO  
FL  
STREG  
LDREG  
MVBT  
LDO  
STOP  
FINIT  
FLINK  
RTER  
WRCH  
BDASC  
BASC

COUT  
LDSTB  
MOVEF  
CPYAR  
NAD  
FPZER  
FPTRS  
DUMMY  
ARDUM  
HMPYD  
TMIN

NMAX 023223  
ZMAX 000231  
CSZE 021526  
EST 000000  
SST 000000

.FREA 000050  
.FWRI 000051  
.BRD 000052  
.BWR 000053  
.ALLO 000074  
.THRE 000075  
.RDFL 000076  
.RDFC 000077  
.WRIT 000100  
.READ 000101  
.WRTS 000102  
.REDS 000103  
.FOPE 000104  
.FRED 000105  
.FALU 000106  
.ARYS 000107  
.FSUB 000110  
.FSBR 000111  
.CGT 000112  
IA.S 000113  
.SMPY 000114  
.SDVD 000115  
XI.X 000116  
IF.X 000116  
.ABS 000117  
.IPWR 000120  
.FLIP 000127  
.FARG 000131  
.FRGL 000132  
.FRGO 000133  
.FRG1 000134  
SN.L 000135  
DB.E 000135  
FL.AT 000143  
.MVBC 000160  
.MVBT 000161  
.LD0 000162  
.LD1 000163  
.LD2 000164  
.ST0 000165  
.ST1 000166  
.ST2 000167  
.STOP 000170  
.PAUS 000171



.FINI 000172  
.FCAL 000173  
.FSAV 000175  
.FRET 000176  
.RTER 000200  
.RTE0 000201  
.RTES 000202  
.WRCH 000203  
.COUT 000204  
.CIN 000205  
.LDBT 000206  
.STBT 000207  
.MOVE 000210  
.CPYA 000211  
.CPYL 000212  
.MAD 000213  
.MADU 000214  
.IOCA 000215  
SUCOM 000216  
.NDSP 000217  
TVR 000217  
AFSE 000220  
SP 000221  
.OVFL 000222  
.SV0 000223  
QSP 000224  
NSP 000225  
FLSP 000225  
USTAD 000400  
C IGD 000444 000036  
.MAIN 000602  
EMG 001710  
CINPU 001716  
DCAL 001727  
TINPU 001730  
GRAPH 001753  
FFTPS 001753  
POUT 002016  
FSAV 002175  
FRET 002176  
FQRET 002177  
FORCE 002462  
FIND 003147  
WINDO 003225  
RFFT 003450  
OUT60 003612  
FFT 005333  
XAS. 006117  
DABS. 006117  
ABS. 006117  
FIPR1 006121  
COS. 006122  
SIN. 006123  
SQRT. 006124  
FPLY1 006125  
FBRK1 006126  
FLIP2 006130  
FLIP1 006130  
FFLD1 006135  
FFST1 006136  
FAD1 006137  
FSB1 006140



2

FDV1 006142  
FXFL1 006143  
FLFX1 006144  
FSGN1 006145  
FNEG1 006146  
FCLE1 006147  
FCLT1 006150  
FCGE1 006151  
FCGT1 006152  
FCEQ1 006153  
FRST1 006154  
FRST2 006155  
FRLD2 006156  
FRLO1 006157  
FCALL 006173  
FRCAL 006174  
MPY0 006226  
MPY 006227  
DVD 006230  
PSDSA 011007  
.I 011656  
IOPTR 012024  
FERT0 012052  
FERT1 012055  
FERTN 012061  
BRD 012100  
BWR 012104  
FREAD 012111  
FWRIT 012115  
ALLOC 015761  
THREA 016010  
RDFLD 016025  
RDFCH 016031  
READL 016142  
WRITL 016152  
REDS 016366  
WRITS 016377  
MTOPD 016460  
APPEN 016463  
OVOPN 016466  
OPEN 016471  
OVLOD 016740  
FOVLO 016740  
FOPEN 017007  
DFILW 017205  
DELET 017205  
CLOSE 017237  
FCLOS 017237  
FREDI 017307  
FALOC 017420  
ARYSZ 017454  
FSUBA 017472  
FSBR 017635  
CGT 017671  
.IABS 017714  
SMPY 017723  
SDVD 017747  
.IFIX 020017  
ABS 020030  
RLSE 020035  
INIT 020053  
MTDIO 020072  
IPWK 020146

RIPWR 020223  
CS 020273  
SN 020277  
SWR 020440  
PLY1 020600  
BRK 020633  
FLP0 020650  
FLP 020653  
FARGU 020676  
FRGLD 020732  
FRG1 020742  
FRG0 020747  
FL 020767  
FS 021026  
FB 021107  
FA 021110  
FM 021240  
FXL 021377  
FLX 021435  
FSG 021502  
FNG 021523  
FLE 021557  
FLT 021561  
FGE 021563  
FGT 021565  
FEQ 021567  
FD 021637  
ST1 021742  
ST2 021744  
LDR1 021767  
LDR2 021771  
MVB1 022023  
MVB0 022027  
LD0 022060  
LD1 022062  
LD2 022064  
EXIT 022100  
STOP 022105  
PAUSE 022112  
FINIT 022163  
SAV0 022205  
SAV2 022255  
SAV3 022262  
RSTR 022275  
QRSTR 022311  
.QFLO 022313  
RTER 022355  
RTESP 022357  
RTE0 022371  
WRCH 022571  
.BDAS 022606  
.BASC 022670  
COUT 022725  
CIN 022737  
LDB 022746  
ST6 022760  
MOVEP 022776  
CPYAR 023023  
CPYLS 023054  
MAD 023063  
MADO 023064  
.FLSP 023110  
INT0 023127

FMMA	177777
FRISK	177777
.FLSL	177777
QICK	177777
OV7	000,000
OV2	000,001
OV3	000,002
OV4	000,003
OV5	000,004
OV6	000,005
OV1	000,006



APPENDIX H  
CF-16/A  
GENERAL A/D PROGRAM  
DOCUMENTATION

Author: Winston Blackmon  
Engineering Systems Branch  
Institutional Data Systems Division  
Data Systems and Analysis Directorate

#### 4.2.0 LBNP ONLINE WITH SIGMA 3 REQUIREMENT DELETED

##### 4.3.0.0 General Data Acquisition System.

The DBAS (General Data Acquisition System) acquires, digitizes, and records test data from an analog playback. The user is allowed up to 3 different rates and sixteen channels of input. The only restriction are that the rates must be integer multiples of each other and that the A/D channels be grouped in the order of highest rate, mid-rate, and slow rate. The GDAS is loaded from the nine-track tape and will operate independently of the Sigma 3 Biomedical System. It will, however, require the FM recorder and Systron Donner used by the Sigma 3 unless the user can supply his own recorder and Systron Donner. The GDAS data storage and output, operator interaction, and subroutines are described in the following sections.

##### 4.3.1.0 Data Storage

There are two types of data storage - that which is used by the executive for control purposes; and that which is used as buffers for digitized data accumulation. The data used by the executive for control is stored in locations  $80_{16}$  to  $A5_{16}$  of scratch pad (locations 0-255<sub>10</sub>). The following table defines these data.

LOCATION	VARIABLE	DESCRIPTION
$80_{16}$ ( $128_{10}$ )	CLCNT	User input number of microseconds to elapse between each A/D read of the highest speed data.
$81_{16}$ ( $129_{10}$ )	NDRATE	User input number of different rates (1, 2, or 3)
$82_{16}$ ( $130_{10}$ )	NØLØ	Minus total number of channels (computed from user input)
$83_{16}$ ( $131_{10}$ )	LØCNT	Number of reads of the high speed channels before the lowest speed channel is read.
$84_{16}$ ( $132_{10}$ )	NØMED	User input minus total number of medium and high speed channels.
$85_{16}$ ( $133_{10}$ )	MEDCNT	Number of reads of the high speed channels before the medium speed channels are read (user input)
$86_{16}$ ( $134_{10}$ )	NØHS	Number of high speed channels to read - user input

LOCATION	VARIABLE	DESCRIPTION
87 <sub>16</sub> (135 <sub>10</sub> )	BUFZER	Buffer zero starting address minus 1. Fixed at 77C <sub>16</sub>
88 <sub>16</sub> (136 <sub>10</sub> )	BUFZND	Buffer zero end address. Set to B64 <sub>16</sub> , but will be changed if the user sets the data buffer length to less than 1000.
89 <sub>16</sub> (137 <sub>10</sub> )	BUFØNE	Buffer one starting address minus 1. Fixed at B7D <sub>16</sub> .
8A <sub>16</sub> (138 <sub>10</sub> )	NDBUF1	End address for Buffer one. Set at F65 <sub>16</sub> , but will be changed if the user selects a data buffer length less than 1000.
8B <sub>16</sub> (139 <sub>10</sub> )	BUFUSE	Buffer in use pointer. Initially set at zero but is changed when data buffer 1 is being filled with digitized data.
8C <sub>16</sub> (140 <sub>10</sub> )	MAXBUF	Maximum buffer length including header information. Fixed at 401 <sub>16</sub> (1025 <sub>10</sub> )
8D <sub>16</sub> (141 <sub>10</sub> )	TPØUT	Tape output flag. Initially set to zero, but is set to 1 when a data buf- fer has been filled and is ready for output to tape.
8E <sub>16</sub> (142 <sub>10</sub> )	DØNCLK	Flag, initially zero, set to one when- ever the Systron Donner clock timer has been read and are ready to be unpacked.
8F <sub>16</sub> (143 <sub>10</sub> )	BGNCHN	Address of the first high speed A/D convertor channel. Default value is zero.
90 <sub>16</sub> (144 <sub>10</sub> )	ADCRD	Flag which indicates an A/D convertor read is in progress (= 1)
91 <sub>16</sub> (145 <sub>10</sub> )	DØNMS	} Timer read from the Systron Donner when data buffers are switched. These times are in a pseudo BCD format, Complemented.
92 <sub>16</sub> (146 <sub>10</sub> )	DNSEL	
93 <sub>16</sub> (147 <sub>10</sub> )	DNMIN	
94 <sub>16</sub> (148 <sub>10</sub> )	DNØUR	
95 <sub>16</sub> (149 <sub>10</sub> )	STHR	} Systron Donner timer at which the user wants to begin digitizing data. Default values are zero.
96 <sub>16</sub> (150 <sub>10</sub> )	STMIN	
97 <sub>16</sub> (151 <sub>10</sub> )	STSEL	
98 <sub>16</sub> (152 <sub>10</sub> )	STRMS	



LOCATION	VARIABLE	DESCRIPTION
99 <sub>16</sub> (153 <sub>10</sub> )	} STPHR STPMIN STPSEL STØPMS	Syston Donner timer at which the user wants to stop digitizing data. Default values are 1E00 <sub>16</sub> .
9A <sub>16</sub> (154 <sub>10</sub> )		
9B <sub>16</sub> (155 <sub>10</sub> )		
9C <sub>16</sub> (156 <sub>10</sub> )		
9D <sub>16</sub> (157 <sub>10</sub> )	BUFHD1	Buffer one header address minus one. Used for tape output.
9E <sub>16</sub> (158 <sub>10</sub> )	BUFHDO	Buffer zero header address minus one. Used for tape output.
9D <sub>16</sub> (159 <sub>10</sub> )	MBUFLG	Negative of Maximum buffer length including header. Set by user with a default value of -1025 <sub>10</sub> .
A0 <sub>16</sub> (160 <sub>10</sub> )	MXDBUF	Maximum data buffer length. Set at 1000 <sub>10</sub> .
A1 <sub>16</sub> (161 <sub>10</sub> )	PAR	Number of parity errors detected during tape output of the digitized data.
A2 <sub>16</sub> (162 <sub>10</sub> )	TIMG	Number of timing errors detected during tape output of the digitized data.
A3 <sub>16</sub> (163 <sub>10</sub> )	HDBFLG	Length of the header portion of each data buffer. Set at 25 <sub>10</sub> .
A4 <sub>16</sub> (164 <sub>10</sub> )	BUFTTY	Address of the teletype data buffer. Address is 100 <sub>16</sub> (256 <sub>10</sub> ).
A5 <sub>16</sub> (165 <sub>10</sub> )	STATUS	Temporary storage for the status word of the tape write operation.
A6 <sub>16</sub> (166 <sub>10</sub> )	TAPEND	Flag signifying end of operation on tape. Writes during test.

Also used for control purposes, are the addresses of subroutines which are general in nature or called by more than one other subroutine. These are defined below

70 <sub>16</sub> (112 <sub>10</sub> )	FETCH	Used to coordinate data input, e.g., reads a buffer from teletype and checks for data errors.
71 <sub>16</sub> (113 <sub>10</sub> )	CRLF	Outputs a carriage return-line feed to the teletype.
72 <sub>16</sub> (114 <sub>10</sub> )	ØTL	Outputs a line to the teletype.

73 <sub>16</sub> (115 <sub>10</sub> )	UNPACK	Unpacks the System Donner Time from pseudo BCD to binary.
74 <sub>16</sub> (116 <sub>10</sub> )	CØNVR	Converts ASCII Input number to a signed binary integer and checks for overflow (number too large).
75 <sub>16</sub> (117 <sub>10</sub> )	TDUMP	Starts tape I/O during a test and writes EOF Marks after test.
76 <sub>16</sub> (118 <sub>10</sub> )	ØDEL	Converts a binary integer to a string of ASCII digits for output.
77 <sub>16</sub> (119 <sub>10</sub> )	ITT	Inputs one ASCII Character from the teletype.
78 <sub>16</sub> (120 <sub>10</sub> )	ØTT	Outputs one ASCII character to the teletype.
79 <sub>16</sub> (121 <sub>10</sub> )	TYPIN	Routine called to fill the teletype input buffer.
7A <sub>16</sub> (122 <sub>10</sub> )	RESTØR	Restores the executive data table (locations 80 <sub>16</sub> to A5 <sub>16</sub> ) to it's initial value, and zero the digitized data buffer areas.
7B <sub>16</sub> (123 <sub>10</sub> )	STØPIT	Routine to stop the data acquisition and restore the acquisition routine to its initial state.

The buffers which are used to accumulate the digitized data consist of two blocks of core, each 1025<sub>10</sub> words in length. The first 25<sub>10</sub> words of each block are leader data, 20<sub>10</sub> of which may be set by the user. The blocks are described below.

LOCATION	BUFFER WORD	PURPOSE
764 <sub>16</sub> } B65 <sub>16</sub>	1	Buffer record counter incremental by the executive.
765-768 <sub>16</sub> } B66-B69 <sub>16</sub>	2 - 5	Storage for Systron Donner time at which the first data for this record were acquired.
769-77C <sub>16</sub> } B6A-B7D <sub>16</sub>	6-25	Header information stored by the user. The 20 <sub>10</sub> words are constant during a given acquisition.

77D-B64<sub>16</sub> }  
B7E-F65<sub>16</sub> }

26  
1025

Digitized data. Note that these end addresses are only valid for the maximum data buffer length, 1000<sub>10</sub>. Also the data will be stored in the order that it is read, i.e. there is no attempt to group the data by channels, rates, etc.

#### 4.3.2.0

#### Operator Interaction with the DGAS

The GDAS acts in a conversational manner with the user via the teletype. In response to questions and for all data (numeric) input, a slash (/) is used to terminate the input. For correction of input prior to the slash, the last character input can be deleted by an up arrow (↑) and the entire line deleted by a back (←). Of the particular input is rejected by the GDAS, the back arrow and a carriage return - line feed will be printed. The quantity will be requested again.

The following inputs are requested by the GDAS.

6 DAS.

No. Rates  $\left\{ \begin{array}{l} 1 \\ 2 \\ 3 \end{array} \right\}$

Number of different data rates, the user must input, 1, 2, or 3.

MICROSEC PER HI SPD READ.

Number of microseconds to elapse between reads of the highest data rate. For example, 1000 samples/second would be 1000 microseconds; 100 samples/second would be 10000 microseconds. The lowest rate which can be input for a high speed rate is 62500 microseconds/read or 16 samples/second. This does not preclude lower rates if more than one rate is used.

DATA BUF LENGTH, MAX I S 1000

Length of the data buffer to be filled before output to tape. This number should be chosen with care-especially with multiple data rates-so that the data from each rate will appear at the same place in each buffer. For example, with four channels of data at 320 samples/second and eight channels at 2 samples/second, a good number for buffer length would be 648. (4x160+8). The buffer would be dumped each  $\frac{1}{2}$  second.

ADC CHANNEL

Address of the first A/D convertor channel to be read, i.e. the first high speed channel.



GAIN	Gain to be applied to each A/D conversion. Gains available are times 1, 2, 4, or 8.
NO. HI SPD CHANS.	Number of high speed channels to read. For one rate, this number will be the total number of channels.
NO. MED SPED CHANS.	Number of medium speed channels to read. The output applies only when there are two or three rates. When there are two rates, this is the number of lower speed channels to read.
NO. LO SPED CHANS.	Number of low speed channels. This output applies only when there are three rates. This input is the number of channels for the lowest rate.
NO. HI SPD READS/MEDIUM READ.	Input the number of times to read the fastest rate before reading the mid-rate (or in the case of two rates, the low rate). e.g., Suppose the high rate is 100 samples/second and the mid rate is 25 samples/second. The number to be input would be $4 = 100/25$ .
NO. HISPD READS/LOW READ.	In the case of three data rates, how many times are the highest speed channels read before the lowest speed channels are read? Eg. Suppose the rates are 500, 250, and 25 samples per second, respectively. The response to this input would be $500/25 = 20$ .
INPUT DATA FOR TAPE HEADER?	Yes or no response. This question asks if the user desires to store constant numeric data in words 6 through 25 in the header of each buffer. This data is output to tape with the digitized data. See 4.3.1.0 for further data description. If a yes response is given, the GDAS asks for the input to word 6, then 7, then eight, etc. If the user desires to terminate the header input, he may enter letter A followed by a slash.
WRITE BEGINNING EOF?	This allows the user to write a end of file on the beginning of his tape. One reason to do so is to allow the tape to space past the extra 3 inches at the beginning. This 3 inch gap is especially important during high speed acquisition rates. (See 4.3.2.1 for more detail.)

## START TIMES?

Does the user want to input the Systron Donner times (hr, min, Sec) at which the GDAS is to start digitizing data? Yes or No response required! If yes, the user must input the times as requested. If no, the GDAS uses zero as the default value.

## STOP TIMES?

Same as start times except these times are when to stop acquiring data. The default value is 1E00<sub>16</sub> if a no response is entered.

A  $\longrightarrow$

## S TO START

This is an input which prevents the GDAS from looking at the start times. It is useful if the default is entered or if the user must position an analog tape for playback. The quantity S/ will start the GDAS.

When a test is completed, either by reaching the stop times or by an x input, the GDAS will allow additional tests to be conducted with the same data input. This question is

## CONTINUE TEST WITH NEW TIMES?

If a yes response is entered, then new start/stop times are requested. Then the test may be restarted by using the s/input. If a no response is entered, the GDAS will ask of the tape is to be rewound between tests. If yes, it rewinds the tapes. In either case, the GDAS initializes itself, i.e., the user must enter new inputs for rates, etc.

### 4.3.2.1. Timing Estimates for GDAS

The GDAS will digitize data at various rates, write various length records depending on the user's inputs. For this reason, the user should put some "handle" on whether his specific rates and data buffer lengths are feasible. As mentioned previously, the minimum single rate speed is 16 samples/second. To determine the maximum rate, several tests were conducted. The results of these tests indicate that the tape speed will be the limiting factor. The only control the user has over tape speed is to vary his buffer length. A longer buffer will take longer to fill than to dump in most cases. The limiting case is:

$$\begin{array}{ll} \text{Time to fill buffer} & & 0.046 + 0.0001 \text{ (buffer} \\ \text{with data in seconds} & > & \text{length in words).} \end{array}$$

The above inequality has been tested with 6 channels of data at 1000 samples/second using various buffer lengths. The minimum buffer length which worked was 690. The time required to fill the buffer and to dump the buffer in this case is identical. Note that the above inequality and test are based on the GDAS not being required to write a record from the BOT (beginning of tape) position. The first operation from BOT requires .197 seconds to simply start and stop the tape. Hence, if the user is dealing with high rates, he must not start the tape from BOT. A beginning dummy record or an end of file can be used to position past the BOT.



#### 4.3.3.0.0 Subroutine Description

All routines are documented according to the pattern shown below:

.0	name
.1	purpose
.2	calling sequence
.3	software/software interfaces
.4	input data
.5	output data
.6	storage required
.7	description
.8	flow diagrams

#### 4.3.3.1.0 BEGIN - Pre- and Post Test Control Routine

4.3.3.1.1 The purpose of this program is to control the flow of input data from the user and to insure that all required inputs are within tolerances. BEGIN also starts the test by calling the test routine.

4.3.3.1.2 Calling Sequence -  
JST \*BEGIN  
Where BEGIN is a core location which contains the value  $110_{16}$ .

4.3.3.1.3 Software/Software Interfaces  
BEGIN uses both external and internal subroutines. If no indication is given, the subroutine should be assumed to be external.

RESTOR Subroutine called to restores all the GDAS variables that may be changed by the user to their default values. It also resets all the I/O devices.

CRLF Subroutine called to issue a carriage return - line feed to the teletype.

ØTL Subroutine called to output a message to the teletype.

FETCH Subroutine called to output a message and return the numeric response in the A register. It skips forward one instruction if the numeric input was free from character errors.

TYPIN Subroutine called to fill the 16 character teletype input buffer, BUFR.

CHKIT Internal subroutine used to check for the no change condition in variables whose input is optional. The routine FETCH returns an error condition for this particular case when possibly no input error has been made.

EINS -	Subroutine called to input data pertinent to a one rate data acquisition.
ZWEI -	Subroutine called to input data pertinent to a two rate data acquisition.
DREI -	Subroutine called to input data pertinent to a three rate data acquisition
HEADER	Subroutine called to input and store constant data into the 25 word tape buffer header. Only twenty words of input are allowed: The other five are reserved for record number and time.
TPCHK	Subroutine called to check for the tape ready condition, i.e. tape is on-line and not write protected.
ADCHK -	Subroutine called to check that the A/D convertor is not in manual mode.
READON -	Internal subroutine called to read the Systron Donner time code generator to determine when to start the test.
UNPACK -	Subroutine called to unpack the Systron Donner time from pseudo BCD to binary.
TEST -	Subroutine called to start the data acquisition and remain in control throughout the test.
4.3.3.1.4	Input Data - The following parameters must be available for storage/use by BEGIN
BUFR	A 16 word teletype input buffer filled by TYPIN.
NORATE	Number of A/D rates
CLCNT	Number of microseconds between highest rate reads of the A/D convertor
MAXBUF	Maximum data buffer length (default 1000 <sub>10</sub> )
BUFZER	Starting address minus 1 for data buffer zero. (Set at 77C <sub>16</sub> )
BUFZND	End address of buffer zero (default B64 <sub>16</sub> )
BUFONE	Start address minus 1 for Buffer One (set to B7D <sub>16</sub> )
NDBUF1	End address of buffer one (default F65 <sub>16</sub> )

HDBFLG	Header length on data buffer (set at 25)
MBUFLG	Minus buffer length including header (default=-1025 <sub>10</sub> )
BGNCHN	First A/D channel to read (default = 0)
STHR	} Systron Donner Strat hour, minute, second, and Millisecond (default is 0)
STMIN	
STSEC	
STRMS	
STPHR	} Systron Donner stop hour, minute, second, and millisecond (default is 1E00 <sub>16</sub> )
STPMIN	
STPSEL	
STØPMS	

In addition to the above values, the pointers to subroutines FETCH, CRLF, RESTØR, UNPACK, and TYPIN must be available in Scratch pad (locations 0 to 255). BUFR is stored at 256<sub>10</sub> to 271<sub>16</sub>. The binary numbers for the above values are returned from subroutine FETCH.

The inputs which may be alphabetical are put in the teletype buffer, BUFR, and then checked for the appropriate alphabetical value by BEGIN.

The hours, minute, second, and millisecond are input by internal subroutine READØN.

#### 4.3.3.1.5 Output Data

The specific names in the input section may be changed if there is a default value. In addition to the above outputs, the following subroutines have calling arguments as indicated.

FETCH	-	The address of the message requesting numeric input is transfered through the x-register.
ØTL	-	The terminating character of a message is transferred through the S-register. The address pointing to the first word of the message is stored in the first location after the call to ØTL.
READØN	-	The code for the time to be returned (hours, minutes, seconds, milliseconds) is passed through the x register.
UNPACK	-	The pseudo BCD time code in passed to this subroutine through the A register.
CHKIT	-	The first character of the teletype buffer, BUFR, must be available. Since CHKIT is an internal subroutine within +255 of BUFR, no special action is required.



TYPIN\_ Any single character to be output before filling the teletype buffer with input must be passed through the A register.

4.3.3.1.6. Storage Required - BEGIN requires 183<sub>16</sub> or 387<sub>10</sub> words of memory.

4.3.3.1.7 Description

The tasks of BEGIN are to: obtain all the required and optional inputs from the user; insure that the tape and A/D converter are ready; and to start the acquisition at the specified Systron Donner time (If the default values are used, the acquisition starts immediately). BEGIN Accomplish. The first task by calling the appropriate support subroutines to get a message to the user requesting input and return a binary number for those numeric inputs. In the case of a question, a yes or no response is required. If an error is detected in either the optional or mandatory input, the request for data is typed again. In the specific case of the number of microseconds per high speed read, BEGIN cannot determine if the error was legitimate. It requests that the user verify that the input is okay. The only correct input which can trigger the verification is if the number of microseconds is greater than 32767 and less than 65535.

The second task is accomplished by calling the A/D convertor and tape check routines. The A/D convertor checks only for the convertor being in manual mode. If so, the routine outputs diagnostic message and continues checking until the convertor is placed in automatic mode. The tape check routine checks for tape on-line and not write protected. If either/or both conditions are not true, a diagnostic is output. The condition(s) must be connected before BEGIN will continue.

After the tape and A/D convertor are checked, BEGIN will then start the test sequence by asking about writing a beginning end of file. The user can write an end of file to position past the first three inches of tape (recommended for high data rates). BEGIN will then request an S to begin checking the Systron Donner times against the start times. This allows the user a chance to review his input data before actually starting the test.

After the test acquisition is completed, the user is asked about continuing the acquisition with new times. This allows the user to have a short calibration run and then a long test without changing his fixed data. The tape can be rewound if the user desires. The GDAS then initializes itself and starts the whole sequence again.

- 4.3.3.2.0 Subroutine FETCH
- 4.3.3.2.1 The purpose of this subroutine is to output a message requesting numeric input and convert the subsequent input to a binary integer.
- 4.3.3.2.2 The calling sequence is LDX Address of message  
JST FETCH  
Returns here if input in error  
Returns here if input is okay  
Binary integer is in the A register.
- 4.3.3.2.3 Software/Software Interfaces.  
FETCH Interfaces with the following subroutines -  
CRLF Subroutine which issues a carriage return - line feed to the teletype  
ØTL Subroutine which writes a line to the teletype  
TYPIN Subroutine which fills the teletype input buffer.  
CONVR Subroutine which converts ASCII integers to a binary number.
- 4.3.3.2.4 Input Data - The address of the message to be output is transferred to FETCH through the x register. A binary number is returned from CONVR in the A register, and if an error was detected, the overflow bit will be set. The address of the teletype buffer must be available in location A4<sub>16</sub> of scratchpad.
- 4.3.3.2.5 Output Data - Subroutine FETCH passes the terminating character of the message requesting input through the A register. The address of the first word of the message is stored in the first location after the call to the message writer, ØTL.  
  
FETCH transfers a blank character to the teletype input routine, TYPIN, through the A register.  
  
The address of the teletype buffer is transferred to the ASCII to binary conversion routine, CONVR, through the x register.  
  
The binary number that FETCH was requested to input is transferred to the calling program through the A register.
- 4.3.3.2.6 Storage Required - Subroutine FETCH requires D<sub>16</sub> or 13<sub>10</sub> locations.

## 4.3.3.2.7

## Description

When FETCH is called, the address of the message requesting an input number is in the x register. FETCH stores that address in the first location after the call to the message writer, ØTL. FETCH then issues a carriage return-line feed by calling the routine, CRLF. After returning from CRLF, FETCH loads the terminating character into the A register before calling, ØTL. The terminating character is always a period (.) for any message from FETCH. After the message is printed by ØTL, FETCH calls the routine TYPIN to fill the teletype input buffer. Upon return from TYPIN, FETCH loads the address of the teletype buffer into the x register and calls the conversion routine, CONVR. When CONVR returns control to FETCH, a binary number will be in the A register. The overflow bit will be set (i.e., 1) if CONVR was unable to convert the number. If the overflow is set, FETCH returns control to the first location after the calling program. If the overflow is reset (i.e., 0) the return is to the second location after the call.



4.3.3.3.0

Subroutine CLØCK and STOPIT

4.3.3.3.1

The purpose of this subroutine is to respond to the real time clock interrupt and to start the data acquisition. The purpose STOPIT is to reset the clock subroutine.

4.3.3.3.2

There is no calling sequence since this is an interrupt subroutine.

4.3.3.3.3

Software/Software Interfaces.

CLOCK calls internal subroutine.  
STRTAD to reset and start the A/D Converter.  
Internal subroutine STOPIT is called by external programs to reset the CLØCK routine variables to their initial state.

4.3.3.3.4

Input Data-There are no data inputs through calling sequences. However, the following variables must be available in scratchpad (locations 0-256<sub>10</sub>)

LOCATION	VARIABLE	DESCRIPTION
81 <sub>16</sub>	NØRATE	Number of Rates (1, 2, or 3)
82 <sub>16</sub>	NØLØ	Number of Low Speed A/D channels
83 <sub>16</sub>	LOCNT	Number of High speed reads per low speed read
84 <sub>16</sub>	NØMED	Number of medium speed channels.
85 <sub>16</sub>	MEDCNT	Number of high speed reads per medium speed read
86 <sub>16</sub>	NØHS	Number of high speed channels
87 <sub>16</sub>	BUFZER	Starting address for storing data into buffer zero.
88 <sub>16</sub>	BUFZND	Ending address for storing data into buffer zero.
89 <sub>16</sub>	BUFØNE	Address for buffer one.
8A <sub>16</sub>	NDBUF1	
8C <sub>16</sub>	MAXBUF	Data buffer length
8D <sub>16</sub>	TPØUT	Tape output flag

LOCATION	VARIABLE	DESCRIPTION
8E <sub>16</sub>	DONCLK	Systron Donner time conversion flag.
8F <sub>16</sub>	BGNCHN	Beginning A/D channel number
91 <sub>16</sub>	DØNMS	Temporary storage for Systron Donner times (Hours, minutes, seconds, and milliseconds)
92 <sub>16</sub>	DNSEC	
93 <sub>16</sub>	DNMIN	
94 <sub>16</sub>	DNHRS	
8B <sub>16</sub>	BUFUSE	Flag indicating which data buffer is currently being filled (buffer zero or one).

Subroutine STOPIT sets the following variables to their initial state.

BUFEND    End of Buffer sentinel set to 0.

DØNCLK    Systron Donner convert flag set to 0

CNTMED    Medium A/D rate internal counter to -1

CNTLO    Low A/D rate internal counter to -1

ADCNUM    Number of A/D channels to read to -1

4.3.3.3.5    Output Data - CLOCK outputs the start codes to the A/D convertor and sets up the automatic input locations for the digital input values.

4.3.3.3.6    Storage Required - The CLOCK program requires 66<sub>16</sub> or 102<sub>10</sub> locations. This includes 14<sub>10</sub> locations required by STØPIT.

4.3.3.3.7    Description - When a clock interrupt occurs, the CLOCK subroutine responds by clearing the interrupt. Since the time period between high speed A/D reads was used as a clock count, each interrupt signals the start time for high speed data acquisition. The low and medium speed counters are checked (and incremented) to determine if the low and medium speed channels are to be read on this particular interrupt. After the number of channels to be read is determined, CLOCK checks the buffer to insure that all the data to be converted will fit. If so, the acquisition is started. If the data will not fit, the buffers are switched and the A/D convertor is started. The systron donner times are read and the donner time convert and tape output flags are set.

In either case, control is returned to the interrupted program.

When the test is complete, STOPIT is called by an external program. This subroutine stops the A/D convertor and the clock to prevent future conversions. STOPIT also resets counters internal to CLOCK.



- 4.3.3.4.0 Subroutine TEST
- 4.3.3.4.1 The purpose of this subroutine is to control the background tasks during a data acquisition.
- 4.3.3.4.2 The calling sequence is JST \* TEST
- 4.3.3.4.3 Software/Software interfaces.

Subroutine TEST calls the following programs.

- CLOCK Called for the first data acquisition  
All additional acquisitions are done with clock interrupts.
- TDUMP Called to start the data transfer to magnetic tape and to write end of files on the tape. The A register contains the minus number of words to be transferred and the X register contains the buffer address minus one. The A register is loaded with -1 for writing end of files.
- UNPACK Called to convert the Systron Donner times from packed BCD to binary numbers. The A register contains the number to be converted.
- STØPIT Called to stop the data acquisition and reset the CLOCK routine.
- OTL Called to output messages to the teletype
- ODEC Called to convert a binary number to decimal and output it to the teletype
- CRLF Called to output a carriage return-line feed.

4.3.3.4.4 Input Data - There are no data passed as calling arguments. However, the following variables must be available in scratchpad (locations 0-255<sub>16</sub>)

LOCATION	VARIABLE	DESCRIPTION
8B <sub>16</sub>	BUFUSE	Buffer in use flag.
8E <sub>16</sub>	DØNCLK	Donner Clock time conversion flag.
80 <sub>16</sub>	CLCNT	Clock count (time between interrupts) in microseconds.

#### 4.3.3.4.5

Output Data - The following messages are output from TEST.

- 1 - DONE - ERRORS -- PARITY XXXX TIMING. XXXX
- 2 - XXXX RECORDS WRITTEN.
- 3 - RATE TOO FAST!

Message one and two give numeric values for the parity and timing errors and the number of records written.

Message three indicates that one buffer is full before the previous buffer has been completely output. This is an abort condition.

#### 4.3.3.4.6

Storage required. TEST Requires  $AB_{16}$  or  $171_{10}$  locations.

#### 4.3.3.4.7

Description. When TEST is called, the appropriate counters (parity error, timing, tape output, etc.) are set to zero and the appropriate I/O devices are set up. The clock is started and subroutine CLOCK is called to start the data acquisition. TEST then enters a loop. The first step of the loop is to check for a buffer ready for output and tape not busy. If a buffer is ready, it is output with a call to TDUMP. If not, then the second step of the loop is executed. The second step checks for Systron Donner times to be converted. If there are times to be converted, then the routine UNPACK is called. After the times are unpacked and stored for output, the test stop times are checked to see if the test has exceeded the time. If so, the test is stopped. If not, the third step of the loop is executed. This step looks for an 'X' character from the teletype. If an X has been input, the test is stopped. If an X has not been input, the loop is repeated.

After the test has been completed, subroutine TEST writes two end of files on the tape and backspaces over the second. It then informs the user of any tape errors and of the total number of records written on tape. Control is then returned to the calling program.

8D <sub>16</sub>	TPØUT	Tape output flag.
91 <sub>16</sub>	DONMS	Temporary storage for Systron Donner times.
92 <sub>16</sub>	DNSEC	
93 <sub>16</sub>	DNMIN	
94 <sub>16</sub>	DNHRS	
99 <sub>16</sub>	STPHR	Test stop times in hours, minutes, seconds, and milliseconds.
9A <sub>16</sub>	STPMIN	
9B <sub>16</sub>	STPSEC	
9C <sub>16</sub>	STØPMS	
9D <sub>16</sub>	BUFHDI	Address of the data buffers.
9E <sub>16</sub>	BUFHDO	
9F <sub>16</sub>	MBUFLG	Minus data buffer length.
A1 <sub>16</sub>	PAR	Tape parity error counter.
A2 <sub>16</sub>	TIMG	Tape timing error counter.
A6 <sub>16</sub>	TAPEND	Tape end of operation flag.

In addition, the addresses of all the external subroutines are stored in scratchpad.

LOCATION	ROUTINE
3B <sub>16</sub>	CLØCK
71 <sub>16</sub>	CRLF
72 <sub>16</sub>	ØTL
73 <sub>16</sub>	UNPACK
75 <sub>16</sub>	TDUMP
76 <sub>16</sub>	ØDEC
7B <sub>16</sub>	STØPIT



#### 4.3.3.5.0 Subroutine RESTØR

4.3.3.5.1 The purpose of this subroutine is to reset all I/O devices and to restore all scratchpad variables to their initial state.

4.3.3.5.2 Calling Sequence. The calling sequence is JST\*RESTØR

4.3.3.5.3 Software/Software Interfaces - none

4.3.3.5.4 Input Data - none

4.3.3.5.5 Output Data - The following list of scratchpad variables are restored to the indicated values.

LOCATION	VARIABLES	VALUE
80 <sub>16</sub>	CLCNT	0
81 <sub>16</sub>	NORATE	0
82 <sub>16</sub>	NØLØ	0
83 <sub>16</sub>	LØCNT	0
84 <sub>16</sub>	NØMED	0
85 <sub>16</sub>	MEDCNT	0
86 <sub>16</sub>	NØHS	0
87 <sub>16</sub>	BUFZER	77C <sub>16</sub>
88 <sub>16</sub>	BUFZND	B64 <sub>16</sub>
89 <sub>16</sub>	BUFØNE	B7D <sub>16</sub>
8A <sub>16</sub>	NDBUF1	F65 <sub>16</sub>
8C <sub>16</sub>	MAXBUF	1025
8F <sub>16</sub>	BGNCHN	0
90 <sub>16</sub>	ADCRD	0
95 <sub>16</sub>	STHP	0
96 <sub>16</sub>	STMIN	0
97 <sub>16</sub>	STSEC	0

LOCATION	VARIABLE	VALUE
98 <sub>16</sub>	STRMS	0
99 <sub>16</sub>	STPHR	0
9A <sub>16</sub>	STPMIN	0
9B <sub>16</sub>	STPSEC	0
9C <sub>16</sub>	STOPMS	0
9D <sub>16</sub>	BUFHDI	B64 <sub>16</sub>
9E <sub>16</sub>	BUFHDO	763 <sub>16</sub>
9F <sub>16</sub>	MBUFLG	-1025
A0 <sub>16</sub>	MXBDUF	1000
A3 <sub>16</sub>	HDBFLG	25
A6 <sub>16</sub>	TAPEND	1

The twenty five header words in each data buffer are also filled with zeros.

4.3.3.5.6                      Storage Required - This routine requires 41<sub>16</sub> or 65<sub>10</sub> locations.

4.3.3.5.7                      Description - The I/O devices are reset first. Then the appropriate values are restored by using constants defined in RESTOR or by constants defined by load A register immediate instructions.

4.3.3.6.0 Subroutine EINS, ZWEI, and DREI

4.3.3.6.1 The purpose of these subroutines is to input Data generic to the number of A/D conversion rates.

4.3.3.6.2 Calling Sequence - The calling sequences are -

```
I      JST * DREI
I + 1  NOP
I + 2  NOP
I + 3  RETURNS HERE
```

```
I      JST * ZWEI
I + 1  RETURNS HERE
```

```
I      JST * EINS
I + 1  NOP
I + 2  RETURNS HERE
```

4.3.3.6.3 Software/Software Interfaces - Subroutine FETCH is called to type the data request and read the input value.

4.3.3.6.4 Input Data - The only data input is through the teletype and is returned in the A register by FETCH.

The following input values are required by the indicated entry point.

#### ENTRY POINTS

#### VARIABLES REQUIRED

EINS, ZWEI, DREI

NØHS number of High speed A/D channels.

ZWEI, DREI

NØMED number of medium speed a/D channels

DREI

NØLØ number of low speed A/D channels

ZWEI, DREI

MEDCNT number of times to read high speed channels before reading medium speed channels

DREI

LØCNT number of times to read high speed channels before reading low speed channels.

4.3.3.6.5 Output Data - The Data which are teletype inputs are stored in scratchpad (locations 0-255). In addition, descriptive messages requesting input are passed to FETCH to be output to the teletype.



#### 4.3.3.6.6

Storage Required - The total storage required by the three entry points is  $71_{16}$  or  $113_{10}$  locations.

#### 4.3.3.6.7

Description - Each subroutine only types out that information which is generic to its specific rate. The X register is loaded with the address of the message to be output. Subroutine FETCH is called to type out the message and read the input number. If the number is in error (i.e., couldn't be converted from ASCII), FETCH returns to the next instruction past the call. The message loop is then repeated. If the value is valid, it is stored and the next call to FETCH is executed. All the entry points return control to the calling program.

#### 4.3.3.7.0 Subroutine HEADER

4.3.3.7.1 The purpose of this subroutine is to accept user input data for the data buffer header words 6-25.

4.3.3.7.2 The Calling sequence is JST \* HEADER

4.3.3.7.3 Software/Software Interfaces -  
Subroutine HEADER calls: CRLF to output carriage return-line feeds; OTT to output a single character from the A register; ØTL to output a message; and FETCH to output a message and to read a number for storage.

4.3.3.7.4 Input Data - The following data must be available in scratchpad -

9D <sub>16</sub>	BUFHDI	Address of the buffer header for
9E <sub>16</sub>	BUFHDO	buffer one and zero.

A4 <sub>16</sub>	BUFTTY	Address of the teletype input buffer.
------------------	--------	---------------------------------------

The only other input data is from the teletype for storage in the headers.

4.3.3.7.5 Output Data - up to 20 words of information may be stored in the buffer headers. Output message addresses are passed to FETCH and ØTL. A single character to be output is transferred to ØTT via the A register.

4.3.3.7.6 Storage Required - HEADER requires 57 or 87 locations.

4.3.3.7.7 Description - HEADER calls ØTL to output the message:

HEADER DATA INPUT TO TERMINATE, INPUT 'A/'. The subroutine FETCH is called to write the word number of the header in to which the next input value will go. If the value input is in error, the FETCH is called to print the message again. If the value is valid, it is stored in the headers. The header word pointers are incremented and FETCH is called again with the next word number. The process is repeated until either 20 values are input or the quantity "A/" is input. HEADER then returns control to the calling program.

- 4.3.3.8.0 Subroutine PFAIL
- 4.3.3.8.1 The purpose of the subroutine is to provide an orderly shutdown of the CF16A during a power failure.
- 4.3.3.8.2 There is no calling sequence since this is an interrupt subroutine.
- 4.3.3.8.3 Software/Software Interfaces  
This subroutine doesn't call any other subroutines. However, it does store a jump to the restart routine at the interrupt location,  $00_{16}$ .
- 4.3.3.8.4 Input Data - None
- 4.3.3.8.5 Output Data - An instruction is stored at the restart interrupt location
- 4.3.3.8.6 Storage Required - This routine requires  $7_{16}$  or  $7_{10}$  locations.
- 4.3.3.8.7 Description - When a power failure is detected, an interrupt is generated. This subroutine responds to that interrupt by halting an tape I/O. It also stores a jump to the restart routine at location 0. The interrupts are enable and the CF16A is halted.



- 4.3.3.9.0 Subroutine PWRØN
- 4.3.3.9.1 The purpose of PWRØN is to insure an orderly restart after a power failure.
- 4.3.3.9.2 There is no calling sequence since this routine responds to an interrupt at location 0.
- 4.3.3.9.3 Software/Software Interfaces - This subroutine calls RESTØR to reset the system variables to their original state. It also calls STØPIT to reset the CLOCK data acquisition routine to its original state.
- 4.3.3.9.4 Input Data - None
- 4.3.3.9.5 Output Data - This subroutine stores a jump to the spurious interrupt handler at location 0.
- 4.3.3.9.6 Storage Required - This routine requires  $9_{16}$  locations.
- 4.3.3.9.7 Description - When a restart interrupt occurs at location 0, this routine calls RESTØR, STØPIT, and then restores the spurious interrupt subroutine call at location 0. (When a device requests an interrupt and does not respond with an address, the CF16A interrupts to location 0.) PWRØN then begins execution at location  $111_{16}$ , the start of GDAS.

- 4.3.3.10.0 Subroutine BADRPT
- 4.3.3.10.1 The purpose of this subroutine is to alert the user that a spurious interrupt has occurred.
- 4.3.3.10.2 There is no calling sequence since this is an interrupt subroutine.
- 4.3.3.10.3 Software/Software Interfaces - This subroutine calls ØTL to output the message -  
"SPURIOUS INTERRUPT! CYCLE TO CONTINUE."
- 4.3.3.10.4 Input Data - None
- 4.3.3.10.5 Output Data - The ASCII character '.' is loaded into the A register before calling ØTL.
- 4.3.3.10.6 Storage Required - This subroutine requires 1C<sub>16</sub> or 27<sub>10</sub> locations.
- 4.3.3.10.7 Description - When a device requests an interrupt and does not respond with an interrupt address, the CF16A interrupts to location 0. This subroutine issues a halt I/O to the tape controller and then outputs a message to alert the user of the interrupt and halts. Although the user can restart at the point of interruption, this interrupt usually signifies a hardware fault and is likely to recur until the hardware is repaired.

- 4.3.3.11.0 Subroutine TPCHK
- 4.3.3.11.1 The checks the tape unit to insure that a tape is mounted and not write protected.
- 4.3.3.11.2 Calling sequence - the calling sequence is JST \* TPCHK.
- 4.3.3.11.3 Software/Software Interfaces - TPCHK calls CRLF to issue a carriage return-line feed to the teletype. Subroutine ØTL is called to output the message.
- "TAPE NOT ON-LINE."
- " TAPE IS WRITE PROTECTED."
- 4.3.3.11.4 Input Data - None
- 4.3.3.11.5 Output Data - The message termination character is transferred to ØTL via the A register. The address of the message to be output is transferred as an argument immediately after the call.
- LDA termination Character  
JST \* ØTL  
DATA message address.
- 4.3.3.11.6 Storage Required - This subroutine requires 30<sub>16</sub> or 48<sub>10</sub> locations.
- 4.3.3.11.7 Description - Using the sense instruction, TPCHK checks to insure that the tape is on-line and not write protected. If either of these conditions are not met, a message is output. TPCHK then enters a loop until the error condition is corrected. After the error conditions are corrected, TPCHK returns to the calling program.



- 4.3.3.12.0 Subroutine ADCHK
- 4.3.3.12.1 Subroutine ADCHK checks the A/D convertor for automatic mode.
- 4.3.3.12.2 Calling Sequence - The calling sequence is  
JST \*ADCHK
- 4.3.3.12.3 Software/Software Interfaces - This subroutine calls CRLF to output a carriage return-line feed and ØTL to output the message  
  
"ADC IN MANUAL MØDE."
- 4.3.3.12.4 Input Data - None
- 4.3.3.12.5 Output Data - The termination character of the message is transferred to ØTL via the A register. The address of the message is transferred as a calling argument.  
  
LDA Termination Character  
JST \* ØTL  
DATA Message Address
- 4.3.3.12.6 Storage Required - ADCHK required 16<sub>16</sub> or 22<sub>10</sub> locations.
- 4.3.3.12.7 Description - Using the sense instruction, ADCHK determines whether or not the A/D convertor is in manual mode. If so, ADCHK outputs a message to the user and waits for the A/D converter to be put in automatic mode. ADCHK returns to the calling program.